



Advance

Advancing the Quality of
Dental Care Products



Advance Dental Products from **MidwestDental**
EQUIPMENT • SERVICES • SUPPLIES

(800) 766-2025 | MidwestDental.com

Log in to
see your
pricing!

Advance Reference Catalog

Table of Contents

Acrylics & Anesthetics	6
Clear Comfort Aligners	4
Burs	8
Diamonds	10
Cements & Cosmetic Restoratives	15
Disposables	17
Endodontics & Rubber Dam	19
Evacuation	25
Impression Materials	26
Instruments	29
Infection Control & Sterilization	30
Preventive & Prophy	34
Surgical	39

*Prices are subject to change without notice.



Advance


(800) 766-2025

MidwestDental.com

06



11



16



19



33



26



36



Advancing the Quality of
Dental Care Products

Invisible
Comfortable
Convenient




HOW DOES IT WORK?

Clear Comfort Aligners is an innovative technology that allows you to straighten your teeth without the appearance of braces. It utilizes a digital 3D treatment plan and a series of custom-made clear aligners to gently and gradually adjust your smile.


WHY IS IT BETTER?

Clear Comfort Aligners are virtually invisible and allow you to smile confidently and comfortably throughout your treatment.


Clear Comfort Aligners apply less force on your teeth than braces. Because of their smooth




First a digital scan or dental impression is created of your teeth.



A series of custom-fitted Clear Comfort Aligners are made to your specifications.



The Clear Comfort Aligner system gradually straightens your teeth and corrects your bite as you wear them.



Once you have completed the Clear Comfort Aligner System, you will be left with a brilliantly corrected smile.

Advancing the Brilliance of Your Smile

surface, they tend to be more comfortable on the inside of your cheeks.

Clear Comfort Aligners are removable so you can take them out to eat, drink, brush, & floss as you normally would. You spend less time & have fewer visits at the doctor's office.



IS IT RIGHT FOR YOU?

You are a candidate for **Clear Comfort Aligners** if you are an adult or teen and have any of the following conditions:

- ✓ Crowded Teeth
- ✓ Unattractive Gaps
- ✓ Overbite
- ✓ Underbite
- ✓ Crossbite

Clear Comfort
ALIGNERS

For more information, go to:
midwestdental.com/clearcomfort





Lidocaine HCl 2% with Epinephrine 1:100,000
(Pictured above.)

- No-methylparaben formulation to avoid allergic reactions on sensitive patients.
- 100 percent latex-free components to help reduce allergic responses.
- Terminal sterilization of product/cartridges for safe use.

Box/50(032-1700)
Ea \$36.50
10+ \$35.50

Topical Anesthetic Gel
(Pictured below.)

Contains 20% benzocaine for maximum strength pain relief while fast working to quickly relieve patient discomfort. Gluten free. Available in six great tasting flavors.

1oz Bottle

- Cherry(036-1600)
- Mint(036-1605)
- Piña Colada (Dye-Free).....(036-1610)
- Bubble Gum(036-1615)
- Strawberry(036-1620)
- Raspberry.....(036-1625)

Ea \$5.50
6+ \$4.75

Temporary Crown & Bridge Material

10:1 Ratio (Pictured above.)

Syringeable bis-acryl composite for chairside provisional restorations. Available in 3 shades: A1, A2, A3. All are delivered in 76g (5ml) cartridges designed to fit on a 10:1 style automix gun. Ten waste saver tips are included per kit.


50ml Cartridge with 10 tips

- A1(168-5500)
- A2(168-5505)
- A3(168-5510)

Ea \$88.95
4+ \$85.95



Raspberry, Piña Colada, Bubble Gum, Strawberry, Cherry, Mint



Mixing Tips for Temporary Crown & Bridge

Blue Mixing Tips 10:1 & 4:1 Ratio
 Pkg/25(520-2030)

Mixing Tips 1:1 Ratio
 Pkg/25(520-2035)

Ea \$21.50

Single-Use Plastic Hub Injection Needles

NEW
Product!

- 25ga (Short, Long)
- 27ga (Short, Long)
- 30ga (X-Short, Short)



Disposable Injection Needles

Single - Use Plastic Hub (Pictured above.)

Sterile, pre-threaded plastic hub needles.
Bevel indicator on the hub for quick orientation of the needle.
Color-coded caps indicate gauge size.

- Single-use only
- Plastic hub
- For dental use only

25ga (Red) - Box/100

Short.....(600-1351)
Long.....(600-1350)

27ga (Yellow) - Box/100

Short.....(600-1353)
Long.....(600-1352)

30ga (Blue) - Box/100


X-short.....(600-1355)
Short.....(600-1354)

Ea \$11.95
3+ \$10.95




Sizing Chart


25 Gauge - Red



27 Gauge - Yellow



30 Gauge - Blue



Illustrations are not actual size.



Clinic Pack of Burs

Operative Carbide Burs

Carefully designed blade structure, rake angle, flute depth and spiral angulation combined with specially formulated tungsten carbide results in a powerful cutting performance. These burs are engineered to deliver the most efficient cutting rate and performance for the most popular procedures.

FG, pkg/10 **\$14.50**, 10+ **\$13.50**
 FG, pkg/100 **\$124.95**, 4+ **\$120.95**

RA, pkg/10 **\$16.95**, 10+ **\$15.50**
 RA, pkg/100 **\$145.50**, 4+ **\$133.50**

FGSS, pkg/10 **\$14.50**, 10+ **\$13.50**
 FGSS, pkg/100 **\$124.95**, 4+ **\$120.95**

	Round					Inverted Cone			Straight/ Flat End Plain		Taper/Flat End Plain		Amalgam Prep	Pear			
	1/2	2	4	6	8	33 1/2	35	37	56	57	169	170	245	330	332	2P	4P
FG (pkg/10)		100-2115	100-2120	100-2125	100-2130	100-2135	100-2140			100-2145			100-2150	100-2155			
FG (pkg/100)	102-1865	102-1845	102-1850	102-1855	102-1860	102-1870	102-1875	102-1880	102-1885	102-1890	102-1825	102-1945	102-1895	102-1900	102-1905	102-1815	102-1820
RA (pkg/10)		120-0300	120-0305	120-0310	120-0315												
RA (pkg/100)		120-0320	120-0325	120-0330	120-0335												
FGSS (pkg/10)														100-2110			
FGSS (pkg/100)													102-1835	102-1840			

	Straight/Flat End Fissure			Taper/Flat End Cross Plain	Straight/ Round End	Straight/ Round End Cross	
	556	557	558	701	1157	1557	1558
FG (pkg/10)	100-2160	100-2165		100-2170	100-2175		100-2180
FG (pkg/100)	102-1910	102-1915	102-1830	102-1920	102-1925	102-1935	102-1940

Turbo Carbide Burs

Pkg/10 **\$25.50**
 Pkg/100 **\$214.95**

	Pear	Straight/ Round End Cross Cut	Straight/ Flat End Cross Cut	Straight Fissure
	330	2	557	1557
FG (pkg/10)		100-2100	100-2105	
FG (pkg/100)	102-1805	102-1800		102-1810

Surgical & Endodontic Burs

Pkg/5 **\$10.50**

	Round			Straight/Flat End Cross Cut		Taper/Flat End Cross Plain	
	2	4	6	557	558	701	702
FGSL	104-1600	104-1605	104-1610	104-1615	104-1620	104-1625	104-1630

Trimming & Finishing Burs

Advance Burs are manufactured with strict quality standards. This collection of burs is created with tungsten carbide and high-quality steel, which delivers dependable strength and powerful cutting.

	Round				Egg		Football			Bullet		
	7006	7008	9006	9008	7379	9406	7404	7406	7408	8404	7802	7803
FG	136-0210	136-0220	136-0295	136-0300	136-0230	136-0305	136-0235	136-0240	136-0245	136-0290	136-0260	136-0265
Pricing	Pkg/5 21.95 4+ 20.95		Pkg/5 30.95 4+ 29.95		Pkg/5 21.95 4+ 20.95		Pkg/5 30.95 4+ 29.95		Pkg/5 21.95 4+ 20.95		Pkg/5 23.95	

	Tapered		Flat End Taper	Needle			
	7642	7664	7204	7901	7902	7903	9903
FG	136-0250	136-0255	136-0225	136-0275	136-0280	136-0285	136-0310
Pricing	Pkg/5 21.95 4+ 20.95						Pkg/5 30.95 4+ 29.95

Advance Safe Tip Series Trimming & Finishing Burs

Ideal for contouring and finishing subgingival and interproximal composite margins. Non-cutting rounded tip prevents "pink composites" due to damage to hard and soft tissue. The 10-blade trims and contours composites. Choose the 20-blade for a smooth pre-polish on composites.

(136-0200) #SE9 (10-blade), pkg/5 **25.95**
 (136-0205) #SE9 (20-blade), pkg/5 **36.95**



SE9



ADVANCE NITRILE GLOVES

Advance 2.7 Nitrile Powder-free Gloves

Hole detection: ASTM D 5151; blood-borne pathogen resistance: ASTM F 1671. 2.0mil palm thickness, 2.3mil finger thickness. Cobalt color.

Box/300
 X-Small (468-7700)
 Small (468-7705)
 Medium (468-7710)
 Large (468-7715)
 X-Large (468-7720)
Ea \$19.95
10+ (case) \$18.95

Advance 3.2 Nitrile Powder-free Gloves

Hole detection: ASTM D 5151; blood-borne pathogen resistance: ASTM F 1671. 2.5mil palm thickness, 3.3mil finger thickness. Dark lavender blue color.

Box/200
 X-Small (468-7721)
 Small (468-7722)
 Medium (468-7723)
 Large (468-7724)
 X-Large (468-7725)
Ea \$14.95
10+ (case) \$13.95

Diamonds

NEW
Multi-Use
Diamonds

Midwest Dental Multi-Use Diamonds utilize natural diamond grit for efficient cutting and less heat. Diamond particles are precisely placed to provide the most cutting surface possible, including the tips of the diamonds. The stainless steel shank is carefully designed to operate smoothly and efficiently. Made in the U.S.

100% Money Back Guarantee if not Completely Satisfied.



	Inverted Cone	Flat End Cylinder				Round End Cylinder				Round End Taper				Flat End Taper
Shape No.	805	835	835	837	830	830	850	850	855	856	856	856	856	847
Size (1/10mm)	016	010	012	012	010	014	014	018	018	014	016	018	021	012
Supercoarse		296-5561		296-5539				296-5550	296-5554		296-5553	296-5557	296-5558	
Coarse	296-5531		296-5562	296-5538	296-5532	296-5534	296-5546	296-5549		296-5547	296-5551	296-5555		
Medium					296-5533	296-5535								296-5540
Fine										296-5548	296-5552	296-5556		
Extra Fine														
Pricing	Pkg/5 . . . 18.95 , 3+ . . . 16.95													

	Flat End Taper, continued				Flame	Curette	Round				Occlusal Reduction/Finishing			
Shape No.	847	847	848	848	863	878K	801	801	801	801	379	368	368	368
Size (1/10mm)	014	016	016	018	012	016	010	014	016	023	023	016	018	023
Supercoarse					296-5524						296-5527			
Coarse	296-5542	296-5544			296-5522	296-5521		296-5504	296-5507		296-5525	296-5529	296-5528	296-5501
Medium		296-5545	296-5541				296-5503	296-5506	296-5509	296-5515				
Fine				296-5543	296-5523		296-5502	296-5505	296-5508		296-5526	296-5530		
Extra Fine														
Pricing	Pkg/5 . . . 18.95 , 3+ . . . 16.95													

	Occlusal Reduction/Finishing, continued								Mosquito	Modified Shoulder		Depth Cutter 0.5mm
Shape No.	811	811	909	909	858	858	859	852	392	837KR	850KR	834
Size (1/10mm)	033	037	035	043	008	014	014	023	016	012	018	021
Supercoarse	296-5518	296-5520								296-5565	296-5563	
Coarse	296-5517	296-5519	296-5566	296-5567	296-5510	296-5514	296-5511			296-5564		
Medium												296-5560
Fine							296-5512	296-5516	296-5536			
Extra Fine							296-5513		296-5537			
Pricing	Pkg/5 . . . 18.95 , 3+ . . . 16.95											

Single-Use Sterile Diamonds

NEW
Single-Use
Diamonds

Unique bonding technology, uniform diamond coating, larger exposure of each diamond crystal. More than 245 shapes and grits, sold in boxes of 25, compare to all leading brands.

100% Money Back Guarantee if not Completely Satisfied.



	Football		Egg		Mosquito		Round				Inverted Cone			
Shape No.	368	368	379	392	392	801	801	801	801	801	805	805	807	807
Size (1/10mm)	016	020	019	016	018	009	012	014	016	018	012	017	016	018
Supercoarse									296-5328	296-5329				
Coarse	296-5316	296-5235	296-5320			296-5325	296-5327		296-5258	296-5280		296-5281	296-5335	296-5336
Medium		296-5319			296-5228	296-5326	296-5257	296-5303	296-5279	296-5304	296-5333	296-5334		
Fine	296-5266	296-5318	296-5321	296-5256	296-5324									
Extra Fine										296-5330				
Pricing	Pkg/25 . . . 39.50 , 4+ . . . 38.50													

	Double Inverted Cone	Flat End Cylinder			Flat End Taper					Round End Taper				
Shape No.	813	835	835	837	847	847	847	848	848	850	850	856	856	856
Size (1/10mm)	014	010	012	012	014	016	018	016	018	016	018	012	014	016
Supercoarse														296-5242
Coarse	296-5305	296-5208	296-5312	296-5294	296-5211	296-5214	296-5252	296-5284	296-5314	296-5204	296-5369	296-5272	296-5230	296-5206
Medium		296-5247	296-5340											296-5232
Fine														296-5229
Pricing	Pkg/25 . . . 39.50 , 4+ . . . 38.50													

	Round End Taper				Pointed Cone					Pointed Cone Short			Flame	
Shape No.	856	856L	856L	856L	858	859	859L	883	883	858S	859S	883S	861	862
Size (1/10mm)	018	012	016	018	014	014	012	007	009	014	014	010	012	010
Supercoarse	296-5251											296-5390		
Coarse	296-5236	296-5223	296-5273	296-5380	296-5245	296-5224	296-5267			296-5388	296-5395		296-5401	
Medium					296-5295		296-5391		296-5291	296-5389	296-5263			
Fine	296-5276					296-5207	296-5285	296-5311		296-5296	296-5396	296-5274		296-5238
Extra Fine						296-5265	296-5392	296-5433		296-5255	296-5264		296-5402	
Pricing	Pkg/25 . . . 39.50 , 4+ . . . 38.50													

Advance Disposable Diamonds continue on the next page.

*Prices are subject to change without notice.

To order call **800.766.2025** or visit **MidwestDental.com** • 11

Advance Disposable Sterile Diamonds

Unique bonding technology, uniform diamond coating, larger exposure of each diamond crystal. More than 245 shapes and grits, sold in boxes of 25, compare to all leading brands.

100% Money Back Guarantee if not Completely Satisfied.



	Flame, continued			Flame Beveled Margin	Beveled Cylinder	Beveled Cyl. Short	Pear		Pear Short		Flat End Cylinder Short		
Shape No.	862	862	863	883	885	886S	830R	830R	830RS	830RS	835S	836S	837S
Size (1/10mm)	012	014	012	010	012	012	009	012	009	012	010	012	016
Supercoarse													
Coarse	296-5217	296-5309	296-5212		296-5233	296-5244		296-5282			296-5341	296-5218	296-5249
Medium	296-5297			296-5435	296-5292	296-5437	296-5260	296-5240	296-5306	296-5293	296-5342	296-5343	
Fine	296-5239		296-5259	296-5434	296-5436								
Extra Fine	296-5403												
Pricing	Pkg/25 . . . 39.50 , 4+ . . . 38.50												

	Flat End Cylinder Short Modified		Flat End Taper Short Modified					Flat End Taper Short					Torpedo Short	
Shape No.	837KRS	837KRS	845KRS	846KRS	847KR	847KRS	847KRS	847S	847S	847S	847S	848S	848S	878S
Size (1/10mm)	014	016	019	017	016	017	020	013	015	017	019	016	018	012
Supercoarse				296-5307		296-5355	296-5356		296-5360	296-5362	296-5364	296-5366	296-5368	
Coarse	296-5344		296-5347		296-5262			296-5357						296-5300
Medium			296-5348	296-5352		296-5354			296-5359	296-5283	296-5363	296-5365	296-5367	
Fine		296-5345				296-5353			296-5358	296-5361				296-5427
Pricing	Pkg/25 . . . 39.50 , 4+ . . . 38.50													

	Torpedo Short	Round End Taper Short										Flame Short		
Shape No.	878S	850S	850S	855S	855S	855S	855S	855SH	856S	856S	856S	859S	861S	862S
Size (1/10mm)	015	014	012	012	016	017	018	018	014	016	019	016	014	010
Supercoarse					296-5373	296-5241	296-5231	296-5250	296-5209	296-5205	296-5210	296-5399		
Coarse	296-5428	296-5213	296-5393	296-5308									296-5287	296-5237
Medium		296-5371	296-5394	296-5372		296-5374	296-5376	296-5378	296-5382	296-5384	296-5387	296-5398	296-5409	296-5315
Fine		296-5370					296-5375		296-5381	296-5383	296-5386	296-5397		296-5404
Extra Fine										296-5385		296-5400		
Pricing	Pkg/25 . . . 39.50 , 4+ . . . 38.50													

Advance Disposable Diamonds continue on the next page.

Advance Disposable Sterile Diamonds

Unique bonding technology, uniform diamond coating, larger exposure of each diamond crystal.
More than 245 shapes and grits, sold in boxes of 25, compare to all leading brands.

100% Money Back Guarantee if not Completely Satisfied.

	Flame Short						Pointed Taper Short							
Shape No. Size (1/10mm)	862S 012	862S 014	862S 015	862S 016	863S 012	878KS 019	877KS 017	878KS 012	878KS 014	878KS 016	878KS 018	879KS 012	879KS 016	879KS 020
Supercoarse											296-5271			296-5430
Coarse	296-5405	296-5288	296-5298	296-5268	296-5414	296-5424	296-5299	296-5269	296-5310	296-5234		296-5289	296-5290	
Medium	296-5407	296-5410		296-5412	296-5416		296-5418	296-5270	296-5420	296-5422		296-5301	296-5429	
Fine	296-5406			296-5411	296-5415	296-5425			296-5419	296-5421	296-5423			
Extra Fine	296-5277			296-5413	296-5417									
Very Fine	296-5408													
Pricing	Pkg/25 . . . 39.50 , 4+ . . . 38.50													

	Round End Cylinder KS-0 Short	Round End Cylinder KS-1 Short	Round End Cylinder KS-1	Round End Cylinder KS-2	Gross Reduction Short	
Shape No. Size (1/10mm)	880S 010	880S 010	880 010	881S 014	RC8 017	RT8S 018
Supercoarse			296-5216	296-5219		
Coarse	296-5222				296-5445	296-5278
Medium		296-5432				
Fine	296-5431					
Pricing	Pkg/25 . . . 39.50 , 4+ . . . 38.50					

	Round End Taper	Round End Taper Short	Pointed Taper Short
Shape No. Size (1/10mm)	856 022	856S 023	878KS 023
Supercoarse	296-5379	296-5254	
Coarse			296-5426
Medium			
Fine			
Extra Fine			
Pricing	Pkg/25 . . . 39.50 , 4+ . . . 38.50		

	End Cutter Short			Occlusal Reduction Short
Shape No. Size (1/10mm)	839S 012	839S 014	839S 016	811S 033
Supercoarse				296-5221
Coarse				
Medium	296-5346	296-5313	296-5243	
Fine				
Extra Fine				
Pricing	Pkg/25 . . . 48.95 , 4+ . . . 47.95			

Advance Disposable Diamonds continue on the next page.

*Prices are subject to change without notice.

Advance Disposable Sterile Diamonds

Unique bonding technology, uniform diamond coating, larger exposure of each diamond crystal. More than 245 shapes and grits, sold in boxes of 25, compare to all leading brands.


100% Money Back Guarantee if not Completely Satisfied.



Shape No. Size (1/10mm)	368 016	368 023	379 023	379 025	801 023	855 023	801S 023	845KRS 026	855S 024
Supercoarse		296-5253	296-5248			296-5226			296-5377
Coarse		296-5202	296-5201		296-5261		296-5220		
Medium		296-5302	296-5275				296-5331	296-5350	
Fine		296-5203	296-5200	296-5323			296-5246	296-5349	
Extra Fine	296-5317	296-5286	296-5322					296-5351	
Pricing	Pkg/25 . . . 42.95 , 4+ . . . 41.95								



Shape No. Size (1/10mm)	811 037	909 035	909 042	909S 040
Supercoarse	296-5215			
Coarse		296-5227	296-5225	296-5438
Pricing	Pkg/25 . . . 61.50 , 4+ . . . 60.50			



Shape No. Size (1/10mm)	815 027	834 017	834 022	998 016	998 021	998 023
Medium	296-5337	296-5338	296-5339	296-5440	296-5442	296-5444
Fine				296-5439	296-5441	296-5443
Pricing	Pkg/25 . . . 93.95 , 4+ . . . 92.95					

Advance Disposable Diamonds continued from the previous page.

*Prices are subject to change without notice.

Clear Comfort ALIGNERS

Clear Comfort Aligners are clear, custom-fitted, removable aligners that gradually straighten your teeth and correct your bite as you wear them. They are an effective and attractive alternative to braces for adults and teens with crowded teeth, excessive spacing (gaps), underbites, or overbites.

See more on page 4!




Invisible
Custom Fitted
Removable
Easy To Use
Comfortable



Call us at (800) 766-2025 or visit MidwestDental.com/clearcomfort to learn more!

Cements & Cosmetic Restoratives



Temporary Luting Cement

Compare to VOCO Provicol QM Plus!

High-retentive, eugenol-free, zinc-oxide based temporary luting cement containing calcium hydroxide. Indicated for luting temporary provisional restorations and temporary definitive restorations.

Temporary Luting Cement..... (199-1455)
Contains: 1 5ml syringe and 10 brown mixing tips.

Ea \$41.50

Mixing tips for Temporary Cement

Mixing Tip 1:1 Ratio

Pkg/20 (520-2040)
Ea \$19.50



Universal Bond Plus & Universal Bond Plus DCA

Compare to VOCO Futurabond M+!

Light-curing universal adhesive. Can be applied with or without phosphoric acid etching. Suitable for direct or indirect restorations. Provides secure adhesion, without the need for additional primer, for various materials such as metals, zirconium, or aluminum oxide, as well as silicate ceramics. Dual-cure activator for self-curing or dual curing. Applied in one layer: only 35 seconds total working time.

Universal Bond Plus 5ml (199-1453)
Contains: 1 1.5ml bond and 2 12-well mixing wells.

Ea \$116.95

Universal Bond Plus DCA
2ml (199-1454)
Ea \$51.95

LC Nano Hybrid Flowable Composite

Compare to VOCO GrandioSO Flow!

Light-cured, flowable, 80% nano-hybrid composite with a 3.5mm depth of cure. Natural aesthetics with lower shrinkage, the wear resistance of a hybrid, and excellent color retention. Radiopaque.*

*Application requires layering technique. Light-cure layer by layer in 2mm increments.

Flowable Composite - 1.8g

A1 (199-1450)
A2 (199-1451)
A3 (199-1452)
A3.5 (199-1456)
B1 (199-1457)
Contains: 1 1.8g syringe and 5 20g black applicator tips.

Ea \$21.95

Pre-Bent Needles (black, 20g)

Pkg/100 (520-2610)

Ea \$12.50

*Prices are subject to change without notice.

To order call **800.766.2025** or visit **MidwestDental.com** • 15



Pre-bent Multipurpose Dispensing Tips

Dispensing Tips are pre-bent to the proper angle and have a blunt end for patient comfort. Great for placing sealants, etchants, cements and other flowable materials. Secure Luer lock, universal color coding for easy identification.

Tips Pkg/100

18ga - Pink (520-2600)

Cavity liners, Bases & Composites

20ga - Yellow (520-2605)

Composites & Sealants

20ga - Black..... (520-2610)

Composites & Cements

22ga - Gray (520-2615)

Sealants, Etchants & Cements

25ga - Blue (520-2620)

Ea \$10.50



NEW!

Advance Articulating Paper

Advance Articulating Paper features soft, smooth, non-stick surfaces that are extra sensitive. Made in the U.S.

Pkg/12 12-sheet books

Red/blue (.003", 71 microns) (040-0002)

Thin blue (.003", 71 microns) (040-0003)

Thick blue (.005", 127 microns) (040-0004)

Ea \$12.50, 6+ \$11.50

Horseshoe Full Arch

Pkg/6 12-sheet books

Red/blue (.004", 89 microns) (040-0005)

Ea \$18.95



NEW!

Advance Matrix Strips

Advance Matrix Strips are strong, economical, and resistant to solvents. Made of polyester film, the strips will not stretch, curl, or break! Sized at 4" long x 3/8" wide x .002" thick. Made in the U.S.

Pkg/100 (564-0003)

Ea \$3.95, 6+ \$3.50

Pkg/500 (564-0004)

Ea \$7.25, 6+ \$6.75

Etching Gel 37% Phosphoric



Formulated for maximum working and handling characteristics.



Syringe Kit, each (320-4400)

Contains: 4 1.2ml syringes and 20 tips.

Ea \$20.50, 3+ \$19.95

Economy Kit, each (320-4405)

Contains: 12 1.2ml syringes. Tips sold separately.

Ea \$48.95, 3+ \$47.95

Jumbo Syringes, each..... (320-4415)

Contains: 2 25ml syringes, 60 tips, 5 3cc (3ml) empty syringes and 1 syringe connector.

Ea \$35.95

Pre-bent Tips (23ga) - Blue

Pkg/20 (320-4420)

Ea \$5.95

3cc Empty Refillable Syringes

Pkg/5 (320-4410)

Ea \$9.95



NEW!

Advance Core Build-up Material

Dual-cure flowable microhybrid composite for core buildups and cementing of posts. Polymer-based restorative material. Universal shade.

Kit (160-0002)

Contains: 2 5ml syringes, 10 mixing cannulas and 10 intra-oral tips.

Ea \$86.95



NEW!

Advance Mixing Pads

Advance Mixing Pads feature a high-gloss, poly-coated finish and a non-skid foam base. Extra heavy duty for mixing cements, composites, silicates and rubber-based materials.

3" x 3", pkg/100 (584-0000)

Ea \$6.50, 6+ \$5.95

3" x 6", pkg/100 (584-0001)

Ea \$7.95, 6+ \$7.50

6" x 6", pkg/100 (584-0002)

Ea \$10.95, 6+ \$10.50

Advance Hemostatic Brush Tip

Compare to Ultradent Metal Dento-Infusor Tips!

- Single use
- Soft bristle end allows for non-invasive brushing of solutions
- Universal luer-lock style hub
- Polyester bristle allows for spreading of viscous solutions while stiff enough to avoid metal cannula contact with soft tissue



Pkg/20 (092-2200)

Ea \$10.50

Pkg/100 (092-2201)

Ea \$55.95

Pkg/500 (092-2202)

Ea \$209.95

Advance Black Brush Tip

Compare to Ultradent Black Mini Brush Tips!

Intended for the precise and smooth delivery of aqueous solutions.

- Minimizes bubbles
- Luer Lock Hub



Brush Tips

Pkg/100 (092-2205)

Ea \$52.95

Advance Metal Brush Tips

Compare to Ultradent Black Mini Brush Tips!

- 19ga Pre-bent (Black)
- Single use
- Curved cannula for optimal entry
- Nylon brush fiber for controlled delivery of viscous material, etchants & resins
- Luer Lock Hub
- Minimizes bubbles
- Adjustable brush tip fiber length



19ga Pre-bent - Opaque Black

Pkg/100 (092-2204)

Ea \$45.95

Advance Spiral Brush Tip

Compare to Ultradent Inspiral Brush Tips!

These syringe tips are ideal for controlled placement of materials in restorative procedures.

- Luer Lock Hub
- Minimizes bubbles
- Adjustable brush tip fiber length
- Curved cannula for optimal entry
- Internal spiral channel
- Nylon brush fiber for controlled delivery of viscous material, etchants & resins



Brush Tips

Pkg/100 (092-2203)

Ea \$68.50



PREMIUM Non-Woven Sponges

- Unique rayon/polyester blend
- Provides maximum absorbency
- Less adhesion to wounds for faster healing
- Soft, Smooth, Strong, Virtually lint free
- Ideal for sustained patient use, disinfecting, general cleaning
- 4 ply non-sterile



2" x 2"
Case/5000 (456-1300)
Ea \$77.95

4" x 4"
Case/2000 (456-1310)
Ea \$109.95



BASIC Non-Woven Sponges

- Soft, viscose/poly blend, non-woven material for increased patient comfort
- Economical alternative to cotton
- No chemical additives to avoid allergic reactions
- Non-sterile

Non-Sterile - 2" x 2"
Sleeve/200 (456-1205)
Ea \$1.95

Case/5000 (456-1200)
Ea \$34.50

Non-Sterile - 4" x 4"
Sleeve/200 (456-1255)
Ea \$5.95

Case/2000 (456-1250)
Ea \$55.95



BASIC Cotton Rolls

Pliable, non-linting super absorbent rolls made from 100% cotton. They do not adhere to mucous membranes and will retain their shape and strength when wet.

Non-Sterile #2 Medium
Case/2000 (212-1100)
Ea \$14.50


Parotid Dry Angle Tips

Soft, comfortable and highly efficient, Advance Parotid Dry Angle Tips replace cotton rolls and enhance saliva absorption with their thin yet highly absorbent pads.



NEW!

Box/50
Small (312-0003)
Large (312-0002)
Ea \$16.50



Both Econoback and Deluxe Polyback Bibs come in a wide variety of colors!

Econoback Patient Bibs

Highly absorbent tissue specially treated to enhance wet strength and tear resistance, and a bonded layer of poly film. Reinforced edge to provide strength for towel clips.

2-ply tissue and 1-ply poly

13" x 19" - Case/500

Aqua	(068-2300)
Beige	(068-2305)
Blue	(068-2310)
Dusty Rose	(068-2315)
Green	(068-2320)
Lavender	(068-2325)
Peach	(068-2330)
Gray	(068-2335)
White	(068-2340)
Yellow	(068-2345)
Ea	\$32.50
4+	\$31.50

Deluxe Polyback Patient Bibs

Highly absorbent tissue specially treated to enhance wet strength and tear resistance, and a bonded later of poly film. Reinforced edge to provide strength for towel clips.

3-Ply tissue and 1-ply poly

13" x 19" - Case/500

Aqua	(068-2400)
Beige	(068-2405)
Blue	(068-2410)
Dusty Rose	(068-2415)
Green	(068-2420)
Lavender	(068-2425)
Peach	(068-2430)
Gray	(068-2435)
White	(068-2440)
Yellow	(068-2445)
Ea	\$39.95
4+	\$38.95



HEAVY weight paper!



Bracket Tray Covers

Heavy weight paper impedes the flow of moisture, minimizing liquid seep-through, protecting the covers and unit beneath the tray.

Ritter - "Size B" (8.5" x 12.25")

Box/1000

White	(804-1020)
Ea	\$29.50
3+	\$28.50

Box/1000

Blue	(804-1000)
Dusty Rose	(804-1005)
Green	(804-1010)
Silver	(804-1015)
Ea	\$40.95
3+	\$39.95



Plastic Cups

Horizontal embossed grip for easier dispensing and holding. Double rolled lip for easy separation. Can be used with any 5oz dispenser.

5oz Cups Case/1000

Blue	(260-1100)
Dusty Rose	(260-1105)
Silver Gray	(260-1110)
White	(260-1115)
Ea	\$35.50
4+	\$34.50



BASIC ProBack 1-Ply Bibs

Absorbent tissue bonded to a layer of poly film.

19" x 13" Case/500

1-ply tissue with 1-ply poly (2ply bibs)	
Blue	(068-2460)
Ea	\$35.95



BASIC Bracket Tray Covers

High quality paper tray cover.

Ritter - Size "B"

(8.5" x 12.25") Case/1000	
White	(804-1050)
Ea	\$24.95
3+	\$23.95



Compare to ProTaper® Gold and EdgeTaper Platinum™!

- Same sizes, lengths and tapers enable you to use the same protocol or technique you use with ProTaper® Gold and Edge Taper Platinum™ file systems for shaping and finishing
- Shaping files cut dentin throughout each canal zone with greater efficiency, while reducing the potential for over-preparation
- Multiple tapers and the Heat Matrix technology of the Advance Taper HT ensure added flexibility while easily cutting through dentin in specific canal zones
- Shaping files have multiple tapers that ensure flexibility and cut dentin in specific canal zones with greater efficiency



Compare to ProTaper® Universal and EdgeTaper™!

- Same sizes, lengths and tapers enable you to use the same protocol or technique you use with ProTaper® Universal and Edge Taper™ file systems for shaping and finishing
- Shaping files cut dentin throughout each canal zone with greater efficiency, while reducing the potential for over-preparation
- Shaping files have multiple tapers that ensure flexibility and cut dentin in specific canal zones with greater efficiency



Pk/6 **\$28.25**

Pk/6 **\$28.25**

Heat Treated	Order#	Size	Taper	Length	Pk Size
	328-0550	SX	04	19MM	6
	328-0551	ASST	X	21MM	6
	328-0552	F1 20	07	21MM	6
	328-0553	F2 25	08	21MM	6
	328-0554	F3 30	09	21MM	6
	328-0555	F4 40	06	21MM	6
	328-0556	F5 50	05	21MM	6
	328-0557	S1 18	02	21MM	6
	328-0558	S2 20	04	21MM	6
Better Flexibility	Order#	Size	Taper	Length	Pk Size
	328-0559	ASST	X	25MM	6
	328-0560	F1 20	07	25MM	6
	328-0561	F2 25	08	25MM	6
	328-0562	F3 30	09	25MM	6
	328-0563	F4 40	06	25MM	6
	328-0564	F5 50	05	25MM	6
	328-0565	S1 18	02	25MM	6
	328-0566	S2 20	04	25MM	6
	328-0567	ASST	X	31MM	6
328-0568	F1 20	07	31MM	6	
328-0569	F2 25	08	31MM	6	
328-0570	F3 30	09	31MM	6	
328-0571	F4 40	06	31MM	6	
328-0572	F5 50	05	31MM	6	
328-0573	S1 18	02	31MM	6	
328-0574	S2 20	04	31MM	6	

Order#	Size	Taper	Length	Pk Size
328-0500	SX	04	19MM	6
328-0501	ASST	X	21MM	6
328-0502	F1 20	07	21MM	6
328-0503	F2 25	08	21MM	6
328-0504	F3 30	09	21MM	6
328-0505	F4 40	06	21MM	6
328-0506	F5 50	05	21MM	6
328-0507	S1 18	05	21MM	6
328-0508	S2 20	04	21MM	6
328-0509	ASST	X	25MM	6
328-0510	F1 20	07	25MM	6
328-0511	F2 25	08	25MM	6
328-0512	F3 30	09	25MM	6
328-0513	F4 40	06	25MM	6
328-0514	F5 50	05	25MM	6
328-0515	S1 18	02	25MM	6
328-0516	S2 20	04	25MM	6
328-0517	ASST	X	31MM	6
328-0518	F1 20	07	31MM	6
328-0519	F2 25	08	31MM	6
328-0520	F3 30	09	31MM	6
328-0521	F4 40	06	31MM	6
328-0522	F5 50	05	31MM	6
328-0523	S1 18	02	31MM	6
328-0524	S2 20	04	31MM	6

*Assortment contains: One each of S1, S2, F1, F2, F3 and SX 19mm.

*Assortment contains: One each of S1, S2, F1, F2, F3 and SX 19mm.



Compare to WaveOne® Gold and EdgeOne Fire™!

- Same sizes, lengths and tapers enable you to use the same protocol or technique you use with Wave One® Gold and Edge One Fire™ file systems
- Shaping files cut dentin throughout each canal zone with greater efficiency, while reducing the potential for over-preparation
- Multiple tapers and the Heat Matrix technology of the Advance One HT ensure added flexibility while the optimized tip promotes safe canal navigation
- Finishing files ensure the proper apical fit for the matching Advance One or your existing Gutta Percha points



Asst. Pk/4

21mm..... (328-0654)
 25mm..... (328-0663)
 Contains: 1 each of small, primary, medium and large.

Ea \$24.50

Compare to WaveOne® and EdgeOne Platinum™!

- Same sizes, lengths and tapers enable you to use the same protocol or technique you use with Wave One® file system
- Shaping files cut dentin throughout each canal zone with greater efficiency, while reducing the potential for over-preparation
- Finishing files ensure the proper apical fit for the matching Advance One or your existing Gutta Percha points



Asst. Pk/4

21mm..... (328-0604)
 25mm..... (328-0613)
 Contains: 1 each of small, primary, medium and large.

Ea \$24.50

Pk/3 \$18.50

**Heat Treated
Better Flexibility**

Order#	Size	Length	Pk Size
328-0650	SM	21MM	3
328-0651	PR	21MM	3
328-0652	MD	21MM	3
328-0653	LG	21MM	3
<hr/>			
328-0655	SM	25MM	3
328-0656	PR	25MM	3
328-0657	MD	25MM	3
328-0658	LG	25MM	3
<hr/>			
328-0659	SM	31MM	3
328-0660	PR	31MM	3
328-0661	MD	31MM	3
328-0662	LG	31MM	3

Pk/3 \$18.50

Order#	Size	Length	Pk Size
328-0600	SM	21MM	3
328-0601	PR	21MM	3
328-0602	MD	21MM	3
328-0603	LG	21MM	3
<hr/>			
328-0605	SM	25MM	3
328-0606	PR	25MM	3
328-0607	MD	25MM	3
328-0608	LG	25MM	3
<hr/>			
328-0609	SM	31MM	3
328-0610	PR	31MM	3
328-0611	MD	31MM	3
328-0612	LG	31MM	3




100% guarantee on all Advance products!



Compare to Vortex® Blue and EdgeSequel™ Sapphire!

- Same sizes, lengths and tapers enable you to use the same protocol or technique you use with Vortex™ Blue and Edge Sequel™ Sapphire file systems
- Creates reproducible access and predictable canal patency
- Increased torque strength for safe and predictable canal navigation while the unique design enhances tactile experience
- Guide files polish dentin surfaces for easier identification of hidden canals



Pk/4 \$18.50

Order#	Size	Length	PK Size
328-0715	15/.04	21MM	4
328-0733	15/.04	25MM	4
328-0749	15/.04	30MM	4
<hr/>			
328-0713	20/.04	21MM	4
328-0731	20/.04	25MM	4
328-0747	20/.04	30MM	4
<hr/>			
328-0711	25/.04	21MM	4
328-0729	25/.04	25MM	4
328-0745	25/.04	30MM	4
<hr/>			
328-0709	30/.04	21MM	4
328-0727	30/.04	25MM	4
328-0743	30/.04	30MM	4
<hr/>			
328-0707	35/.04	21MM	4
328-0725	35/.04	25MM	4
328-0741	35/.04	30MM	4
<hr/>			
328-0405	40/.04	21MM	4
328-0723	40/.04	25MM	4
328-0739	40/.04	30MM	4
<hr/>			
328-0703	45/.04	21MM	4
328-0721	45/.04	25MM	4
328-0737	45/.04	30MM	4
<hr/>			
328-0701	50/.04	21MM	4
328-0719	50/.04	25MM	4
329-0735	50/.04	30MM	4
<hr/>			
328-0716	.04 (15-40)	25MM	4

Order#	Size	Length	PK Size
328-0714	15/.06	21MM	4
328-0732	15/.06	25MM	4
328-0748	15/.06	30MM	4
<hr/>			
328-0712	20/.06	21MM	4
328-0730	20/.06	25MM	4
328-0746	20/.06	30MM	4
<hr/>			
328-0710	25/.06	21MM	4
328-0728	25/.06	25MM	4
328-0744	25/.06	30MM	4
<hr/>			
328-0708	30/.04	21MM	4
328-0726	30/.04	25MM	4
328-0742	30/.04	30MM	4
<hr/>			
328-0706	35/.06	21MM	4
328-0724	35/.06	25MM	4
328-0740	35/.06	30MM	4
<hr/>			
328-0704	40/.06	21MM	4
328-0722	40/.06	25MM	4
328-0738	40/.06	30MM	4
<hr/>			
328-0702	45/.06	21MM	4
328-0720	45/.06	25MM	4
328-0736	45/.06	30MM	4
<hr/>			
328-0700	50/.06	21MM	4
328-0718	50/.06	25MM	4
328-0734	50/.06	30MM	4
<hr/>			
328-0717	.06 (15-40)	25MM	4

*Prices are subject to change without notice.

Advance TAPER
ABSORBENT PAPER POINTS



Advance TAPER
GUTTA PERCHA POINTS



Each Advance Taper NiTi File system includes matching:
- Paper Points
- Gutta Percha

Compare to ProTaper® and EdgeTaper™!

- Size-matched to the Advance Taper & Advance Taper HT File System
- Superior absorbency

Compare to ProTaper® and EdgeTaper™!

- Size-matched to the Advance Taper & Advance Taper HT File System
- Excellent radiopacity

Box/60 \$13.95

Paper Points	Size	Order #	PK Size
SX-F3 PP ASST	SX-F3/28	322-1500	60
F1 PP 28 MM	F1/28	322-1501	60
F2 PP 28 MM	F2/28	322-1502	60
F3 PP 28 MM	F3/28	322-1503	60
F4 PP 28 MM	F4/28	322-1504	60
F5 PP 28 MM	F5/28	322-1505	60

Box/60 \$13.95

Gutta Percha	Size	Order #	PK Size
SX-F3 GP ASST	SX-F3/28	332-9500	60
F1 GP 28 MM	F1/28	332-9501	60
F2 GP 28 MM	F2/28	332-9502	60
F3 GP 28 MM	F3/28	332-9503	60
F4 GP 28 MM	F4/28	332-9504	60
F5 GP 28 MM	F5/28	332-9505	60

Advance ONE
ABSORBENT PAPER POINTS



Advance ONE
GUTTA PERCHA POINTS



Each Advance One NiTi File system includes matching:
- Paper Points
- Gutta Percha

Compare to WaveOne® and EdgeOne™!

- Size-matched to the Advance One & Advance One HT File System
- Millimeter-marked length indicator lines
- Superior absorbency

Compare to WaveOne® and EdgeOne™!

- Size-matched to the Advance One & Advance One HT File System
- Millimeter-marked length indicator lines
- Excellent radiopacity
- Color coded size identification

Box/60 \$13.95

Paper Points	Size	Order #	PK Size
ADV ONE PP	ASST	322-1509	60
ADV ONE PP	SMALL	322-1506	60
ADV ONE PP	PRIMARY	322-1507	60
ADV ONE PP	MEDIUM	322-1508	60
ADV ONE PP	LARGE	322-1509	60

Box/60 \$13.95

Gutta Percha	Size	Order #	PK Size
ADV ONE GP	ASST	332-9510	60
ADV ONE GP	SMALL	332-9506	60
ADV ONE GP	PRIMARY	332-9507	60
ADV ONE GP	MEDIUM	332-9508	60
ADV ONE GP	LARGE	332-9509	60



Compare to Vortex® Blue and EdgeSequel™ Sapphire!

- Size-matched to the Advance Cobalt
- File System
- Superior absorbency

Box/60 \$13.95



Compare to Vortex® Blue and EdgeSequel™ Sapphire!

- Size-matched to the Advance Cobalt
- File System
- Excellent radiopacity

Box/60 \$13.95



Each Advance Cobalt Twirl NiTi File system includes matching:
 - Paper Points
 - Gutta Percha

Paper Points	Gutta Percha	Size	PK Size
322-1715	332-9715	15/.04	60
322-1714	332-9714	15/.06	60
322-1713	332-9713	20/.04	60
322-1712	332-9712	20/.06	60
322-1711	332-9711	25/.04	60
322-1710	332-9710	25/.06	60
322-1709	332-9709	30/.04	60
322-1708	332-9708	30/.06	60
322-1707	332-9707	35/.04	60
322-1706	332-9706	35/.06	60
322-1705	332-9705	40/.04	60
322-1704	332-9704	40/.06	60
322-1703	332-9703	45/.04	60
322-1702	332-9702	45/.06	60
322-1701	332-9701	50/.04	60
322-1700	332-9700	50/.06	60
322-1716	332-9716	ASST .04	60
322-1717	332-9717	ASST .06	60

*Prices are subject to change without notice.

Optimize Tissue Recovery Time with an Anti-Bacterial Endodontic Sealer

Compare to EdgeBioCeramic™ Sealer and ProRoot™ MTA Sealer!

Advance Dental Products Bioceramic Endodontic Sealer combines the proven advantages of MTA (Mineral Trioxide Aggregate) with a superior canal obturation product. Our bioceramic sealer is as biological as the actual pulp tissue, and is engineered to penetrate, fill, and seal main and lateral root canals. The 13% MTA based sealer maintains an optimal pH, which elicits anti-bacterial effects and encourages new tissue formation and recovery.

- Provides rapid recovery of tissue without causing inflammatory reaction
- Avoids tooth stain and discoloration; eugenol and bismuth oxide free
- Easy to use dual syringe automix
- Easy removal for retreatment
- No need for low temperature gutta percha cone use standard heat settings for warm vertical
- Working Time: 23 minutes, Set Time: 130 minutes

Up to 20%-40% savings over comparable leading brands per gram!



Package (362-1000)
 Contains: 1 4g syringe and 15 tips. **Ea \$59.95**

10 Automix Tips..... (362-1001)
Ea \$15.25



Advance Irrigating Tips
Slotted-end,
side-vented

Compare to **Ultradent® Endo-Eze®** and **Monoject® Notched Tip!**

- Bendable 1" needle tip designed for irrigation of canals, pockets and fistulas
- Blunt end, side port delivery needle tip
- Color coded for easier identification
- Universal Luer-style design



Pkg/100 Tips

23ga Blue	(602-0800)
27ga Yellow	(602-0801)
30ga Purple	(602-0802)

Ea \$31.50

Advance Probe Irrigating Tips
Closed-end,
side-port

Compare to **Dentsply® Max-I-Probe®!**

- Bendable 1" needle tip designed to enhance the clinician's ability to irrigate subgingival surgical sites, and sulcus
- Closed end, side port delivery needle
- For use in apical proximity
- Universal Luer-style design



Pkg/100

23ga Blue	(602-0803)
27ga Yellow	(602-0804)
30ga Purple	(602-0805)

Ea \$49.95

Advance Syringes with Needle Tips Attached
Slot-end,
side-vented
needle

Compare to **Monoject®!**

Endodontic syringe and irrigation needle are the ideal tools for irrigation of pockets, fistulas and canals.



3cc Syringe with Tip - Pkg/100

23ga x 1.25" 3cc Syringe	(744-1100)
27ga x 1.25" 3cc Syringe	(744-1101)
30ga x 1.25" 3cc Syringe	(744-1102)

Ea \$37.50

Advance Sodium Hypochlorite 6%
Compare to **Ultradent® ChlorCid®!**

For the debridement and irrigation of root canals during and after instrumentation.



Bottle - 16oz.....(570-4700)
Ea \$13.95

Advance EDTA 17% Solution
Compare to **Ultradent® EDTA®!**

EDTA is a 17% aqueous solution used for endodontic irrigation, removal of the smear layer and dentin mud from root canals.



Bottle - 16oz.....(570-4701)
Ea \$43.95

Advance Calcium Hydroxide 35%

Compare to **Ultradent® UltraCal® XS!**

- Contains no more than 35% calcium hydroxide paste
- 12.5 pH aqueous solution
- Radiopaque
- Antibacterial agent
- Stimulates healing of bone



Syringe 20-pack.....(570-4703)
Contains: 20 1.2ml syringes.

Ea \$93.50

Syringe 4-pack.....(570-4702)
Contains: 4 1.2ml syringes and 10 tips.

Ea \$36.95



Advance Ferric Sulfate Liquid

Compare to **Ultradent® Astringent™!**
15.5% Ferric Sulfate Liquid. Available in Pre-filled syringe or bottle.

30ml Bottle.....(488-1002)
Ea \$16.50

30ml Syringe.....(488-1003)
Ea \$16.50

Advance Ferric Sulfate Gel Standard Kit

Compare to **Ultradent® ViscoStat®!**



- 15.5% Ferric Sulfate Hemostatic Gel
- Provides effective tissue management when completing impressions, restorations, crown & bridge procedures
- Will not flow onto surrounding tissue

Syringe 4-pack.....(488-1000)
Contains: 4 1.2ml prefilled syringes and 8 brush tips.

Ea \$16.95

30ml Syringe.....(488-1001)
Contains: 1 30ml prefilled syringe and 1 docking port.

Ea \$26.50

Hemostatic Brush Tip Refills

Pkg/20.....(092-2200)

Ea \$10.50

Pkg/100.....(092-2201)

Ea \$55.95

Pkg/500.....(092-2202)

Ea \$209.95



Advance Aluminum Clear Hemostatic Chloride Gel
25% Aluminum Chloride Clear Hemostatic Gel

Compare to **Ultradent® ViscoStat®!**

- Surfactant lowers surface tension for better wound penetration
- Controls bleeding fast – WILL NOT STAIN tissue
- Viscous gel will not flow onto surrounding tissue

Syringe 4-pack.....(488-1004)
Contains: 4 1.2ml prefilled syringes and 8 brush tips.

Ea \$16.95

30ml Syringe.....(488-1005)
Contains: 1 30ml prefilled syringe and 1 docking port.

Ea \$26.50

Hemostatic Brush Tip Refills

Pkg/20.....(092-2200)

Ea \$10.50

Pkg/100.....(092-2201)

Ea \$55.95

Pkg/500.....(092-2202)

Ea \$209.95



Pulp Vitality Refrigerant Spray

A refrigerant spray used in dentistry for the purpose of thermo testing the vitality (existing blood flow) of the tooth prior to beginning a procedure. Peppermint flavor.

6oz Spray..... (368-0600)
Ea **\$14.95**



Advance Micro Suction Tip

Compare to Ultradent® Capillary Tips!

This narrow, flexible tip easily navigates curved canals. Intended for endodontic evacuation, though not limited to, this luer lock tip is easily attached to a luer HVE adapter. This allows for connection to high or low speed suction valves.

Pkg/50..... (092-2206)
Ea **\$19.95**

Multi-Enzyme Evacuation Cleaner

NEW!

Neutral pH formula (7.3-9.3 level)!

Bring your evacuation system and hoses to peak performance with Advance Multi-Enzyme Evacuation Cleaner. This non-foaming formula contains three powerful enzymes (protease, amylase and lipase) that dissolve debris while deodorizing your evacuation system lines with a fresh scent. Mix 1oz of Advance's concentrated cleaner with warm water to make 32oz of enzymatic solution. Also contains a corrosion inhibitor.



33.8oz Bottle.....(412-0231)
Ea **\$23.95**
4+ **\$21.95**



NEW!

Latex Dental Dams

Created from low-protein natural rubber latex that can be stretched up to 800% elongation while maintaining its strength and resisting tears. Non-toxic blue pigment. Powder free. Unscented.

Medium Gauge
(5" x 5") Box/52 (700-0003)
Ea **\$10.50**
(6" x 6") Box/36 (700-0004)
Ea **\$10.50**



NEW!

Non-Latex Dental Dams

Created from an incredibly elastic synthetic polyisoprene material that can be stretched up to 1200% elongation – all while virtually eliminating tearing even with a punched dam. Non-toxic dark blue pigment. Powder free. Unscented.

Medium Gauge
(5" x 5") Box/20 (700-0001)
Ea **\$21.50**
(6" x 6") Box/15 (700-0002)
Ea **\$21.50**

Advance Evacuation System Cleaner

NEW!

Neutral pH formula (7 level)!

All-natural, eco-friendly microbial vacuum line cleaner that breaks down organic waste and eliminates the need for daily trap cleaning. No harmful ingredients. Non-foaming formula is safe for dry pump systems, amalgam separators, and is EPA pH 6-8 compliant. Each bottle contains 64 treatments. May be used at the end of each day, or twice per week. Lemon/ginger scent. Made in the U.S.



32oz Bottle..... (412-0002)
Ea **\$57.95**

BASIC Saliva Ejectors

Vented, smooth non-removable tip. Wire-reinforced tube maintains desired shape. Phthalate free. Nonsterile. Made in Italy.




NEW!

Bag/100
White with White Tips (416-2415)
Clear with Blue Tips (416-2416)
Ea **\$4.75**
10+ **\$4.50**

Saliva Ejectors



Provides optimal suction without aspirating tissue using a unique wire-reinforced tube, which maintains any desired configuration. Latex free.




Pkg/100

White / White tip (416-2400)
Clear / Clear tip..... (416-2405)
Clear / Blue Tip (416-2410)
Ea **\$5.75**
10+ **\$5.50**

Combo Tip Evacuators (HVE)

Double-ended design with both a vented and non-vented end on the same tip. Latex free.



Pkg/100..... (416-2350)
Ea **\$4.95**
10+ **\$4.50**


Surgical Aspirator Tips

- Designed to fit all high-volume evacuation systems
- 30° angle for maximum patient comfort
- Made of rigid opaque plastic



Bag/25
1/8" x 6 1/4" White.....(416-2300)
1/4" x 6 1/4" Green.....(416-2305)
1/16" x 6 1/4" Blue(416-2310)
Ea **\$3.95**
10+ **\$3.50**

VPS Impression Material



VPS flavors:
Mint
Vanilla
Unflavored



VPS Impression Material

A hydrophilic impression material with super wetting agents that improve wettability resulting in a highly accurate impression. Excellent tear strength and dimensional stability also help to reduce the chances for distortion. Available in multiple viscosities, flavored and non-flavored in both fast and regular set. Mixing tips sold separately.

Light Body - 50mL - Pkg/4

- Mint (Fast Set).....(528-4200)
- Vanilla (Fast Set).....(528-4201)
- Unflavored (Fast Set).....(528-4202)
- Unflavored (Regular Set)(528-4203)

Ea \$40.95

Regular Body - 50mL - Pkg/4

- Mint (Fast Set).....(528-4204)
- Vanilla (Fast Set).....(528-4205)
- Unflavored (Fast Set).....(528-4206)
- Unflavored (Regular Set)(528-4207)

Ea \$40.95

Heavy Body - 50mL - Pkg/4

- Mint (Fast Set).....(528-4208)
- Vanilla (Fast Set).....(528-4209)
- Unflavored (Fast Set).....(528-4210)
- Unflavored (Regular Set)(528-4211)

Ea \$40.95

Heavy Body - 380mL - Pkg/2

- Unflavored (Fast Set).....(528-4212)

Ea \$135.95

6+ \$129.95

FAST SET	Work Time	IO Set Time
Light Body (50mL)	1:00	1:00
Medium (50mL)	1:00	1:00
Heavy Body (50mL)	1:00	1:00
Heavy Body (380mL)	1:30	1:00

REGULAR SET	Work Time	IO Set Time
Light Body (50mL)	2:30	1:00
Medium (50mL)	2:30	1:00
Heavy Body (50mL)	2:30	1:00



Putty Kit Contains (300mL)
- Base Jar
- Catalyst Jar



Alginate Substitute

An A-silicone replacement impression material that is ready to use and requires no messy mixing of powders and liquids. Ideal for preliminary impressions, orthodontic and study models and provisional and temporary impressions. Dimensionally stabilized when cured, eliminates the need to pour immediately and able to do multiple pours without loss of accuracy. Can be disinfected using standard methods. Mixing tips sold separately.

50mL - Pkg/8

- Mint (Fast Set).....(020-1510)

Ea \$59.95

380mL - Pkg/2

- Mint (Fast Set).....(020-1511)

Ea \$94.50

	Work Time	IO Set Time
Fast Set (50mL)	1:00	1:00
Fast Set (380mL)	1:15	1:00

VPS Impression Material - Putty

This hand mixable addition-cure silicone material is an ideal tray material for a 2-step putty/wash technique. This easy to mix, non-oily material when cured hardens to a high durometer which helps resist distortion and provides excellent hydraulic force when combined with a light body material.

Putty Kit contains:

- 300 mL Base Jar
- 300 mL Catalyst Jar

Putty Kit (Fast Set).....(528-4213)

Ea \$43.95

	Work Time	IO Set Time
Fast Set (300mL)	1:00	2:00



- Ideal For:
- Occlusal registration
 - Copying transfers
 - Primary impressions with bite trays



Color Changing and Dust Free



High Performance Bite Registration Material

Formulated to have a light, mousse-like consistency ideal for occlusal registration, copying transfers and primary impressions with bite trays. Our material has excellent durometer or hardness preventing distortion in disposable trays, specifically the triple trays. This unique bite registration approaches zero resistance, thus allowing for unmodified natural movement of the jaw. The non-flowing, thixotropic formula stays in place and will not drip, run, or sag.

Refill Pkg - 50ml

Contains: 2 50ml cartridges & 6 6.5mm mixing tips.

Regular Set

Work Time 2:00, Intraoral Set 1:00
 Unflavored.....(084-2605)
Ea \$22.75
4+ \$21.75

Fast Set

Work Time 0:30, Intraoral Set 1:00
 Unflavored.....(084-2610)
 Vanilla.....(084-2615)
Ea \$22.75
4+ \$21.75

Super Fast Set

Work Time 0:30, Intraoral Set 0:30
 Unflavored.....(084-2600)
 Vanilla.....(084-2620)
Ea \$22.75
4+ \$21.75

Rigid Bite - Fast Set

Unflavored.....(084-2625)
Ea \$22.75
4+ \$21.75

Clear Bite - Fast Set

Unflavored.....(084-2630)
Ea \$29.50
4+ \$28.50

Mixing Tips for VPS Impression Material

Pkg/48

Large - Green (6.5mm).....(520-2000)
 Small - Yellow (4.2mm).....(520-2005)
 Medium - Pink (5.4mm).....(520-2010)
 X-Large - Purple (7.5mm).....(520-2015)
Ea \$27.95

Pkg/100

Intra-oral Clear (fits 6.5mm).....(520-2020)
 Intra-oral Yellow (fits 4.2mm).....(520-2025)
Ea \$17.95

Dust Free Alginate *Mint*

- Color changing and dust free.
- Easy to use, avoiding remakes.
- Shorter duration in the mouth offers maximum patient comfort.
- Pleasing mint flavor.

1 lb. Bag - Color Changing

Fast Set (Mint).....(020-1460)
Ea \$13.50
10+* (case) \$11.50

1 lb. Bag

Regular Set (Mint).....(020-1451)
 Fast Set (Mint).....(020-1450)
Ea \$13.50
10+* (case) \$11.50

*Indicates a case which includes 10 1lb pouches, one canister, one alginate powder scoop and one vial.





NEW & IMPROVED DESIGN



Premium Disposable Bite Trays

- Triple function trays
- Rigid frame with nylon webbing to ensure accurate bite
- Ergonomically designed thumb grip aids in placement and removal
- Patent Pending grooves provide superior material retention

Ea \$39.95
4+ \$38.95

Box/50
Posterior (Blue)..... (524-6020)
Sideless (White)..... (524-6010)

Box/35
Anterior (Purple) (524-6015)
Quadrant (Yellow)..... (524-6005)
Assorted Sizes* (524-6025)
*Assorted contains: 7 of each size.

Box/30
Full Arch (Pink)..... (524-6000)

Multi-color Disposable Impression Trays

- Strategically perforated trays allow easy modification with dental bur to accommodate implants
- Can be heated and molded for a custom fit
- Bright color-coded body for quick and easy identification during procedure
- Easy to apply I.D. labels included in every bag
- Non-sterile

Ea \$6.95
6+ \$5.95

Large Perforated Yellow - Bag/12

#1 - Upper (524-6030)
#2 - Lower (524-6035)

Medium Perforated Blue - Bag/12

#3 - Upper (524-6040)
#4 - Lower (524-6045)

Small Perforated Orange - Bag/12

#5 - Upper (524-6050)
#6 - Lower (524-6055)

Partial Perforated Green - Bag/12

#7 - UL/LR..... (524-6060)
#8 - UR/LL..... (524-6065)

Anterior Perforated Pink - Bag/12

#9 - Upper (524-6070)
#10 - Lower (524-6075)

Assorted Perforated

Bag/12 (524-6080)
Contains: 4 each of medium, and 2 of each other size.

Premium Disposable Impression Trays

- Rigid high-impact styrene with smooth rounded edges
- Interior retention slots and vertical perforation enhance material retention
- Generates a consistently detailed reproduction with dimensionally-stable results

Ea \$10.95
6+ \$9.95

Perforated Upper - Bag/12

#1 - Large (524-6100)
#3 - Medium..... (524-6101)
#5 - Small (524-6102)

Perforated Lower - Bag/12

#2 - Large (524-6103)
#4 - Medium..... (524-6104)
#6 - Small (524-6105)

Perforated - Bag/12

#7S - UR/LL Perforated Partial..... (524-6110)
#8S - UL/LR Perforated Partial..... (524-6111)
#9S - Lower Perforated Anterior..... (524-6112)
#10S - Upper Perforated Anterior.... (524-6113)

Solid - Bag/12

#7 - Partial UR/LL Solid (524-6106)
#8 - Partial UL/LR Solid (524-6107)
#9 - Anterior Lower Solid (524-6108)
#10 - Anterior Upper Solid (524-6109)



Instruments



Explorers

Explorers are used to detect caries and calculus. The knurled handle provides a solid, ergonomic grip for greater tactile sensitivity.

- Single Ended - #23 Shepherds Hook (542-2003)
Ea \$12.50, 6+ \$11.95
- Double Ended - #2 (542-2004)
Ea \$15.50, 6+ \$14.95
- Double Ended - #5 (542-2005)
Ea \$15.50, 6+ \$14.95

Expros

Expros have an explorer on one end and a probe on the other. This ergonomic design allows the practitioner to examine and evaluate the oral cavity with more precision without putting the instrument down.

- #23/CP11 Color Coded - 3,6,8,11 (542-2006)
Ea \$16.50, 6+ \$15.95
- #23/CP12 Color Coded - 3,6,9,12 (542-2007)
Ea \$17.95, 6+ \$17.50
- #23/Williams Color Coded - 1,2,3,5,7,8,9,10 (542-2008)
Ea \$18.50, 6+ \$17.50

Probes

Probes are used to assess pocket depth and measure dental anomalies and lesions during diagnostic examinations.

- #CP12 Color Coded - 3,6,9,12 (542-2009)
Ea \$15.25, 6+ \$14.75
- Williams Round - 1,2,3,5,7,8,9,10 (542-2010)
Ea \$15.25, 6+ \$14.75

College Pliers

Our college pliers can be used to grasp gauze or other materials from the tray or oral cavity. Precision tips and our two-stage guide pin ensure accuracy and prevent glove punctures when closed.

- Serrated 6" (542-2000)
Ea \$9.95, 6+ \$9.50
- Self-Locking X-Fine Serrated 6" (542-2001)
Ea \$15.50, 6+ \$14.95
- Narrow Tip #317B 6" (542-2002)
Ea \$15.50, 6+ \$14.95



Mirrors

Front surface mirrors give a true reflection, eliminating the ghost or double image seen on plane front mirrors. The stainless steel housing on this mirror offers a longer life than typical brass housings.



Pkg/12 Front Surface Cone Socket

- #4 (572-2000)
 - #5 (572-2002)
- Pkg/12 \$23.50, 3+ \$22.50**

Pkg/12 Front Surface Simple Stem

- #4 (572-2001)
 - #5 (572-2003)
- Pkg/12 \$23.50, 3+ \$22.50**

Mirror Handles

Our solid construction knurled mirror handles are the perfect complement to our front surface mirrors, offering a tight fit and tactile connection to the instrument.

- Round Mirror Handle - Simple Stem (572-2007)
- Round Mirror Handle - Cone Socket (572-2006)
- Octagonal Mirror Handle - Simple Stem (572-2005)
- Octagonal Mirror Handle - Cone Socket (572-2004)
Ea \$3.95, 12+ \$2.95

Infection Control and Sterilization.

CSR Sterilization Wrap

Advance Sterilization Wrap brings an affordable balance of drapeability and softness with excellent bacterial barrier and strength.



Cellulose Fibers

Case/500 Sheets

12" x 12"	(546-1400)	Ea \$88.95
15" x 15"	(546-1401)	Ea \$56.95
18" x 18"	(546-1402)	Ea \$113.50
20" x 20"	(546-1403)	Ea \$126.50
24" x 24"	(546-1404)	Ea \$181.95

Advance Biological Mail-in Monitoring Service



Provides third-party documentation of testing and results with a tangible online record when proof must be provided. Compatible with steam, ethylene oxide, chemical vapor and dry heat sterilizers. Steam results are available in 24 hours, chemical vapor in 72 hours, and dry heat and ethylene oxide in 7 days (from start of incubation period at the processing laboratory). Complimentary reminder service notifies practice of missed test date by either fax or email. Made in the U.S.

Single Strip Value Pack

(056-0190) **Ea \$223.95**
Contains: 52 sets of one test strip and one control strip. Postage not included.

2 Strips Tests

(056-0191) **Ea \$112.50**
Contains: 12 sets of two test strips, one control strip and one pre-paid postage.

Enzymatic Ultrasonic Cleaner

NEW!



Effervescent Tablets

Enzymatic ultrasonic cleaning solution that is available in easy-to-use tablets which are also effective for cleaning the evacuation system and pre-soaking instruments.



64 Tablets..... (812-0264)
Ea \$51.95

General Purpose Cleaner

NEW!



Non-ammoniated concentrated liquid formula employs surfactants and anti-corrosive additives to safely and thoroughly clean instruments prior to sterilization. Phosphate-free.

Gallon..... (812-0260)
Ea \$11.95
4+ \$10.95

Tartar & Stain Remover

NEW!



Powerful, ready-to-use ultrasonic solution to remove stains from dentures, partials, appliances. Also removes temporary cement.

Gallon..... (812-0261)
Ea \$11.95
4+ \$10.95

Temporary Cement Remover

NEW!



Strong ready-to-use ultrasonic solution to remove temporary cement.

Gallon..... (812-0262)
Ea \$11.95
4+ \$10.95

Plaster & Stone Remover

NEW!



For removal of plaster and stone from bridges, crowns, dentures, inlays, metal impression trays, spatulas, etc.

Gallon..... (812-0263)
Ea \$11.95
4+ \$10.95

Glutaraldehyde 2.65% Sterilant

NEW!



Ready to use sterilizing and disinfecting solution. No activator required.

Gallon..... (300-3710)
Ea \$35.50
4+ \$30.50



Premium Sterilization Pouches

Autoclavable self-sealing, blue tinted, see-through pouches. Easy-to-read color changing arrows are indicators for chemical vapor, steam or gas sterilization. Utilizes a triple seal plus a fourth "Secure Seal" at the chevron.

- 3.5" x 9", box/200 (048-2400)
Ea \$11.95, 6+ \$10.95
- 2.75" x 9", box/200 (048-2410)
Ea \$10.50, 6+ \$9.95
- 2.25" x 4", box/200 (048-2425)
Ea \$8.95, 6+ \$8.50
- 3.5" x 5.25", box/200 (048-2420)
Ea \$13.95, 6+ \$12.95
- 5.25" x 10", box/200 (048-2415)
Ea \$18.50, 6+ \$17.50
- 3.5" x 9", box/500 (048-2405)
Ea \$32.50

*Quantity rates apply to each size only.



BASIC Self-Sealing Sterilization Pouches

- Dual External Indicators, Built-in Internal Indicator
- Blue-Tinted, Transparent Film
- Triple Seal & Heat Processed
- Wide Self-Seal Adhesive Strip

- 3.5" x 9", box/200 (048-2470)
Ea \$9.95, 6+ \$9.50
- 2.75" x 9", box/200 (048-2455)
Ea \$7.95, 6+ \$7.50
- 2.25" x 4", box/200 (048-2450)
Ea \$5.50, 6+ \$4.95
- 3.5" x 5.25", box/200 (048-2465)
Ea \$7.50, 6+ \$6.95
- 5.25" x 10", box/200 (048-2480)
Ea \$13.95, 6+ \$12.95
- 7.5" x 13", box/200 (048-2485)
Ea \$25.50
- 3.5" x 9", box/1,000 (048-2475)
Ea \$43.95
- 2.75" x 9", box/1,000 (048-2460)
Ea \$37.50

*Quantity rates apply to each size only.

*Prices are subject to change without notice.



Chair Sleeves

Tapered design conforms neatly to patient chairs, covering chair back

- Box/225**
27.5" x 24" (304-3725)
Ea \$33.95
3+ (case) \$33.50



Air/Water Syringe Sleeve

Syringe sleeves have pre-cut openings for easy application. Made of clear poly material that protects instruments from cross-contamination. Clear.

- Box/500**
21/2" x 10" (304-3720)
Ea \$9.25



T-Handle Sleeves

Specially designed barriers effectively cover equipment with near custom fit.

- Box/500 T-Style**
4" x 5 3/4" (304-3705)
Ea \$14.95



Tray Sleeves

Back-flip closure ensures tray is completely covered. Made of a clear poly material that protects tray from cross-contamination.

- Box/500**
10 1/2" x 14", Size B (304-3700)
Ea \$23.50



X-Ray Cover

Specially designed barriers effectively cover equipment with near custom fit.

- Box/250**
15" x 26" (304-3730)
Ea \$28.95

Headrest Cover

Tapered design conforms neatly to patient chairs, covering chair back




- 9.5" x 11", box/250 (484-0700)
Ea \$13.50, 4+ \$12.95
- 9.5" x 14", box/250 (484-0701)
Ea \$14.95, 4+ \$14.50



Barrier Film

4" x 6" perforated sheets with non-adhesive side edges for easy removal. Specially formulated adhesive leaves no residue behind. Helps to prevent cross-contamination. Convenient dispenser box contains 1 roll of 1200 perforated. 4" x 6" sheets.

- Roll/1200**
Clear (304-3710)
Blue (304-3715)
Ea \$18.50



Air/Water Syringe Tips

- Rounded edge for maximum patient comfort
- Convenient dispenser bottle keeps work area clear and assists in preventing cross-contamination
- Rounded non-scratch end to minimize wear and tear of O-ring in hand piece

- Assorted Colors, bag/250 (752-0800)
Ea \$30.95
3+ \$28.95

WE GOT YOU COVERED.



Premium Earloop Mask

Level 1

Procedural mask with soft white facial tissue inner layer, fluid resistant outer facing, and a front (outside) earloop attachment to eliminate potential irritation. Filtration efficiency of 99.5% at 1 micron. Latex Free.

Box/50

Blue	(418-3200)
Pink	(418-3205)
Lavender.....	(418-3210)
Ea	\$9.95
5+	\$8.95



BASIC Earloop Face Mask

Level 1

Level 1 protection for low-risk procedures. Ideal for patient appointments and exams.

- Inner and outer layer constructed for fluid resistance, with ultrasonically sealed, soft, stretchy, latex-free earloops for superior comfort.
- Adjustable nosepiece for custom fit and easy breathing.

Level 1 Mask - Box/50

Blue	(418-2700)
Pink	(418-2705)
Lavender.....	(418-3710)
Ea	\$7.50
5+	\$6.50



BASIC Earloop Face Mask

Level 3

Level 3 protection for a bit higher-risk procedures. Comes in white only.

- Inner and outer layer constructed for fluid resistance, with ultrasonically sealed, soft, stretchy, latex-free earloops for superior comfort.
- Adjustable nosepiece for custom fit and easy breathing.

Level 3 Mask - Box/50

White	(418-2750)
Ea	\$7.95
5+	\$6.95



Advance 2.7 Nitrile Powder-free Gloves

Hole detection: ASTM D 5151; blood-borne pathogen resistance: ASTM F 1671.

Thinner, stronger and far more elastic than normal nitrile gloves. Textured grip. Tested for use with chemo drugs. Gloves are flat-packed for easy dispensing. 2.0mil palm thickness, 2.3mil finger thickness. Cobalt color.

Box/300

- X-Small (468-7700)
- Small (468-7705)
- Medium (468-7710)
- Large (468-7715)
- X-Large (468-7720)

Ea \$19.95
10+ (case) \$18.95



Advance 3.2 Nitrile Powder-free Gloves

Hole detection: ASTM D 5151; blood-borne pathogen resistance: ASTM F 1671.

Thinner, lighter nitrile glove. Textured fingertips. Tested for use with chemo drugs. Gloves are flat-packed for easy dispensing. 2.5mil palm thickness, 3.3mil finger thickness. Dark lavender blue color.

Box/200

- X-Small (468-7721)
- Small (468-7722)
- Medium (468-7723)
- Large (468-7724)
- X-Large (468-7725)

Ea \$14.95
10+ (case) \$13.95



NEW!

Advance 3.0 Ultra Black Nitrile Powder-free Gloves

Hole detection: ASTM D 5151; blood-borne pathogen resistance: ASTM F 1671.

Soft, relaxed fit for added comfort and reduced muscle fatigue during longer exams. Textured finish for a stronger grip. Gloves are flat-packed for easy dispensing. 2.0mil palm thickness, 2.3mil finger thickness. Ultra black color.

Box/200

- X-Small (468-7726)
- Small (468-7727)
- Medium (468-7728)
- Large (468-7729)
- X-Large (468-7730)

Ea \$15.95
10+ (case) \$14.95



Advance Aqua CR Chloroprene Powder-free Gloves

Hole detection: ASTM D 5151; blood-borne pathogen resistance: ASTM F 1671.

Advance Aqua CR Chloroprene Gloves are soft, stretchy, and strong. Made from polychloroprene synthetic polymer, Advance Aqua CR provides a textured finish for excellent wet or dry gripping ability. Gloves are flat-packed for easy dispensing. 2.8mil palm thickness, 3.5mil finger thickness. Aqua color.

Box/200

- X-Small (468-7750)
- Small (468-7755)
- Medium (468-7760)
- Large (468-7765)
- X-Large (468-7770)

Ea \$21.95
10+ (case) \$20.95



COMING 2023!

Advance Pink CR Chloroprene Powder-free Gloves

Hole detection: ASTM D 5151; blood-borne pathogen resistance: ASTM F 1671.

Advance Pink CR Chloroprene Gloves are soft, stretchy, and strong. Made from polychloroprene synthetic polymer, Advance Pink CR provides a textured finish for excellent wet or dry gripping ability. Gloves are flat-packed for easy dispensing. 2.8mil palm thickness, 3.5mil finger thickness. Pink color.

Box/200

- X-Small (468-7775)
- Small (468-7780)
- Medium (468-7785)



NEW!

Advance Latex Powder-free Gloves

Hole detection: ASTM D 5151; blood-borne pathogen resistance: ASTM F 1671.

Advance Latex Gloves are micro-chlorinated for softness and comfort. They're slightly thinner and have a textured finish for a stronger grip. The lower protein content of Advance Latex helps reduce the chances for allergic reactions. Gloves are flat-packed for easy dispensing. 4.0mil palm thickness, 4.5mil finger thickness. White color.

Box/100

- X-Small (468-7740)
- Small (468-7741)
- Medium (468-7742)
- Large (468-7743)
- X-Large (468-7744)

Ea \$10.95
10+ (case) \$9.95

*Prices are subject to change without notice.



Advance Technical Adult

Adult 30 Tuft

Assorted: 4 assorted colors

No Imprint, Box/72.....(776-8741)
Ea **\$42.95**

With Imprint, Box/144*.....(777-0201)
Ea **\$131.95**



Adult 38 Tuft

Compare to Reach® Toothbrush!

Assorted: Yellow, Red, Dark Blue, Dark Green

No Imprint, Box/72.....(776-8732)
Ea **\$38.95**

With Imprint, Box/144*.....(777-0105)
Ea **\$123.95**



Adult 38 Tuft

With indicator.

Assorted: Blue, Dark Green, Orange, Bright Purple

No Imprint, Box/72.....(776-8729)
Ea **\$43.95**

With Imprint, Box/144*.....(777-0035)
Ea **\$118.50**



Adult 32 Tuft

Interdental with Tongue Cleaner and Power Tip

Assorted: Green, Purple, Red, Blue

No Imprint, Box/72.....(776-8738)
Ea **\$43.95**

With Imprint, Box/144*.....(777-0200)
Ea **\$128.95**



Adult 35 Tuft

Compare to Colgate® Advantage Plus Toothbrush!

Interdental Cut Tuft with Super Grip

Assorted: Dark Purple, Dark Blue, Teal, Orange

No Imprint, Box/72.....(776-8734)
Ea **\$38.95**

With Imprint, Box/144*.....(777-0155)
Ea **\$122.95**



Prepasted Disposable

Adult 39 Tuft (Mint Flavor). Cannot be imprinted.

Assorted: Yellow, Red, Green, Blue

No Imprint, Box/100.....(776-8751)
Ea **\$39.50**

10+ **\$38.50**



Adult 33 Tuft

Assorted: Blue, Green, Pink, Purple

No Imprint, Box/72.....(776-8742)
Ea **\$43.95**

With Imprint, Box/144*.....(777-0202)
Ea **\$132.95**



Adult 35 Tuft

Soft Rounded Tuft

Compare to Colgate® Advantage Compact Head Toothbrush!

Assorted: Purple, Dark Blue, Dark Green, Orange

No Imprint, Box/72.....(776-8733)
Ea **\$38.95**

With Imprint, Box/144*.....(777-0140)
Ea **\$109.95**



Adult 40 Tuft

Interdental Cut

Assorted: Dark Blue, Red, Dark Green, Yellow

No Imprint, Box/72.....(776-8731)
Ea **\$38.95**

With Imprint, Box/144*.....(777-0085)
Ea **\$123.50**



Adult 36 Tuft

Gum Stimulator and Tongue Cleaner

Assorted: Dark Pink, Bright Blue, Vibrant Green, Orange

No Imprint, Box/72.....(776-8735)
Ea **\$88.95**

With Imprint, Box/144*.....(777-0170)
Ea **\$128.95**



Denture Brush

Dupont Tynex® 12mil, dual-headed bristles are attached to a dishwasher-safe, angled handle with finger grips. Individually wrapped. Cannot be imprinted.

Assorted: Red, Green, Purple, Blue

No Imprint, Box/12.....(776-8750)
Ea **\$12.95**



Adult Basic Twist

Compare to Colgate® Wave Compact Head Toothbrush!

Compact, Extra-soft Bristles

Assorted: Red, Green, Blue, Purple

No Imprint, Box/72.....(776-8737)
Ea **\$25.95**

With Imprint, Box/144*.....(777-0201)
Ea **\$131.95**



Adult Basic Straight Handle

Adult 39 Rounded Tuft

Assorted: Red, Dark Blue, Dark Purple, Dark Green

No Imprint, Box/72.....(776-8726)
Ea **\$27.95**

With Imprint, Box/144*.....(777-0010)
Ea **\$109.95**



Adult Basic Curved Handle

Adult 38 Soft Tuft Rounded

Assorted: Dark Blue, Dark Purple, Dark Green, Orange

No Imprint, Box/72.....(776-8728)
Ea **\$38.50**

With Imprint, Box/144*.....(777-0025)
Ea **\$109.95**



Youth Basic Compact Head

Straight Handle

Young Adult 26 Soft Tuft

Assorted: Orange, Red, Dark Blue, Aqua

No Imprint, Box/72.....(776-8725)
Ea **\$27.95**

With Imprint, Box/144*.....(777-0005)
Ea **\$109.95**



*Please allow extra time for imprinting and delivery.

YOUTH/CHILD BRUSHES

YOUTH



Advance Color Youth

Compare to Oral-B® Stages 4!

Age 8+ 32 Tuft

Assorted: Two Tone colors

No Imprint, Box/72..... (776-8739)

Ea \$35.25

With Imprint, Box/144..... (777-0060)

Ea \$114.95

CHILD



Ladybug Grip

Compare to Oral-B® Stages 1!

Infant 24 Soft Tuft

Assorted: 4 assorted colors

Box/72..... (776-8740)

Ea \$25.95

With Imprint, Box/144..... (777-0135)

Ea \$106.95

CHILD



Child's Sparkle

Ages 5-7 27 Rounded Tuft

Assorted: 5 assorted two-tone colors with sparkles

No Imprint, Box/72..... (776-8730)

Ea \$27.95

With Imprint, Box/144..... (777-0055)

Ea \$92.50

CHILD



Basic

Child's Curved Handle

Ages 3-5 23 Soft Tuft

Assorted: Red, Green, Blue, Purple, White

No Imprint, Box/72..... (776-8727)

Ea \$21.95

With Imprint, Box/144..... (777-0020)

Ea \$90.95

*Please allow extra time for imprinting and delivery. Prices are subject to change without notice.



Denture Boxes

- One-piece "hinged" lid provides easy access to contents inside
- Sturdy plastic helps protect dentures from damage
- Deep soaking well fits both uppers and lowers at the same time
- Assorted colors 3¾" x 3" x 2¾"

Pkg/12.....(704-0401)

Ea \$15.95

4+ \$14.50



Retainer Boxes

- Easy to open/close "snap-lock"
- Sturdy plastic protects from dust and damage
- Vented holes top and bottom
- Assorted colors 3" x 3½" x 1½"

Pkg/12.....(704-0400)

Ea \$15.95

4+ \$14.50



Advance Thick Flexible Interdental Pick



Advance Thin Flexible Interdental Pick


Advance Flexible Interdental Picks

Cylinder interproximal picks. Perfect for cleaning between teeth and under orthodontic work.

72 packs of 10 (Thick)..... (776-8701)

72 packs of 10 (Thin)..... (776-8700)

Ea \$74.95



Advance Interdental Brushes

Great for removal of food and plaque on bridges, crowns, braces and between teeth.

Wide Cone - Blue Grip

Box/72..... (776-8705)

Moderate Cylinder - Red Grip

Box/72..... (776-8704)

Ultra Tight - Yellow Grip

Box/72..... (776-8702)

Tight Cylinder - Green Grip

Box/72..... (776-8703)

Ea \$26.95



NEW
Product
Coming
2023!

TruGuard 5000

TruGuard 5000 is a sodium fluoride dentifrice formulated for the treatment of patients best suited to receive single step additional fluoride incorporated into their everyday home oral care regimen. TruGuard 5000 also contains 5% potassium nitrate that will stop decay and sensitivity.



1.1% Neutral
Sodium Fluoride
5000 ppm
Toothpaste with
5% Potassium
Nitrate.

**COMING
2023!**



TruGuard 5000

Compare to Colgate® PreviDent® 5000 ppm Booster Plus Toothpaste!

1.1% Neutral Sodium Fluoride 5000 ppm Toothpaste with 5% Potassium Nitrate. Sweetened with Xylitol. Does not Contain Sodium. Fresh mint flavor.

- Reduces Caries
- Relieves Hypersensitivity
- Reduces Plaque Formation
- Strengthens Tooth Enamel
- Control root surface caries
- Reduces Dry Mouth
- Provides Gentle Cleaning

4 oz Tube. COMING SOON!



SonicTouch

Rechargeable Sonic Power Toothbrush

Compare to Phillips® Sonicare®!

- Produces 40,000 brush strokes per minute
- Built-in 2 minute timer
- Inductive charging
- Protection from direct contact with water
- Built-in 30 second interval timer
- 4 brushing modes
- 3 brush heads included

SonicTouch.....(776-1700)
Ea **\$54.95**

Contains: 1 toothbrush and 3 brush heads.



SonicTouch Replacement Brushes

- The rubberized brush head is designed to protect tooth enamel.
- The brush conforms to the contours of your teeth for effective, deep cleaning.
- High quality DuPont bristles promote naturally whiter teeth by gently polishing off surface stains, end-rounded.

Replacement Brushes

Pkg/2.....(776-1701)
Ea **\$9.95**

Advance SonicTouch

Rechargeable Sonic Power Toothbrush



Advanced Sonic Cleaning Technology

Produces 40,000 brush strokes per minute for a powerful dynamic cleaning.



2 Minute Auto-Timer

Built-in timer to ensure dentist recommended brushing time.



Base Charging Indicator

Inductive charging. Extra long usage time, Li-ion rechargeable battery inside.



Waterproof Design

Protection from direct contact from water, spills or splashes.



30 Second Interval


Built-in interval timer helps you brush each area of your mouth evenly.



4 Brushing Modes:

Clean, Soft, Whiten, Massage

Choose from Clean, Soft, Whiten, or Massage modes to cater all of your oral care needs.



Advance SonicTouch Replacement Brushes

(776-1701)

Ea **\$9.95**

Contains: 2 brush heads.

Compare to Phillips® Sonicare®!

Advance SonicTouch Rechargeable Toothbrush

(776-1700)

Ea **\$54.95**

Contains: 1 toothbrush and 3 brush heads.



Preventive & Prophylaxis



Varnish Comes in two popular flavors – Mint & Toffee!



Advance Varnish

5% Sodium Fluoride

Compare to VOCO Profluorid Varnish!
Easy non-messy single dose delivery system. Tooth-shaded color. Contains Xylitol. Tasty and transparent 5% sodium fluoride.

Box/50

0.4ml unidose (Toffee).....(452-6110)
0.4ml unidose (Mint).....(452-6111)

Ea \$53.50
4+ \$52.50



COMING EARLY 2023!



Basic Disposable Propyl Angles

Latex Free

Comfortable, ergonomic design with non-slip grip for reduced hand fatigue. Runs smoothly and quietly. Latex-free cups are webbed and ribbed. Made in the U.S.

Box/100

Soft, gray.....(670-4020)
Firm, white.....(670-4021)

Ea \$32.95
5+ \$31.95



Propyl Paste with Fluoride

Gluten-free, non-splatter formula. Each cup contains 1.23% fluoride ion.

Mint - Box/200

X-Coarse.....(684-7035)
Coarse.....(684-7030)
Medium.....(684-7025)

Bubblegum - Box/200

Coarse.....(684-7020)
Medium.....(684-7015)
Fine.....(684-7010)

Raspberry - Box/200


Coarse.....(684-7045)
Medium.....(684-7040)

Assorted - Box/200*

Coarse.....(684-7005)
Medium.....(684-7000)

Ea \$19.95
3+ \$18.95

*Assorted contains bubblegum, mint and raspberry.



Fluoride Trays

Fluoride trays are color-coded for easy identification and are disposable for your convenience. The Single Arch trays are made of soft, light foam for patient comfort and have a convenient handle for easy placement and removal. The Dual Arch trays have a fill line that helps reduce gel overflow and locking handles for easy placement and removal.

Single Arch Pkg/100

Small, white.....(448-0510)
Medium, yellow.....(448-0515)
Large, blue.....(448-0520)

Ea \$10.50
3+ \$9.95

Dual Arch Pkg/100

Small, white.....(448-0400)
Small/Medium, teal.....(448-0405)
Medium, yellow.....(448-0410)
Large, blue.....(448-0415)

Ea \$27.95
3+ \$26.95




Dual Arch Pkg/50

Small, white.....(448-0460)
Small/Medium, teal.....(448-0455)
Medium, yellow.....(448-0450)
Large, blue.....(448-0500)

Ea \$15.95
3+ \$14.95



18" Sutures

18" Plain Gut Sutures with Needles

Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C-6	592B	Absorbable	Plain Gut	3/8 circle -19mm reverse cut	3/0	18"	740-3040
	C-6	591B	Absorbable	Plain Gut	3/8 circle -19mm reverse cut	4/0	18"	740-3020
	C-31	594B	Absorbable	Plain Gut	1/2 circle -24mm reverse cut	3/0	18"	740-3010





Pkg/12 **\$34.50**

18" Chromic Gut Sutures with Needles

Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C-6	525B	Absorbable	Chromic Gut	3/8 circle -19mm reverse cut	3/0	18"	740-3030
	C-6	558B	Absorbable	Chromic Gut	3/8 circle -19mm reverse cut	4/0	18"	740-3035

Pkg/12 **\$34.50**



18" Silk Braided Sutures with Needles

Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C-6	781B	Non - Absorbable	Silk - Black Braided	3/8 circle -19mm reverse cut	4/0	18"	740-3015
	C-6	784B	Non - Absorbable	Silk - Black Braided	3/8 circle -19mm reverse cut	3/0	18"	740-3025
	C26	778B	Non - Absorbable	Silk - Black Braided	1/2 circle -15mm reverse cut	3/0	18"	740-3000
	C-31	785B	Non - Absorbable	Silk - Black Braided	1/2 circle -24mm reverse cut	3/0	18"	740-3005

Pkg/12 **\$26.50**



27" Sutures

27" Plain Gut Sutures with Needles

Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C-6	554B	Absorbable	Plain Gut	3/8 circle -19mm reverse cut	3/0	27"	740-3060
	C-7	553B	Absorbable	Plain Gut	3/8 circle -19mm reverse cut	3/0	27"	740-3055

Pkg/12 **\$33.50**

27" Chromic Gut Sutures with Needles




Needle Style	Needle Code	Suture Code	Usage	Material	Needle Description	Suture Size	Suture Length	Item #
	C-6	559B	Absorbable	Chromic Gut	3/8 circle -19mm reverse cut	3/0	27"	740-3045
	C-6	560B	Absorbable	Chromic Gut	3/8 circle -19mm reverse cut	4/0	27"	740-3050

Pkg/12 **\$33.50**



Sutures

Advance Sutures are manufactured using precise methods, which provides consistent strength and durability to each needle. Super-sharp needles are siliconized for repeated passage through tissue. In addition, the smooth transition from needle to material maximizes the needle-attachment strength. Sterile packed in color-coded suture boxes with a front dispenser.

-  Silk = Light Blue
-  Chromic Gut = Beige
-  Plain Gut = Yellow



Advancing the Quality of Dental Care Products



Advance Dental Products from **MidwestDental**[®]
EQUIPMENT · SERVICES · SUPPLIES

(800) 766-2025 | MidwestDental.com

Not responsible for printing and typographical errors. Prices are subject to change without notice. Version 1_2023.