



# Safety Data Sheet

Issue Date: 12-Oct-2020

Revision Date: 12-Oct-2020

Version 1

## 1. IDENTIFICATION

### Product identifier

**Product Name** Advance Dental Bite Registration

### Other means of identification

**SDS #** ADV-006

**Product Code** 39000, 39001, 39002, 39003, 39004, 39005, 39006

### Recommended use of the chemical and restrictions on use

**Recommended Use** For industrial use.

### Details of the supplier of the safety data sheet

#### **Supplier Address**

Advance Dental Products  
P.O. Box 29279  
Atlanta, GA 30359

### Emergency telephone number

**Company Phone Number** (877) 896-7341 (Monday-Friday 8am - 4:45pm EST)  
**Emergency Telephone** INFOTRAC 1-352-323-3500 (International)  
1-800-535-5053 (North America)

## 2. HAZARDS IDENTIFICATION

**Appearance** viscous Light blue base  
White catalyst

**Physical state** Liquid

**Odor** Odorless unless flavored as specified

### Classification

This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200).

## 3. COMPOSITION/INFORMATION ON INGREDIENTS

The product contains no substances which, at their given concentration, are considered to be hazardous to health.

## 4. FIRST AID MEASURES

### Description of first aid measures

**Eye Contact** Flush immediately with plenty of water. If irritation persists, see eye doctor.

**Skin Contact** Wash off immediately with plenty of water for at least 15 minutes.

**Inhalation** Remove to fresh air.

**Ingestion** Small quantities will pass through system. Consult a physician.

### Most important symptoms and effects, both acute and delayed

**Symptoms** Slight irritant for eyes.

**Indication of any immediate medical attention and special treatment needed**

**Notes to Physician**                      Treat symptomatically.

**5. FIRE-FIGHTING MEASURES**

**Suitable Extinguishing Media**

Water spray (fog). Foam. Powder. Carbon dioxide (CO2).

**Unsuitable Extinguishing Media** None.

**Specific Hazards Arising from the Chemical**

Material can burn generating irritating gases.

**Protective equipment and precautions for firefighters**

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

**6. ACCIDENTAL RELEASE MEASURES**

**Personal precautions, protective equipment and emergency procedures**

**Personal Precautions**                      Use personal protective equipment as required.

**Environmental precautions**

**Environmental precautions**              See Section 12 for additional Ecological Information.

**Methods and material for containment and cleaning up**

**Methods for Containment**                  Prevent further leakage or spillage if safe to do so.

**Methods for Clean-Up**                      Collect in a suitable container for disposal. Wipe up affected area with absorbent material. Dispose of material in accordance with local, state and federal requirements.

**7. HANDLING AND STORAGE**

**Precautions for safe handling**

**Advice on Safe Handling**                      Does not require specific technical measures. Not that handling with latex or rubber gloves may inhibit set. Use only vinyl gloves.

**Conditions for safe storage, including any incompatibilities**

**Storage Conditions**                          Store in supplied packaging. Store away from heat sources. Temperature will affect set time.

**Incompatible Materials**                      Strong oxidizing agents.

**8. EXPOSURE CONTROLS/PERSONAL PROTECTION**

**Exposure Guidelines**

This product, as supplied, does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies

**Appropriate engineering controls**

**Engineering Controls**                          Apply technical measures to comply with the occupational exposure limits.

**Individual protection measures, such as personal protective equipment**

- Eye/Face Protection**                      None required under normal use.
- Skin and Body Protection**                None required under normal use.
- Respiratory Protection**                    Refer to 29 CFR 1910.134 for respiratory protection requirements.
- General Hygiene Considerations** Handle in accordance with good industrial hygiene and safety practice.

**9. PHYSICAL AND CHEMICAL PROPERTIES**

**Information on basic physical and chemical properties**

<b>Physical state</b>	Liquid		
<b>Appearance</b>	viscous Light blue base White catalyst	<b>Odor</b>	Odorless unless flavored as specified
<b>Color</b>	Light blue/white	<b>Odor Threshold</b>	Not determined
<b>Property</b>	<b>Values</b>	<b>Remarks • Method</b>	
pH	Not determined		
Melting point / freezing point	Not determined		
Boiling point / boiling range	>65.6 °C / >150 °F		
Flash point	>232.2 °C / >450 °F		
Evaporation Rate	Not determined		
Flammability (Solid, Gas)	Liquid-Not applicable		
Flammability Limit in Air			
Upper flammability or explosive limits	Not determined		
Lower flammability or explosive limits	Not determined		
Vapor Pressure	<5 mm Hg	@ 25°C (77°F)	
Vapor Density	Not determined		
Relative Density	>1		
Water Solubility	Insoluble in water		
Solubility in other solvents	Not determined		
Partition Coefficient	Not determined		
Autoignition temperature	Not determined		
Decomposition temperature	Not determined		
Kinematic viscosity	Not determined		
Dynamic Viscosity	Not determined		
Explosive Properties	Not determined		
Oxidizing Properties	Not determined		

**10. STABILITY AND REACTIVITY**

**Reactivity**  
Not reactive under normal conditions.

**Chemical stability**  
Stable under recommended storage conditions.

**Possibility of hazardous reactions**  
None under normal processing.

**Conditions to Avoid**  
Keep out of reach of children.

**Incompatible materials**  
Strong oxidizing agents.

**Hazardous decomposition products**

None known based on information supplied.

**11. TOXICOLOGICAL INFORMATION****Information on likely routes of exposure****Product Information****Eye Contact** Avoid contact with eyes.**Skin Contact** Avoid contact with skin.**Inhalation** Do not inhale.**Ingestion** Do not ingest.**Component Information**

Not available

**Symptoms related to the physical, chemical and toxicological characteristics****Symptoms** Please see section 4 of this SDS for symptoms.**Delayed and immediate effects as well as chronic effects from short and long-term exposure****Carcinogenicity** Based on the information provided, this product does not contain any carcinogens or potential carcinogens as listed by OSHA, IARC or NTP.**Numerical measures of toxicity**

Not determined.

**12. ECOLOGICAL INFORMATION****Ecotoxicity**

The product is not classified as environmentally hazardous. However, this does not exclude the possibility that large or frequent spills can have a harmful or damaging effect on the environment.

**Component Information**

Not available

**Persistence/Degradability**

Not determined.

**Bioaccumulation**

There is no data for this product.

**Mobility**

Not determined

**Other Adverse Effects**

Not determined

**13. DISPOSAL CONSIDERATIONS****Waste Treatment Methods****Disposal of Wastes** Disposal should be in accordance with applicable regional, national and local laws and regulations.

**Contaminated Packaging**

Disposal should be in accordance with applicable regional, national and local laws and regulations.

**14. TRANSPORT INFORMATION****Note**

Please see current shipping paper for most up to date shipping information, including exemptions and special circumstances.

**DOT**

Not regulated

**IATA**

Not regulated

**IMDG**

Not regulated

**15. REGULATORY INFORMATION****International Inventories**

All ingredients are listed or exempt from listing on Chemical Substance Inventory

**US Federal Regulations****CERCLA**

This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355).

**SARA 313**

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372

**CWA (Clean Water Act)**

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

**US State Regulations****California Proposition 65**

This product does not contain any Proposition 65 chemicals.

**U.S. State Right-to-Know Regulations**

This product does not contain any substances regulated under applicable state right-to-know regulations

**16. OTHER INFORMATION**

<b><u>NFPA</u></b>	<b>Health Hazards</b> Not determined	<b>Flammability</b> Not determined	<b>Instability</b> Not determined	<b>Special Hazards</b> Not determined
<b><u>HMIS</u></b>	<b>Health Hazards</b> Not determined	<b>Flammability</b> Not determined	<b>Physical hazards</b> Not determined	<b>Personal Protection</b> Not determined

Issue Date: 12-Oct-2020  
 Revision Date: 12-Oct-2020  
 Revision Note: New format

**Disclaimer**

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

**End of Safety Data Sheet**