

SAFETY DATA SHEET

Section 1. Product And Company Identification

Product Name: SonicFill 3 Product Use: Dental product: Composite

Manufacturer: Kerr Corporation 1717 W. Collins Ave. Orange, CA 92867-5422 U.S.A.

Information Phone Number: 1-800-841-1428 (Customer Service)

Chemical Emergency Phone Number (Chemical Spills, Leaks, Fire, Exposure or Accident only): CHEMTREC 1-800-424-9300 (in the US) 1-703-527-3887 (Outside the US)

SDS Date of Preparation/Revision: August 14, 2019

Section 2. Hazards Identification

GHS Classification:

Skin Irritation Category 2 Eye Irritation Category 2A

Label Elements:



Hazard Phrases Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation.

Precautionary Phrases:

Wash hands thoroughly after handling. Wear eye protection, protective gloves. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. Take off contaminated clothing and wash it before reuse.



Section 3. Composition/Information on Ingredients

Component	CAS No.	Amount
Silicon dioxide	7631-86-9	30-60%
Glass, oxide, chemicals	65997-17-3	10-30%
Poly(oxy-1,2-ethanediyl), α,α'-[(1-	41637-38-1	10-30%
methylethylidene)di-4,1-phenylene]bis[ω-[(2-		
methyl-1-oxo-2-propen-1-yl)oxy]-		
Ytterbium trifluoride	13760-80-0	1-5%
2,2'-ethylenedioxydiethyl dimethacrylate	109-16-0	0.1-1%

Section 4. First Aid Measures

Inhalation: Immediately remove victim to fresh air. If breathing is difficult, oxygen should be administered by qualified personnel. If breathing has stopped, administer artificial respiration. Get immediate medical attention.

Skin Contact: Flush thoroughly with water. Get medical attention if irritation or symptoms of exposure develop. Remove and launder contaminated clothing before re-use.

Eye Contact: Rinse thoroughly with water. Get medical attention if irritation occurs and persists.

Ingestion: Do NOT induce vomiting. Rinse mouth with water. Never give anything by mouth to an unconscious or convulsing person. Keep the victim calm and warm. Get immediate medical attention.

Most important symptoms and effects, acute and delayed: Not expected to cause adverse effects.

Indication of immediate medical attention and special treatment, if needed: None required under normal conditions of use.

Section 5. Fire Fighting Measures

Suitable (and Unsuitable) Extinguishing Media: Use any media appropriate for the surrounding fire. Cool fire exposed containers with water.

Specific Hazards Arising from the Chemical: Combustion may produce carbon monoxide, carbon dioxide, halogenated compounds, and metal oxides.

Special Protective Equipment and Precautions for Fire-fighters: Firefighters should wear positive pressure self-contained breathing apparatus and full protective clothing for fires in areas where chemicals are used or stored. Cool fire-exposed containers with water. Contain water used in firefighting from entering sewers or natural waterways.

Section 6: Accidental Release Measures

Personal precautions, Protective equipment, and Emergency procedures: Evacuate spill area and keep unprotected personnel away. Avoid contact with eyes, skin and clothing. Wear appropriate protective clothing and equipment.



Environmental Precautions: Avoid releases to the environment. Report spill as required by local and federal regulations.

Methods and Materials for Containment and Cleaning up: Prompt cleanup and removal are necessary. Clean up any spills as soon as possible, using an absorbent material to collect it. Wash off with plenty of water.

Section 7. Handling and Storage

Precautions for Safe Handling: Prevent contact with eyes, skin and clothing. Always wear impervious gloves, chemical safety goggles and protective clothing when handling this material. Wash thoroughly with soap and water after handling. Do not eat, drink or smoke in the work area. Remove and wash contaminated clothing before reuse.

Empty containers retain product residues which can be hazardous. Follow all SDS precautions when handling empty containers.

Conditions for Safe Storage, Including any Incompatibilities: Store in a cool, dry, well-ventilated area away from direct sunlight. Keep in the original container or an approved alternative made from a compatible material, kept tightly closed when not in use. Containers that have been opened must be carefully resealed and kept upright to prevent leakage. Do not store in unlabeled containers.

Section 8. Exposure Controls / Personal Protection

Exposure Limits

Chemical	Exposure Limit
Silicon dioxide	6 mg/m ³ TWA NIOSH REL
Glass, oxide, chemicals	5 mg/m ³ TWA ACGIH TLV
Poly(oxy-1,2-ethanediyl), α,α'-[(1- methylethylidene)di-4,1-phenylene]bis[ω-[(2- methyl-1-oxo-2-propen-1-yl)oxy]-	None Established
Ytterbium trifluoride	2.5 mg/m ³ TWA ACGIH TLV
2,2'-ethylenedioxydiethyl dimethacrylate	None Established

Appropriate Engineering Controls: Use with adequate general or local exhaust ventilation to maintain exposure levels below the occupational exposure limits.

Respiratory Protection: None under normal use conditions with adequate ventilation. For operations where the occupational exposure limits are exceeded, an approved respirator with particulate cartridges is recommended. Equipment selection depends on contaminant type and concentration. Select in accordance with applicable regulations and good industrial hygiene practice. For firefighting, use self-contained breathing apparatus.

Hand protection: Impervious gloves are suggested to prevent skin contact. Contact your glove supplier for selection assistance.

Eye Protection: Chemical safety goggles are recommended if contact is possible.



Skin Protection: Wear protective clothing as needed to avoid skin contact and contamination of personal clothing.

Hygiene measures: Suitable eye and skin washing facilities should be available in the work area.

Section 9. Physical and Chemical Properties					
Appearance:	Various colored paste	Odor:	Fruity		
Odor Threshold:	Not available	pH:	Not available		
Melting/Freezing	Not available	Boiling	Not available		
Point:		Point/Range:			
Flash Point:	Not flammable	Evaporation	Not available		
		Rate:			
Flammability: (Solid,	Not applicable	Flammability	LEL: Not applicable		
Gas)		Limits:	UEL: Not applicable		
Vapor Pressure:	Not available	Vapor	Not available		
•		Density:			
Relative Density:	1.93	Solubilities:	Insoluble in water		
Partition Coefficient:	Not available	Autoignition	Not available		
(N-Octanol/Water)		Temperature:			
Decomposition	Not available	Viscosity:	Not available		
Temperature:		•			

Section 10. Stability and Reactivity

Reactivity: The product is not expected to be reactive.

Chemical Stability: Stable under normal storage and handling conditions.

Possibility of Hazardous Reactions: Hazardous polymerization will not occur.

Conditions to avoid: Avoid heat and direct sunlight.

Incompatible Materials: Oxidizing materials and acids.

Hazardous decomposition products: None if stored normally.

Section 11. Toxicological Information

Potential Health Effects:

Inhalation: May cause irritation of the nose, throat and upper respiratory tract. **Skin Contact:** Causes skin irritation.

Skin Contact: Causes skin imitation.

Eye Contact: Causes serious eye irritation.

Ingestion: Swallowing may cause gastrointestinal irritation and symptoms similar to those described under inhalation.

Chronic Hazards: None expected.

Skin Sensitization: May cause an allergic skin reaction.

Respiratory Sensitization: No data available. This product is not expected to cause respiratory sensitization.

Germ Cell Mutagenicity: None of the components have shown mutagenic activity in animal studies.



Carcinogen:

Silicon dioxide is listed as "The Agent is not Classifiable as to its Carcinogenicity to Humans" (Group 3) by IARC.

Glass, oxide, chemicals is listed as "The Agent is not Classifiable as to its Carcinogenicity to Humans" (Group 3) by IARC.

Ytterbium trifluoride is listed as "The Agent is not Classifiable as to its Carcinogenicity to Humans" (Group 3) by IARC.

None of the other components are listed as a carcinogen or potential carcinogen by IARC, NTP, ACGIH, or OSHA.

Developmental / Reproductive Toxicity: None of the components have been shown to cause reproductive or developmental toxicity.

Specific Target Organ Toxicity (Single Exposure): No data available.

Specific Target Organ Toxicity (Repeated Exposure): No data available.

Aspiration Toxicity: Not an aspiration hazard.

Acute Toxicity Values:

2,2'-ethylenedioxydiethyl dimethacrylate: LD50 Oral rat: 10837 mg/kg

Section 12. Ecological Information

Toxicity:

Silicon dioxide: 72 hr EC50 Pseudokirchneriella subcapitata 55.5 mg/L

Persistence and degradability: Biodegradation is not applicable to inorganic substances.

Bioaccumulative Potential:

Poly(oxy-1,2-ethanediyl), α,α' -[(1-methylethylidene)di-4,1-phenylene]bis[ω -[(2-methyl-1-oxo-2-propen-1-yl)oxy]- has a BCF of 2372, log P_{ow} 3.43 to 5.62, potential of bioaccumulative is high 2,2'-ethylenedioxydiethyl dimethacrylate: log P_{ow} 1.88, potential of bioaccumulative is low

Mobility in Soil: No data available.

Other Adverse Effects: No data available.

Section 13. Disposal Considerations

Disposal: For unused product, dispose of in accordance with Federal and local regulations. **Container Disposal:** Dispose of empty container in accordance with Federal and local regulations.

Section 14. Transport Information

	UN Number	UN Proper Shipping Name	Hazard Class(s)	Packing Group	Environmental Hazards
US DOT	None	Not Regulated			None
Canada TDG	None	Not Regulated			None



IMDG	None	Not Regulated		None
IATA/ICAO	None	Not Regulated		None

Special Precautions for User: None identified

Transport in Bulk According to Annex II MARPOL 73/78 and the IBC Code: Not applicable – product is transported only in packaged form.

Section 15. Regulatory Information

U.S. Federal Regulations:

EPA SARA 311/312 Hazard Classification: Refer to Section 2 for OSHA Hazard Classification.

EPA SARA 313: This Product Contains the Following Chemicals Subject to Annual Release Reporting Requirements Under SARA Title III, Section 313 (40 CFR 372): None

Protection Of Stratospheric Ozone: This product is not known to contain or to have been manufactured with ozone depleting substances as defined in 40 CFR Part 82, Appendix A to Subpart A.

CERCLA SECTION 103: This product is not subject to CERCLA reporting requirements; however, many states have more stringent release reporting requirements. Report spills required under federal, state and local regulations.

US EPA TSCA Inventory: All of the components of this product are listed on the Toxic Substances Control Act (TSCA) Chemical Substances Inventory or exempt.

Canadian Regulations:

Canadian Environmental Protection Act: All of the components in this product are listed on the Domestic Substances List (DSL) or exempt.

National Pollutant Release Inventory (NPRI): This Product Contains the Following Chemicals Subject to Annual Release Reporting Requirements NPRI: None

International Inventories

Australia: All of the components in this product are listed on the Australian Inventory of Chemical Substances (AICS) or exempt.

China: All of the components in this product are listed on the Inventory of Existing Chemical Substances in China (IECSC) or exempt.

European Union: All the components in this product are listed on the EINECS inventory or exempt. **Korea:** All of the components in this product are listed on the Korean Existing Chemicals List (KECL) or exempt.

New Zealand: All of the components in this product are listed on the New Zealand Inventory of Chemicals (NZIoC) or exempt.



Section 16. Other Information

Effective Date: August 14, 2019 Supersedes Date: New SDS Revision Summary: New SDS

The information and recommendations set forth herein are taken from sources believed to be accurate as of the date of preparation, however, KERR Corporation makes no warranty with respect to the accuracy or suitability of the recommendations, and assumes no liability to any use thereof.