

## SAFETY DATA SHEETS

**This SDS packet was issued with item:**

073781242

**The safety data sheets (SDS) in this packet apply to the individual products listed below. Please refer to invoice for specific item number(s).**

070430603 070430645 070431312 073779907 073779915 073779923 073779931 073779949 073779956 073779964  
073779972 073779980 073779998 073780004 073780012 073780020 073780152 073780160 073780178 073780186  
073780194 073780202 073780210 073780228 073780236 073780244 073780251 073780269 073780277 073780285  
073780293 073780533 073780541 073780558 073780566 073780574 073780582 073780590 073780608 073780616  
073780624 073780632 073780640 073780657 073780665 073780673 073780939 073780947 073780954 073780962  
073780970 073780988 073780996 073781002 073781010 073781028 073781036 073781044 073781051 073781069  
073781077 073781085 073781093 073781101 073781119 073781176 073781184 073781192 073781200 073781218  
073781226 073781234 073781259 073781267 073781275 073782117 073782125 073782133 073782141 073782158  
073782166 073782174 073782257 073782323 073782331 073782372 073782380 073782398 073782406 073782414  
073782422 073782430 073782448 073782455 073782463 073782471 073782489 073782497 073782505 073795002  
073795010 073795028 073795036 073795044 073795051 073795069 073795077 073795150 073795168 073795176  
073795184 073795192 073795200 073795218 073795226 073795234 073795242 073795259 073795267 073795275  
073795309 073795317 073795325 073796091 073796653 073796661 073796679 073796687 073796695 073796703  
073796711 073796729 073796737 073796745 073796752 073796760 073796778 073796786 073796794 073802501  
073802519 073802527 073802535 073802543 073802550 073802568 073802576 073802584 073802592 073802600  
073802618 073802626 073802634 073802642 073802659 073802667 073802675 073802683 073802691 073802709  
073802717 077659691 273008664 273008666 273009612 273009613 273009614 273009616 273009641 273009643  
273009644 273009705 273011083 273016507 273016508 273018739 273018963 273020702 273022869 273022871  
273022872

**1. Identification of the Substance/Preparation and of the Company/Undertakings****Product Details:**

Trade Name: Gutta Percha Points

**Manufacturer:****Coltène/Whaledent Inc.**

235 Ascot Parkway

Cuyahoga Falls, OH 44223/USA

Fax +1 330 916 7077

info@coltenewhaledent.com

Authorized EU Representative

**Coltène/Whaledent GmbH+ Co. KG**

Raiffeisenstrasse 30

89129 Langenau/Germany

Tel. +49 (0) 7345 805 0

Fax +49 (0) 7345 805 201

**Emergency Contact Numbers:**

USA &amp; Canada: +1 800 535 5053

Main Phone Number: +1 330 916 8800

**2. Composition/Information on Ingredients****Chemical Characterization:** Paper Points**Hazardous Components:**Chemical NameWeight %Index, R/S Phrases

None

**NON Hazardous Components**

USP Zinc Oxide (CAS# 1314-13-2)

&lt; 70%

Barium Sulphate (CAS# 7727-43-7)

&lt; 15%

FD&amp;C Yellow #5

&lt; 1%

**3. Hazard Identification**

Hazard Designation: Not hazardous.

Information Pertaining to Particular Dangers for Man and Environment: None.

Classification is in Line with BC Directive 1999:45/EC: N.A.

**4. First-Aid Measures**

After Inhalation: N.A.

After Skin Contact: N.A.

After Eye Contact: If scratches or punctures occur, seek medical attention immediately.

After Swallowing: Do not induce vomiting. Material is chemically inert.

**5. Fire-Fighting Measures****Flammability: Nonflammable****Flash Point: N.A.****LEL: N.A.****UEL: N.A.**Extinguishing Media: Water, CO<sub>2</sub>, dry chemical or foam.

Unsuitable Extinguishing Media: None known.

Protective Equipment: Wear NIOSH approved self-contained breathing respirator in positive pressure mode.

**6. Accidental Release Measures**

Personal Precautions: None.

Environmental Precautions: None.

Methods of Cleaning Up: N.A.

**7. Handling and Storage**

Precautions for Safe Handling: Always wear gloves and eye protection. Keep container closed when not in use.

Storage: Keep out of reach of children. Keep container closed.

**8. Exposure Controls/Personal Protection**

Components with Recommended Limits for the Work Place:

CAS No.Chemical NameWeight %Type/Value

None

Engineering Controls: N.A.

Work/Hygienic Protection: N.A.

Hand/Skin Protection: N.A.

Respiratory Protection: Wear safety mask when shortening or removing post from tooth.

Eye Protection: Wear safety goggles.

Date of Issue: December 2008

Coltène/Whaledent Inc.

N.A. = not applicable, N.K. = not known

Gutta Percha Points

**9. Physical and Chemical Properties**

Appearance: Colored Thermoplastic Cones  
Boiling Point: N.A.  
Flammability: Nonflammable  
Vapor Pressure: N.A.  
pH: N.A.

Odor: None  
Flash Point: N.A.  
Autoflammability: N.A.  
Explosive Properties: N.A.  
Viscosity: N.A.

Solubility : N.A.  
Method: N.A.  
Density: N.A.  
Explosive Limits in Air: N.K.  
Method: N.A.

**10. Stability and Reactivity**

Chemical Stability: Stable.  
Materials to Avoid: None.

Conditions to Avoid: N.K.  
Hazardous Decomposition Products: CO.

**11. Toxicological Information**

Acute Toxicity: N.A.  
Primary Effects to Skin: Not a primary dermal irritant.  
Primary Effects to Eyes: Not an ocular irritant.  
Primary Effects to Respiratory System: Not an acute inhalation hazard.  
Sensitization: N.A.  
Additional Toxicological Info: N.A.

**12. Ecological Information**

General Information: N.A.  
Classification of Water Endangerment: N.A.

**13. Disposal Considerations**

Product: Not considered a RCRA hazardous waste as identified in Fed. Reg. 40CFR261. If this product cannot be used as directed by the label, contact the local Environmental Control Agency or Agent.

**14. Transportation Information**

DOT Shipping Name:  
Not regulated  
Sea Transport IMDG/GGV Sea  
Class: N.A.  
Air Transport ICAO-Ti and IATA-DGR  
Class: N.A.

DOT Hazard Class:  
Not regulated  
UN. No.: N.A.  
UN. No.: N.A.

UN-No.: N.A.  
EmS-No.: N.A.  
MFGA-No.: N.A.

**15. Regulatory Information**

TSCA Inventory Status: This product and/or all of its components are included on the TSCA Inventory of Chemical Substances.  
SARA 311/312 Hazard Categories: None.  
SARA 302 Extremely Hazardous Substances: None.  
SARA 313 Toxic Chemicals:

Chemical Name  
None

CAS No.

Weight %

Classification According to EEC Guidelines: N.A.  
Classification According to VBF: N.A.

Canadian Regulations require the inclusion of the following statements under "Regulatory Information";  
"This product has been classified in accordance with the hazard criteria of the controlled Products Regulations and the MSDS contains all the information required by the Controlled Products Regulations."

**16. Other Information**

Revision Summary: 6/14/04 Changed Address and added EU Rep. 12/5/08 Updated to current EU standards in header and section 3.

The information and recommendations set forth herein are presented in good faith and are believed to be correct as of the date hereof. Coltène/Whaledent Inc., however, makes no representations as to the completeness of security of this information and supplies it on the condition that the persons receiving same will make their own determination as to its suitability of their purposes prior to use. In no event will Coltène/Whaledent Inc. be responsible for damages of any nature whatsoever resulting from use of or reliance upon information.