

SAFETY DATA SHEETS

This SDS packet was issued with item:

072347896

The safety data sheets (SDS) in this packet apply to the individual products listed below. Please refer to invoice for specific item number(s).

072347888

CROSSTEX SPARKLE V – VARNISH

Item code: UFVBG

MATERIAL SAFETY DATA SHEET

PRODUCT IDENTIFICATION

Flavor: Bubble Gum

COMPOSITION/INFORMATION ON INGREDIENTS

Product Contents: Colophony Desensitizing Dental Varnish, Tooth Shade A2-Bulk

Hazardous Ingredients:

Rosin 50-70% EINECS: 232-475-7 CAS: 8050-09-7 [Sens.] R43

Ethanol 10-30% EINECS: 200-578-6 CAS: 64-17-5 [F] R11

Sodium Fluoride 1-10% EINECS: 231-667-8 CAS: 7681-49-4 [T] R25; [-] R32; [Xi] R36/38

HAZARDS IDENTIFICATION

Harmful if swallowed. Contact with acids liberates very toxic gas. May cause sensitisation by skin contact.

FIRST AID MEASURES

Symptoms Action

Skin Contact: There may be mild irritation at the site of contact. Wash immediately with plenty of soap and water.

Eye Contact: There may be irritation and redness. Bathe the eye with running water for 15 minutes.

Ingestion: There may be irritation of the throat. Wash out mouth with water.

Inhalation: No symptoms. Consult a doctor.

FIRE-FIGHTING MEASURES

Extinguishing media: Suitable extinguishing media for the surrounding fire should be used.

Use water spray to cool containers.

Exposure hazards: In combustion emits toxic fumes.

Protection of firefighters: Wear self-contained breathing apparatus.

Wear protective clothing to prevent contact with skin and eyes.

ACCIDENTAL RELEASE MEASURES

Personal precautions: Refer to section 8 of SDS for personal protection details.

Turn leaking containers leak-side up to prevent the escape of liquid.

Environmental precautions: Do not discharge into drains or rivers. Contain the spillage using bunding.

Clean-up procedures: Absorb into dry earth or sand. Transfer to a closable, labelled salvage container for disposal by an appropriate method.

HANDLING AND STORAGE

Storage conditions: Store in cool, well ventilated area. Keep container tightly closed.

EXPOSURE CONTROLS/PERSONAL PROTECTION

Hazardous ingredients: ETHANOL, TWA 1920 mg/m³

Respiratory protection: Respiratory protection not required.

Hand protection: Protective gloves.

Eye protection: Safety glasses. Ensure eye bath is at hand.

Skin protection: Protective clothing.

PHYSICAL AND CHEMICAL PROPERTIES

State: Paste

Color: Off-white

Odor: Characteristic odor

Evaporation rate: Moderate

Solubility in water: Slightly soluble

Also soluble in: Ethanol

Viscosity: Viscous

STABILITY AND REACTIVITY

Stability: Stable under normal conditions.

Conditions to avoid: Heat

Materials to avoid: Strong oxidizing agents. Strong acids.

Haz. Decomp. Products: In combustion emits toxic fumes.

PRD/UFVBG-INS

Rev. 0 - 7/10

TOXICOLOGICAL INFORMATION

Hazardous ingredients:

ETHANOL SODIUM FLUORIDE

IVN RAT LD50 1440 mg/kg ORL MUS LD50 57 mg/kg

ORL MUS LD50 3450 mg/kg ORL RAT LD50 52 mg/kg

ORL RAT LD50 7060 mg/kg SCU RAT LD50 175 mg/kg

Routes of exposure: Refer to section 4 of SDS for routes of exposure and corresponding symptoms.

ECOLOGICAL INFORMATION

Mobility: Readily absorbed into soil.

Persistence and degradability: Biodegradable.

Bioaccumulative potential: No bioaccumulative potential.

Other adverse affects: Negligible ecotoxicity.

DISPOSAL CONSIDERATIONS

NB: The user's attention is drawn to the possible existence of regional or national regulations regarding disposal.

TRANSPORT INFORMATION

ADR/RID UN no: -

Shipping name: NOT CLASSIFIED AS DANGEROUS IN THE MEANING OF TRANSPORT REGULATIONS.

IMDG/IMA UN no: -

IATA/ICAO UN no: -

REGULATORY INFORMATION

Hazard symbols: Harmful.

Risk phrases: R22: Harmful if swallowed;

R32: Contact with acids liberates very toxic gas.

R43: May cause sensitisation by skin contact.

Safety phrases: S24: Avoid contact with skin.

S37: Wear suitable gloves.

PRECAUTIONARY PHRASES: Restricted to professional users.

Note: The regulatory information given above only indicates the principal regulations specifically applicable to the product described in the safety data sheet.

The user's attention is drawn to the possible existence of additional provisions which complete these regulations. Refer to all applicable national, international and local regulations or provisions.

OTHER INFORMATION

Risk phrases used in s.2: R43: May cause sensitisation by skin contact.

R11: Highly flammable.

R25: Toxic if swallowed.

R32: Contact with acids liberates very toxic gas.

R36/38: Irritating to eyes and skin.

Legal disclaimer:

CAUTION: PRODUCT FOR PROFESSIONAL USE ONLY.

The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility.

SPECIAL INFORMATION – 24 HR EMERGENCY TELEPHONE NUMBER:

*INFOTRAC = 1-800-535-5035 or local poison control center.

Note: Crosstex Int. has compiled this data based on the research & input of qualified experts in the field of microbiology, chemistry and quality assurance.

The data are not to be taken as a warranty or representation for which Crosstex Int. assumes legal responsibility.

They are offered only for your consideration and verification.