

## **SAFETY DATA SHEETS**

**This SDS packet was issued with item:**

077044365

**The safety data sheets (SDS) in this packet apply to the individual products listed below. Please refer to invoice for specific item number(s).**

077044399 077049190 077049216

# Dux Dental Material Safety Data Sheet

Dux Dental  
600 E. Hueneme Road  
Oxnard, CA, USA

Phone: 800.833.8267  
Fax: 805.488.2266  
www.duxdental.com

**Product** Aluminum Chloride Hemostatic Gingival Retraction Solutions  
**Other Name** Hemogin L (25 % Aluminum Chloride), Styptin (20 %), Hemodettes (20 %)  
**REF** 13052, 13055, 13057, 13058, 13150, 13156, 13157

**MSDS Number** MS110

## Hazardous Ingredients

Principal Hazardous Components	pel / tlv	CAS Number	Percentage
None	NA	NA	NA

## Physical and Chemical Characteristics

**Appearance** Clear, slightly viscous liquids, green or clear, with a bitterseet taste.  
**Boiling Point** > 212°F  
**Melting Point** NA  
**Vapor Pressure** No Data Found  
**Evaporation Rate** Slower than water.  
**Odor** Spearmint, None  
**Vapor Density** NA  
**Specific Gravity** 1.25  
**Water Solubility** Aqueous System

## Fire & Explosion Data

**Explosion** Not considered to be an explosion hazard.  
**Fire** Not considered to be a fire hazard.  
**Fire Ext. Type** A B C  
**Flashpoint** NA  
**Special Fire Procedures** None  
**Unusual Hazards** None Known

## Reactivity Data

**Stability** Stable  
**Incompatibilities** None Known  
**Hzd Decomp Products** Fumes of hydrochloric acid.  
**Hazardous Polymerization** Will not occur.

# Dux Dental Material Safety Data Sheet

Dux Dental  
600 E. Hueneme Road  
Oxnard, CA, USA

Phone: 800.833.8267  
Fax: 805.488.2266  
www.duxdental.com

**Product** Aluminum Chloride Hemostatic Gingival Retraction Solutions  
**Other Name** Hemogin L (25 % Aluminum Chloride), Styptin (20 %), Hemodettes (20 %)  
**REF** 13052, 13055, 13057, 13058, 13150, 13156, 13157

**MSDS Number** MS110

## Health Hazard Information

<b>Ingestion</b>	May cause irritation. Dilute with water if ingested. Consult physician.
<b>Inhalation</b>	NA
<b>Skin/Eye Contact</b>	May cause irritation. Rinse with water for 15 min, consult physician.
<b>Toxicity Data</b>	May cause reversible or irreversible changes to exposed tissue, not permanent injury or death; can cause considerable discomfort. Use in accordance with manufacturers instructions and sound professional
<b>Acute Hazards</b>	None
<b>Chronic Hazards</b>	None
<b>Carcinogenicity</b>	NA, no product ingredient listed by NTP, IARC, or OSHA.

## Precautions for Safe Handling and Use

<b>Cleanup</b>	Collect while minimizing production of dust.
<b>Disposal</b>	Dispose of in accordance with local, state, and federal regulations.
<b>Precautions</b>	As part of good industrial, personal hygiene, and safety practice avoid all unnecessary exposure to the substance and ensure prompt removal from skin, eyes, and clothing.

## Occupational Control Measures

<b>Eye Protection</b>	NA
<b>Gloves</b>	NA
<b>Dust/Particulate Mask</b>	NA
<b>Other</b>	NA

## Additional Information

<b>DOT Hazard Class</b>	Not Regulated
<b>DOT Small Qty Exception</b>	NA
<b>DOT Other Exception</b>	NA
<b>DOT Exception Permitted</b>	NA
<b>Allergenic Ingredients</b>	This product does not contain Latex, Peanut Products, Gluten, or Red Dye #3.

The information contained herein is provided in good faith and is believed to be correct as of the date hereof. However, the company makes no representation as to the comprehensiveness or accuracy of the information. It is expected that individuals receiving the information will exercise their independent judgement in determining its appropriateness for a particular purpose. Accordingly, we will not be responsible for damages of any kind resulting from the use of or reliance upon such information. No representations or warranties, either express or implied, of merchantability, fitness for a particular purpose or of any other nature are made hereunder with respect to the information set forth herein or to the product to which the information refers.