

SAFETY DATA SHEETS

This SDS packet was issued with item:

074353900

N/A



MATERIAL SAFETY DATA SHEET

Dist. By: SybronEndo A Division of Sybron Dental Specialties 1332 S. Lone Hill Ave Glendora, CA 91740	Emergency Phone: 800-346-3636 Information Phone: 800-346-3636	NFPA HAZARD CLASSIFICATION HEALTH: 1 FLAMIBILITY: 0 REACTIVITY: 0 SPECIFIC HAZARD: OXY N/E Not established N/A Not applicable	
Section I. Product Identification			
Label Name: SlickGel ES®			
Part Number: 973-1015			
Revision Date: April 1, 2009			
Section II. Hazardous Ingredients/ Identity Information			
Hazardous Components	CAS#	OSHA PEL	ACGIH TLV
Ethylenediamine	62-33-9	N/A	N/A
Tetraacetic acid			
Urea Peroxide	124-43-06	N/A	N/A
Propylene Glycol	57-55-6	N/A	N/A
Polyethylene Glycol	25322-68-3	15 (mg/m ³)	10 (mg/m ³)
Section III. Physical/Chemical Identity			
Boiling Point	N/E	Specific Gravity	N/E
Vapor Pressure	N/E	Melting Point	N/E
Vapor Density	N/E	Evaporation rate	N/E
Solubility in Water	Soluble		
Appearance and odor	Off White Gel, slight odor		
Section IV. Fire and Explosion Hazard Data			
Flash Point (method used)	N/A	Unusual fire and explosive hazards	None
Flammability	N/E		
LEL: N/A	UEL: N/A		
Special fire fighting procedures	Wear protective self-breathing apparatus		
Extinguishing Media	Water fog, CO ₂ or dry chemical		
Section V. Reactivity Data			
Stability	Stable: yes Unstable:	Hazardous polymerization	May occur: Will not occur: yes
Conditions to Avoid	High Heat	Conditions to avoid	N/A
Incompatibility (materials to avoid)	Alkalies,metals (high carbon steels), ammonia, combustible material, salts, sulfates, phenols		
Hazardous decomposition or byproducts	If burning may liberate Nitrous oxides, CO ₂ , CO, ammonia		

Section VI. Health Hazard Information			
Route or Entry:	Inhalation: No Ingestion: Yes Skin: Yes	Health Hazards (acute or chronic)	May cause mild skin irritation. Avoid Eye contact, can cause irritation
Carcinogenicity:	NTP: no IRAC monographs: no OSHA regulated: no	Signs and symptoms of exposure	Mild irritation
Medical conditions generally aggravated by exposure:	None known	Emergency First Aid Procedures	
		Eye Contact: Flush eye with running water for 15 minutes Skin Contact: Wash with soap and water Swallowed: Take two glasses of water and call physician	
Section VII. Precautions for Safe Handling and Use			
Steps to be taken in case material is spilled or released:	Scoop up, mop up and wash with water and dispose in a DOT approved container		
Waste disposal methods:	Extra contents can be flushed down sanitary sewer. Follow local, state and federal regulations to dispose packaging		
Precautions taken in storage and handling:	Store in cool and dry place. Wash hands with soap and water after using		
Other precautions:	None		
Section VIII. Control Measures			
Respiratory protection:	None under normal anticipated circumstances		
Ventilation:	General: no Mechanical: no	Local: no Special: no	
Protective gloves:	Recommended		
Eye Protection:	Safety glasses		
Other Protective clothing or equipment:	Normal protective garment/apron		
Work/hygienic practices:	Wash hands with soap and water after use		
Note: The information provided herein was collected and reported in good faith by Jordco Inc. While the information is believed to be accurate, Jordco Inc. makes no representation as to its accuracy or sufficiency. Data is provided without any warranty, express or implied. It is the user's responsibility to determine whether the product is suitable for their particular purposes and they assume all risks from the use, safe handling and disposal of the product. User assumes all liability for loss, injury damage or expense resulting from improper use of this product.			