

## Material Safety Data Sheet acc. to OSHA HCS

Printing date 01/10/2014

Reviewed on 01/10/2014

### 1 Identification

· **Product identifier**

· **Trade name:** **Venus Temp 2**

· **Relevant identified uses of the substance or mixture and uses advised against**  
No further relevant information available.

· **Application of the substance / the preparation** Auxiliary for manufacture of dental prothesis

· **Details of the supplier of the safety data sheet**

· **Manufacturer/Supplier:**

Heraeus Kulzer GmbH  
Grüner Weg 11, D-63450 Hanau

Tel.: 0800 4372522

· **Information department:**

Marc Henn  
Tel. +1 (574) 2995444 / Fax: +1 (574) 2912542  
e-mail: marc.henn@kulzer-dental.com

· **Emergency telephone number:**

Emergency CONTACT (24-Hour-Number)  
GBK/Infotrac ID 105860: (domestic) 1 800 535 5053 or international (001) 352 323 3500

### 2 Hazard(s) identification

· **Classification of the substance or mixture**

· **Classification according to Directive 67/548/EEC or Directive 1999/45/EC not applicable**

· **Information concerning particular hazards for human and environment:**

The product does not have to be labelled due to the calculation procedure of the "General Classification guideline for preparations of the EU" in the latest valid version.

· **Label elements**

· **Labelling according to EU guidelines:**

The product is not subject to classification according to the calculation methods of the "General Classification Guideline for Preparations of the EU" as issued in the latest valid version.

The product is not subject to identification regulations under EU Directives and the Ordinance on Hazardous Materials (GefStoffV).

Observe the general safety regulations when handling chemicals

· **Classification system**

· **NFPA ratings for USA (scale 0-4)**



· **HMIS-Ratings (Scale 0-4)**



· **Other hazards -**

### 3 Composition/information on ingredients

· **Chemical characterization: Mixtures**

· **Description:** Composition based on methacrylates

· **Dangerous components:** Void

USA

(Contd. on page 2)

## Material Safety Data Sheet acc. to OSHA HCS

Printing date 01/10/2014

Reviewed on 01/10/2014

Trade name: **Venus Temp 2**

(Contd. of page 1)

### 4 First-aid measures

- **Description of first aid measures**
  - **General information** No special measures required.
  - **After inhalation** Supply fresh air; consult doctor in case of complaints.
  - **After skin contact**  
Immediately wash with water and soap and rinse thoroughly.  
If skin irritation continues, consult a doctor.
  - **After eye contact**  
Rinse opened eye for several minutes under running water. Then consult a doctor.
  - **After swallowing** Immediately call a doctor.
  - **Information for doctor**
    - **Most important symptoms and effects, both acute and delayed**  
No further relevant information available.
    - **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### 5 Fire-fighting measures

- **Extinguishing media**
  - **Suitable extinguishing agents**  
CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.  
Use fire fighting measures that suit the environment.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
  - **Protective equipment:** No special measures required.
- **Additional information** -

### 6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures**  
Wear protective equipment. Keep unprotected persons away.
- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**  
Dispose of the collected material according to regulations.  
Absorb with liquid binding material (diatomite, universal binders, for small amounts tissues).
- **Reference to other sections**  
See Section 7 for information on safe handling  
See Section 13 for disposal information.

### 7 Handling and storage

- **Handling**
  - **Precautions for safe handling**  
Wear protective equipment. Keep unprotected persons away.  
Please observe the additional instructions in the product's instructions for use.  
Store in cool, dry place in tightly closed receptacles.
  - **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
  - **Storage**
    - **Requirements to be met by storerooms and receptacles:**  
Protect from exposure to the light.
    - **Information about storage in one common storage facility:** Store away from foodstuffs.

(Contd. on page 3)

## Material Safety Data Sheet acc. to OSHA HCS

Printing date 01/10/2014

Reviewed on 01/10/2014

**Trade name: Venus Temp 2**

(Contd. of page 2)

- **Further information about storage conditions:** Store in dry conditions.
- **Recommended storage temperature:** <25 °C
- **Specific end use(s)** No further relevant information available.

### 8 Exposure controls/personal protection

- **Control parameters**
  - **Components with limit values that require monitoring at the workplace:**  
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
    - **Additional information:** The lists that were valid during the creation were used as basis.
- **Exposure controls**
  - **Personal protective equipment**
    - **General protective and hygienic measures**  
The usual precautionary measures for handling chemicals should be followed.
    - **Breathing equipment:** Not required.
    - **Protection of hands:** Impervious gloves
      - **Penetration time of glove material**  
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
      - **For the permanent contact of a maximum of 15 minutes gloves made of the following materials are suitable:**  
Nitrile rubber, NBR
    - **Eye protection:** Protective goggles are recommended.
    - **Body protection:** Light weight protective clothing

### 9 Physical and chemical properties

- **Information on basic physical and chemical properties**
  - **General Information**
    - **Appearance:**
      - **Form:** Pasty
      - **Color:** Various colors
    - **Odor:** Nearly odorless
    - **Odour threshold:** Not determined.
  - **pH-value:** Neutral
  - **Change in condition**
    - **Melting point/Melting range:** undetermined
    - **Boiling point/Boiling range:** undetermined
  - **Flash point:** Not applicable
  - **Flammability (solid, gaseous)** Not applicable.
  - **Ignition temperature:**
    - **Decomposition temperature:** Not determined.
  - **Auto igniting:** Product is not selfigniting.
  - **Danger of explosion:** Product does not present an explosion hazard.
  - **Explosion limits:**
    - **Lower:** Not determined.
    - **Upper:** Not determined.

(Contd. on page 4)

## Material Safety Data Sheet acc. to OSHA HCS

Printing date 01/10/2014

Reviewed on 01/10/2014

**Trade name: Venus Temp 2**

(Contd. of page 3)

· <b>Vapor pressure:</b>	Not determined.
· <b>Density at 20 °C (68 °F):</b>	1.52 g/cm <sup>3</sup> (12.684 lbs/gal)
· <b>Relative density</b>	Not determined.
· <b>Vapour density</b>	Not determined.
· <b>Evaporation rate</b>	Not determined.
· <b>Solubility in / Miscibility with</b>	
· <b>Water:</b>	Not miscible or difficult to mix
· <b>Partition coefficient (n-octanol/water):</b>	Not determined.
· <b>Viscosity:</b>	
· <b>dynamic:</b>	Not determined.
· <b>kinematic:</b>	Not determined.
· <b>Other information</b>	No further relevant information available.

### 10 Stability and reactivity

- **Possibility of hazardous reactions** No dangerous reactions known
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** none
- **Additional information:** -

### 11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
  - **Primary irritant effect:**
    - **on the skin:** No irritant effect.
    - **on the eye:** No irritating effect.
  - **Sensitization:** No sensitizing effects known.
- **Additional toxicological information:**  
When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.
- **Carcinogenic categories**

#### · **IARC (International Agency for Research on Cancer)**

None of the ingredients is listed.

#### · **NTP (National Toxicology Program)**

None of the ingredients is listed.

### 12 Ecological information

- **Toxicity**
  - **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behavior in environmental systems:**
  - **Bioaccumulative potential** No further relevant information available.
  - **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
  - **General notes:** Avoid transfer into the environment.

(Contd. on page 5)

USA

## Material Safety Data Sheet acc. to OSHA HCS

Printing date 01/10/2014

Reviewed on 01/10/2014

**Trade name: Venus Temp 2**

(Contd. of page 4)

· **Other adverse effects** No further relevant information available.

**13 Disposal considerations**

- **Waste treatment methods**
  - **Recommendation** Smaller quantities can be disposed of with household waste.
- **Uncleaned packagings:**
  - **Recommendation:** Disposal must be made according to official regulations.

\* **14 Transport information**

· <b>UN-Number</b>	· DOT, ADR, ADN, IMDG, IATA	Void
· <b>UN proper shipping name</b>	· DOT, ADN, IMDG, IATA	Void
	· ADR	Void
· <b>Transport hazard class(es)</b>	· DOT, ADR, ADN, IMDG, IATA	
	· Class	Void
· <b>Packing group</b>	· DOT, ADR, IMDG, IATA	Void
· <b>Environmental hazards:</b>	· <b>Marine pollutant:</b>	No
· <b>Special precautions for user</b>	Not applicable.	
· <b>Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code</b>	Not applicable.	
· <b>Transport/Additional information:</b>	-	

\* **15 Regulatory information**

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**

· <b>SARA Section 355 (extremely hazardous substances)</b>
None of the ingredients is listed.

  - **TSCA (Toxic Substances Control Act)** Medical devices are exempted from TSCA.
- **Carcinogenicity categories**

· <b>TLV (Threshold Limit Value established by ACGIH)</b>
None of the ingredients is listed.
- **OSHA-Ca (Occupational Safety & Health Administration)**

None of the ingredients is listed.
------------------------------------

**16 Other information**

These data are based on our present knowledge. However, they shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- \* **Data compared to the previous version altered.**

USA