*Printing date 09/19/2013* 

Version number 5

Reviewed on 09/19/2013

# 1 Identification

· Product identifier

· Trade name: Tetric EvoCeram

• *Relevant identified uses of the substance or mixture and uses advised against* No further relevant information available.

· Application of the substance / the preparation Light-curing resin-based dental restorative material

· Details of the supplier of the safety data sheet

• Manufacturer/Supplier: Ivoclar Vivadent Inc. 175 Pineview Drive, Amherst, N.Y. 14228 USA Tel. +1 800 533 6825 Fax +1 716 691 2285

Ivoclar Vivadent Inc. 1-6600 Dixie Road Mississauga, Ontario L5T 2Y2 Canada Phone: +1 905 670 8499 Fax: +1 905 670 3102

 Information department: Quality Assurance / Regulatory Affairs
Emergency telephone number: 24 Hour Emergency Assistance: Emergency-Call USA - Infotrac: 1-800-535-5053 Emergency-Call Canada - Canutec: 1-613-996-6666

General MSDS Assistance: US: 1-800-533-6825 Canada: 1-800-263-8182

## 2 Hazard(s) identification

- · Classification of the substance or mixture
- · Classification according to Directive 67/548/EEC or Directive 1999/45/EC
- Irritant
- May cause sensitisation by skin contact.
- Classification system:

The classification was made according to the latest editions of international substances lists, and expanded upon from company and literature data.

#### · Label elements

#### · Labelling according to EU guidelines:

Irritant

Medical devices as defined in Directive 93/42/EEC and which are invasive or used in direct physical contact with the human body, are exempted from the provisions of Dangerous Substance Directive 1999/45/ EC usually if they are in the finished state and intended for the final user.

#### · Code letter and hazard designation of product:

(Contd. on page 2)

Printing date 09/19/2013

Version number 5

Reviewed on 09/19/2013

## Trade name: Tetric EvoCeram

(Contd. of page 1)

• *Hazard-determining components of labeling: urethane dimethacrylate* 

## · Risk phrases:

May cause sensitisation by skin contact.

#### · Safety phrases:

Avoid contact with skin. Wear suitable gloves. This material and its container must be disposed of as hazardous waste.

• Classification system: • NFPA ratings (scale 0 - 4)



#### · HMIS-ratings (scale 0 - 4)

HEALTH①Health = 0FIRE①Fire = 0REACTIVITY0Reactivity = 0

- · Other hazards
- · Results of PBT and vPvB assessment
- **PBT:** Not applicable.
- · vPvB: Not applicable.

### 3 Composition/information on ingredients

· Chemical characterization: Mixtures

• Description: Mixture of the substances listed below with nonhazardous additions.

### · Dangerous components:

72869-86-4 urethane dimethacrylate

1565-94-2 Bis-GMA

### 4 First-aid measures

- · Description of first aid measures
- · After inhalation:
- Supply fresh air and to be sure call for a doctor.
- In case of unconsciousness place patient stably in side position for transportation.
- · After skin contact: Rinse with water.
- · After eye contact:
- Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- · After swallowing: If symptoms persist consult doctor.
- · Information for doctor:
- Most important symptoms and effects, both acute and delayed No further relevant information available.
- · Indication of any immediate medical attention and special treatment needed
- No further relevant information available.

(Contd. on page 3)

USA

2.5-10%

2.5-10%

#### Version number 5

Reviewed on 09/19/2013

## Trade name: Tetric EvoCeram

(Contd. of page 2)

## **5** *Fire-fighting measures*

· Extinguishing media

• Suitable extinguishing agents:

CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

• Special hazards arising from the substance or mixture No further relevant information available.

• Advice for firefighters

· Protective equipment: No special measures required.

### 6 Accidental release measures

- · Personal precautions, protective equipment and emergency procedures Not required.
- · Environmental precautions: Do not allow to enter sewers/ surface or ground water.
- Methods and material for containment and cleaning up: Pick up mechanically.
- · Reference to other sections

See Section 7 for information on safe handling. See Section 8 for information on personal protection equipment. See Section 13 for disposal information.

### 7 Handling and storage

#### · Handling:

• *Precautions for safe handling* Only adequately trained personnel should handle this product. For use in dentistry only.

• Information about protection against explosions and fires: No special measures required.

· Conditions for safe storage, including any incompatibilities

· Storage:

- Requirements to be met by storerooms and receptacles: Store only in the original receptacle.
- · Information about storage in one common storage facility: Not required.
- Further information about storage conditions:
- Protect from exposure to the light.

Protect from heat and direct sunlight.

Keep receptacle tightly sealed.

• Specific end use(s) No further relevant information available.

## 8 Exposure controls/personal protection

• Additional information about design of technical systems: No further data; see item 7.

- · Control parameters
- · Components with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

- Additional information: The lists that were valid during the creation were used as basis.
- · Exposure controls
- · Personal protective equipment:
- General protective and hygienic measures:

Usual hygienic measures for dental practice. Wash hands before breaks and at the end of work.

Avoid contact with the eyes and skin.

(Contd. on page 4)

(Contd. of page 3)

# Safety Data Sheet acc. to OSHA HCS

*Printing date 09/19/2013* 

Version number 5

Reviewed on 09/19/2013

## Trade name: Tetric EvoCeram

Keep away from foodstuffs, beverages and feed. • Breathing equipment: Not required. • Protection of hands:

Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

· Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· Eye protection: Safety glasses

Information on basic physical and General Information Appearance:	chemical properties	
Form:	Pasty	
Color:	According to product specification	
Odor:	Characteristic	
Odour threshold:	Not determined.	
pH-value:	Not applicable.	
Change in condition Melting point/Melting range: Boiling point/Boiling range:	Undetermined. Undetermined.	
Flash point:	Not applicable.	
Flammability (solid, gaseous):	Not determined.	
Auto igniting:	Product is not selfigniting.	
Danger of explosion:	Product does not present an explosion hazard.	
Explosion limits:		
Lower:	Not determined.	
Upper:	Not determined.	
Vapor pressure:	Not applicable.	
Density at 20 °C (68 °F):	2.1 g/cm <sup>3</sup> (17.525 lbs/gal)	
Relative density	Not determined.	
Vapour density	Not applicable.	
Evaporation rate	Not applicable.	

Version number 5

Reviewed on 09/19/2013

## Trade name: Tetric EvoCeram

		(Contd. of page
Solubility in / Miscibility with	h	
Water:	Nearly insoluble.	
Partition coefficient (n-octan	ol/water): Not determined.	
Viscosity:		
Dynamic:	Not applicable.	
Kinematic:	Not applicable.	
Other information	No further relevant information available.	

# 10 Stability and reactivity

· Reactivity No further relevant information available.

#### · Chemical stability

Stable under normal handling and storage conditions.

Photoreactive (Danger of polymerization).

- Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid
- Keep away from heat and direct sunlight.
- Protect from exposure to the light.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: None under normal conditions of storage and use.

## **11** Toxicological information

· Information on toxicological effects

#### • Acute toxicity:

- Primary irritant effect:
- on the skin: No irritant effect.
- on the eye: Irritating effect.
- · Sensitization: Sensitization possible through skin contact.
- Additional toxicological information:
- The product shows the following dangers according to internally approved calculation methods for preparations:

Irritant

- · Carcinogenic categories
- · NTP (National Toxicology Program)

None of the ingredients is listed.

## **12** Ecological information

- · Toxicity
- Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- · Behavior in environmental systems:
- · Bioaccumulative potential No further relevant information available.
- *Mobility in soil* No further relevant information available.

(Contd. on page 6)

Printing date 09/19/2013

Version number 5

Reviewed on 09/19/2013

# Trade name: Tetric EvoCeram

(Contd. of page 5)

#### · Additional ecological information:

· General notes:

Water hazard class 1 (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

- · Results of PBT and vPvB assessment
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- Other adverse effects No further relevant information available.

### **13 Disposal considerations**

#### · Waste treatment methods

#### · Recommendation:

Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.

#### · Uncleaned packagings:

• *Recommendation: Disposal must be made according to official regulations.* 

# 14 Transport information

· UN-Number · DOT, ADR, ADN, IMDG, IATA	Void
· UN proper shipping name · DOT, ADR, ADN, IMDG, IATA	Void
· Transport hazard class(es)	
· DOT, ADR, ADN, IMDG, IATA · Class	Void
· Packing group · DOT, ADR, IMDG, IATA	Void
• Environmental hazards: • Marine pollutant:	No
· Special precautions for user	Not applicable.
• Transport in bulk according to Annex II o MARPOL73/78 and the IBC Code	f Not applicable.

## **15 Regulatory information**

· Safety, health and environmental regulations/legislation specific for the substance or mixture

· Sara

· Section 355 (extremely hazardous substances):

None of the ingredients is listed.

· Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

(Contd. on page 7)

(Contd. of page 6)

A4

# Safety Data Sheet acc. to OSHA HCS

#### Version number 5

Reviewed on 09/19/2013

## Trade name: Tetric EvoCeram

• TSCA (Toxic Substances Control Act):

All ingredients are listed.

#### · Proposition 65

· Chemicals known to cause cancer:

None of the ingredients is listed.

• Chemicals known to cause reproductive toxicity for females: None of the ingredients is listed.

· Chemicals known to cause reproductive toxicity for males:

None of the ingredients is listed.

· Chemicals known to cause developmental toxicity:

None of the ingredients is listed.

#### · Carcinogenic categories

· EPA (Environmental Protection Agency)

None of the ingredients is listed.

• TLV (Threshold Limit Value established by ACGIH)

13760-80-0 ytterbium trifluoride

· NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

· OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

### · Product related hazard informations:

Medical devices as defined in Directive 93/42/EEC and which are invasive or used in direct physical contact with the human body, are exempted from the provisions of Dangerous Substance Directive 1999/45/ EC usually if they are in the finished state and intended for the final user.

#### · Hazard symbols:



Irritant

- *Hazard-determining components of labeling: urethane dimethacrylate*
- · Risk phrases:
- May cause sensitisation by skin contact.
- Safety phrases: Avoid contact with
- Avoid contact with skin. Wear suitable gloves. This material and its container must be disposed of as hazardous waste.
- · National regulations:
- Other regulations, limitations and prohibitive regulations The product is a medical device according to the Directive 93/42/EEC. This product is classified as a medical device under US and Canadian regulations and has been reviewed by the US Food and Drug Administration and Health Canada.
- · Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

USA

Printing date 09/19/2013

Version number 5

Reviewed on 09/19/2013

# Trade name: Tetric EvoCeram

(Contd. of page 7)

# **16 Other information**

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.