



# ISODAM<sup>®</sup>

## Polyisoprene Dental Dam



### Product Description

With a calibre of experience in producing dental dams for leading dentistry institutions, Four D Rubber has developed the most advanced dental dam on the market; Isodam<sup>®</sup>

The primary function of Isodam<sup>®</sup> is to isolate specific areas of the mouth during procedures, reduce the risk of infection, and provide comfort to the patient.

The Isodam<sup>®</sup> provides exceptional stretch, strength, and tear-resistance while eliminating the potential for allergic reactions associated with latex products.

### Key Features

#### Superior Elasticity

The polyisoprene material used to create the Isodam<sup>®</sup> allows for an incredible up to 1200% elongation that virtually eliminates tearing even with a punched dam.

#### Size and Gauge Options

The Isodam<sup>®</sup> is available in medium and heavy gauges and 5" x 5" or 6" x 6" squares, ensuring the best non-latex dental dam option is available for any procedure.

#### Economy Package

For those that use Isodam<sup>®</sup> on a regular basis, our economy package holds 5 times more dams than our standard package, creating less waste and greater savings.

Visit [www.fourdrubber.com](http://www.fourdrubber.com) to learn more or request a sample



# ISODAM<sup>®</sup>

## Polyisoprene Dental Dam

### Product Specifications



- **Synthetic Polyisoprene, not made with natural rubber latex**
- **Premium Light Blue Colour**
- **100% Non-Toxic Pigment**
- **Low Powdered**
- **Non-Sterile/Single Use**

### Physical Properties

Tested under ISO37 Standard

<b>Shape</b>	Square cut
<b>Tensile Strength (MPa)</b>	≥20
<b>Elongation at Break (%)</b>	≥1050%
<b>Modulus at 500% Extension (MPa)</b>	≥1.6
<b>Tear Strength (N/mm)</b>	≥50

### Dimensions and Packaging

#### Standard Package

Code	Size	Colour	Weight	Thickness	Pcs per Box	Boxes per Inner Carton	Boxes per Outer Carton
ISO01800520	5"x 5" (127mm x 127mm)	Blue	Medium	(0.18mm-0.23mm)	20	10	60
ISO02300520	5"x 5" (127mm x 127mm)	Blue	Heavy	(0.23mm-0.28mm)	20	10	60
ISO01800615	6"x 6" (152mm x 152mm)	Blue	Medium	(0.18mm-0.23mm)	15	10	60
ISO02300615	6"x 6" (152mm x 152mm)	Blue	Heavy	(0.23mm-0.28mm)	15	10	60

#### Economy Package

Code	Size	Colour	Weight	Thickness	Pcs per Box	Boxes per Inner Carton	Boxes per Outer Carton
ISO01800510	5"x 5" (127mm x 127mm)	Blue	Medium	(0.18mm-0.23mm)	100	5	30
ISO02300510	5"x 5" (127mm x 127mm)	Blue	Heavy	(0.23mm-0.28mm)	100	5	30
ISO01800675	6"x 6" (152mm x 152mm)	Blue	Medium	(0.18mm-0.23mm)	75	5	30
ISO02300675	6"x 6" (152mm x 152mm)	Blue	Heavy	(0.23mm-0.28mm)	75	5	30