



Advancing the Quality
of Dental Care Products



Advance Orthodontic Products from

MidwestDental[®]

EQUIPMENT • SERVICES • SUPPLIES

(800) 766-2025 • MidwestDental.com

Log in to
see your
pricing!



Orthodontics Reference Catalog

4

Advance Precision Series
Metal Brackets



6

Advance Premium Precision Series
Metal Brackets



8

Advance Pure Ice
Ceramic Brackets



10

Advance Pure Ice-S
Aesthetic Brackets



12

SLB Precision Series
Self-Ligating Metal Brackets



14

Orthodontic Wires



16

Tubes



17

Orthodontic Accessories






18

Elastics



Guide of Orthodontic Specifications



Precision Series

TECHNOLOGY, PRECISION AND EFFICIENCY

Advance Precision Series Metal Brackets are produced using the MIM (metal-injection molding) process that combines all quality and technology requirements in a single piece.

Excellent retention, prescription precision, and an anatomical design provide patients with a more efficient and comfortable treatment.

The Advance Precision Series brackets feature a patented laser texturing process on the base to increase retention, thus ensuring greater adhesion efficiency.

LOW PROFILE DESIGN

Superior aesthetics and patient comfort.

BETTER FINISHING AND HIGH PRECISION

The MIM manufacturing process results in an extremely precise prescription and an exceptional surface finish.

IDENTIFICATION COLOR

Easily identify brackets with the color-coded system.

BETTER FITTING

Advance Precision's anatomical design provides better retention during bonding phase.

BETTER RETENTION

Bases are laser textured and include retaining pins to help increase retention and adhesion.

EASIER LIGATURE PLACEMENT

Twin wings design for easier placement of ligatures, better rotational control, and less friction on the slot edges.

SMOOTH ROUNDED HOOKS AND EDGES

Help ensure that patients have a more comfortable experience.

DIFFERENT PRESCRIPTIONS

Available in Roth 0.022" & MBT 0.022".



ROTH 0.022"

Precision Series



PKG/10 **\$24.95**
5+ **\$22.95**

CASE



(605-2660) Roth (.022")

PKG/20 **\$49.95**
5+ **\$44.95**

ROTH 0.022" | SLOT 0.022"x 0.028" (0.56 x 0.71mm)

	Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
	Angulation	0°	0°	11°	9°	5°	5°	9°	11°	0°	0°
	Hook	✓	✓	✓	-	-	-	-	✓	✓	✓
	Width (mm)	2.8	2.8	2.8	2.8	3.0	3.0	2.8	2.8	2.8	2.8
	MWD # ▶	605-2664	605-2664	605-2663	605-2662	605-2661	605-2665	605-2666	605-2667	605-2668	605-2668
	Advance # ▶	ARO106H	ARO106H	ARO104H	ARO102	ARO100	ARO101	ARO103	ARO105H	ARO107H	ARO107H
		5 15	4 14	3 13	2 12	1 11	1 21	2 22	3 23	4 24	5 25
		U45R	U45R	U3R	U2R	U1R	U1L	U2L	U3L	U45L	U45L
Maxillary											
	Upper Right										
	Lower Right										
Mandibular											
	Lower Right										
	Lower Left										
		45 5	44 4	43 3	42 2	41 1	31 1	32 2	33 3	34 4	35 5
		L5R	L4R	L3R	L12RL	L12RL	L12RL	L12RL	L3L	L4L	L5L
	Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
	Angulation	0°	0°	5°	0°	0°	0°	0°	5°	0°	0°
	Hook	✓	✓	✓	-	-	-	-	✓	✓	✓
	Width (mm)	2.8	2.8	2.8	2.4	2.4	2.4	2.4	2.8	2.8	2.8
	MWD # ▶	605-2672	605-2671	605-2670	605-2669	605-2669	605-2669	605-2669	605-2673	605-2674	605-2675
	Advance # ▶	ARO216H	ARO214H	ARO212H	ARO210	ARO210	ARO210	ARO210	ARO213H	ARO215H	ARO217H

MBT 0.022"

MBT 0.022" | SLOT 0.022"x 0.028" (0.56 x 0.71mm)



PKG/10 **\$24.95**
5+ **\$22.95**

CASE



(605-2680) MBT (.022")

PKG/20 **\$49.95**
5+ **\$44.95**

	Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
	Angulation	0°	0°	0°	8°	4°	4°	8°	0°	0°	0°
	Hook	✓	✓	✓	-	-	-	-	✓	✓	✓
	Width (mm)	2.8	2.8	2.7	2.8	2.8	2.8	2.8	2.7	2.8	2.8
	MWD # ▶	605-2684	605-2684	605-2683	605-2682	605-2681	605-2685	605-2686	605-2687	605-2688	605-2688
	Advance # ▶	AMO306H	AMO306H	AMO304H	AMO302	AMO300	AMO301	AMO303	AMO305H	AMO307H	AMO307H
		5 15	4 14	3 13	2 12	1 11	1 21	2 22	3 23	4 24	5 25
		U45R	U45R	U3R	U2R	U1R	U1L	U2L	U3L	U45L	U45L
Maxillary											
	Upper Right										
	Lower Right										
Mandibular											
	Lower Right										
	Lower Left										
		45 5	44 4	43 3	42 2	41 1	31 1	32 2	33 3	34 4	35 5
		L5R	L4R	L3R	L12RL	L12RL	L12RL	L12RL	L3L	L4L	L5L
	Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
	Angulation	2°	2°	3°	0°	0°	0°	0°	3°	2°	2°
	Hook	✓	✓	✓	-	-	-	-	✓	✓	✓
	Width (mm)	2.8	2.8	2.7	2.4	2.4	2.4	2.4	2.7	2.8	2.8
	MWD # ▶	605-2692	605-2691	605-2690	605-2689	605-2689	605-2689	605-2689	605-2693	605-2694	605-2695
	Advance # ▶	AMO416H	AMO414H	AMO412H	AMO410	AMO410	AMO410	AMO410	AMO413H	AMO415H	AMO417H

Premium Precision Series

PRECISION, EFFICIENCY, AND GREATER COMFORT!

Advance Premium Precision Series is a classic, 2-piece design, welded by an exclusive technique using 17/4 PH surgical steel to meet the highest quality standards.

Produced by a fully automated and high-technology CNC machining process that provides 100% torque and angulation accuracy.

An 80 gauge micro-etched mesh base ensures maximum bond strength.



DIFFERENTIALS:

- MADE OF 17/4 STEEL
- 80 GAUGE MICRO-ETCHED MESH BASE
- INNOVATIVE LASER WELDING PROCESS
- LOW PROFILE DESIGN PROVIDES GREATER PATIENT COMFORT
- ANATOMICALLY DESIGNED
- PRECISE TORQUE AND ANGULATION
- LASER MARKING FOR SIMPLIFIED BRACKET IDENTIFICATION
- MIDLINE MARKING FOR EASY POSITIONING
- AVAILABLE IN 2 DIFFERENT PRESCRIPTIONS: ROTH 0.022" AND MBT 0.022"



ROTH 0.022"

Premium Precision Series



PKG/5 **\$17.95**
5+ **\$16.95**

CASE



(605-2700) Roth (.022")

PKG/20 **\$69.95**
5+ **\$62.95**

ROTH | STANDARD 0.022" | SLOT 0.022"x 0.028" (0.56 x 0.71mm)

		-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
		0°	0°	11°	9°	5°	5°	9°	11°	0°	0°
		✓	✓	✓	-	-	-	-	✓	✓	✓
		2.8	2.8	2.8	2.8	3.1	3.1	2.8	2.8	2.8	2.8
	MWD # ▶	605-2704	605-2704	605-2703	605-2702	605-2701	605-2705	605-2706	605-2707	605-2708	605-2708
	Advance # ▶	PR0106H	PR0106H	PR0104H	PR0102	PR0100	PR0101	PR0103	PR0105H	PR0107H	PR0107H
		5 15	4 14	3 13	2 12	1 11	1 21	2 22	3 23	4 24	5 25
		U45R	U45R	U3R	U2R	U1R	U1L	U2L	U3L	U45L	U45L
	Maxillary										
		Upper Right									Upper Left
		Lower Right									Lower Left
	Mandibular										
		L5R	L4R	L3R	L12RL	L12RL	L12RL	L12RL	L3L	L4L	L5L
		45 5	44 4	43 3	42 2	41 1	31 1	32 2	33 3	34 4	35 5
		-22°	-17°	-11°	-1°	-1°	-1°	-1°	-11°	-17°	22°
		0°	0°	5°	0°	0°	0°	0°	5°	0°	0°
		✓	✓	✓	-	-	-	-	✓	✓	✓
		2.8	2.8	2.8	2.3	2.3	2.3	2.3	2.8	2.8	2.8
	MWD # ▶	605-2712	605-2711	605-2710	605-2709	605-2709	605-2709	605-2709	605-2713	605-2714	605-2715
	Advance # ▶	PR0216H	PR0214H	PR0212H	PR0210	PR0210	PR0210	PR0210	PR0213H	PR0215H	PR0217H

MBT 0.022"

MBT 0.022" | SLOT 0.022"x 0.028" (0.56 x 0.71mm)



PKG/5 **\$17.95**
5+ **\$16.95**

CASE



(605-2720) MBT (.022")

PKG/20 **\$69.95**
5+ **\$62.95**

		-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
		0°	0°	0°	8°	4°	4°	8°	0°	0°	0°
		✓	✓	✓	-	-	-	-	✓	✓	✓
		2.8	2.8	2.8	2.8	3.1	3.1	2.8	2.8	2.8	2.8
	MWD # ▶	605-2724	605-2724	605-2723	605-2722	605-2721	605-2725	605-2726	605-2727	605-2728	605-2728
	Advance # ▶	PM0306H	PM0306H	PM0304H	PM0302	PM0300	PM0301	PM0303	PM0305H	PM0307H	PM0307H
		5 15	4 14	3 13	2 12	1 11	1 21	2 22	3 23	4 24	5 25
		U45R	U45R	U3R	U2R	U1R	U1L	U2L	U3L	U45L	U45L
	Maxillary										
		Upper Right									Upper Left
		Lower Right									Lower Left
	Mandibular										
		L5R	L4R	L3R	L12RL	L12RL	L12RL	L12RL	L3L	L4L	L5L
		45 5	44 4	43 3	42 2	41 1	31 1	32 2	33 3	34 4	35 5
		-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
		2°	2°	3°	0°	0°	0°	0°	3°	2°	2°
		✓	✓	✓	-	-	-	-	✓	✓	✓
		2.8	2.8	2.8	2.5	2.5	2.5	2.5	2.8	2.8	2.8
	MWD # ▶	605-2732	605-2731	605-2730	605-2729	605-2729	605-2729	605-2729	605-2733	605-2734	605-2735
	Advance # ▶	PM0416H	PM0414H	PM0412H	PM0420	PM0420	PM0420	PM0420	PM0413H	PM0415H	PM0417H

Pure Ice

Polycrystalline

COMFORT, AESTHETICS, AND OPAQUE

Advance Pure Ice Brackets are produced using CIM (ceramic injection molding) technology and made with the purest polycrystalline aluminum oxide, which ensures excellent aesthetic results from the beginning to the end of treatment. It's the bracket of choice when you desire an aesthetic treatment with clinical performance, comfort and safety.

COMFORT AND AESTHETICS

- Low profile, rounded edges and ball hooks for patient comfort
- Resistant to pigmentation
- Made of 100% polycrystalline aluminum oxide.

COLOR-CODED DESIGN

- Color coded for easy identification and bracket placement

RETENTION AND FIT

- Mesiodistal and occluso-gingival curvatures allow for precise fitting to the tooth
- Exclusive laser texturing provides up to 30% higher retention when compared to traditional blasting
- Prevents adhesive leakage and facilitates removal of excess adhesive

ACCURACY AND HIGH PERFORMANCE:

- Strong ball hooks in the canines and premolars ensure both comfort and better fitting of elastics and accessories
- Twin wings design for easier placement of ligatures, better rotational control, and less friction on the slot edges
- Smooth slot design brings better sliding mechanics with limited friction
- Designed with torque in the base for greater accuracy during treatment



ROTH 0.022"

Pure Ice

ROTH 0.022" | SLOT 0.022"x 0.028" (0.56 x 0.71mm)



PKG/5 **\$52.95**
10+ **\$48.95**

CASE



(605-2590) Roth (.022")

PKG/20 **\$179.95**
10+ **\$165.95**

	Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
	Angulation	0°	0°	11°	9°	5°	5°	9°	11°	0°	0°
	Hook	✓	✓	✓	-	-	-	-	✓	✓	✓
	Width (mm)	3.2	3.2	3.7	3.2	3.7	3.7	3.2	3.7	3.2	3.2
	MWD # ▶	605-2594	605-2594	605-2593	605-2592	605-2591	605-2595	605-2596	605-2597	605-2598	605-2598
	Advance # ▶	APIRUR4/SH	APIRUR4/SH	APIRUR3H	APIRUR2	APIRUR1	APIRUL1	APIRUL2	APIRUL3H	APIRUL4/SH	APIRUL4/SH
		5 15 U4SR	4 14 U4SR	3 13 U3R	2 12 U2R	1 11 U1R	1 21 U1L	2 22 U2L	3 23 U3L	4 24 U4SL	5 25 U4SL
Maxillary											
		Upper Right Lower Right								Upper Left Lower Left	
Mandibular											
		L5R 45 5	L4R 44 4	L3R 43 3	L12RL 42 2	L12RL 41 1	L12RL 31 1	L12RL 32 2	L3L 33 3	L4L 34 4	L5L 35 5
	Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
	Angulation	0°	0°	5°	0°	0°	0°	0°	5°	0°	0°
	Hook	✓	✓	✓	-	-	-	-	✓	✓	✓
	Width (mm)	3.2	3.2	3.7	2.9	2.9	2.9	2.9	3.7	3.2	3.2
	MWD # ▶	605-2605	605-2604	605-2603	605-2599	605-2599	605-2599	605-2599	605-2600	605-2601	605-2602
	Advance # ▶	APIRLR5H	APIRLR4H	APIRLR3H	APIRL12UN	APIRL12UN	APIRL12UN	APIRL12UN	APIRL13H	APIRL14H	APIRL15H

MBT 0.022"

MBT 0.022" | SLOT 0.022"x 0.028" (0.56 x 0.71mm)



PKG/5 **\$52.95**
10+ **\$48.95**

CASE



(605-2570) MBT (.022")

PKG/20 **\$179.95**
10+ **\$165.95**

	Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
	Angulation	0°	0°	8°	8°	4°	4°	8°	8°	0°	0°
	Hook	✓	✓	✓	-	-	-	-	✓	✓	✓
	Width (mm)	3.2	3.2	3.7	3.2	3.7	3.7	3.2	3.7	3.2	3.2
	MWD # ▶	605-2574	605-2574	605-2573	605-2572	605-2571	605-2575	605-2576	605-2577	605-2578	605-2578
	Advance # ▶	APIMUR4/SH	APIMUR4/SH	APIMUR3H	APIMUR2	APIMUR1	APIMUL1	APIMUL2	APIMUL3H	APIMUL4/SH	APIMUL4/SH
		5 15 U4SR	4 14 U4SR	3 13 U3R	2 12 U2R	1 11 U1R	1 21 U1L	2 22 U2L	3 23 U3L	4 24 U4SL	5 25 U4SL
Maxillary											
		Upper Right Lower Right								Upper Left Lower Left	
Mandibular											
		L5R 45 5	L4R 44 4	L3R 43 3	L12RL 42 2	L12RL 41 1	L12RL 31 1	L12RL 32 2	L3L 33 3	L4L 34 4	L5L 35 5
	Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
	Angulation	2°	2°	3°	0°	0°	0°	0°	3°	2°	2°
	Hook	✓	✓	✓	-	-	-	-	✓	✓	✓
	Width (mm)	3.2	3.2	3.7	2.9	2.9	2.9	2.9	3.7	3.2	3.2
	MWD # ▶	605-2585	605-2584	605-2583	605-2579	605-2579	605-2579	605-2579	605-2580	605-2581	605-2582
	Advance # ▶	APIMLR5H	APIMLR4H	APIMLR3H	APIML12UN	APIML12UN	APIML12UN	APIML12UN	APIML13H	APIML14H	APIML15H

Pure Ice-S

Poly-Sapphire

COMFORT, AESTHETICS, AND TRANSPARENT

Advance Pure Ice-S Brackets are made with polycrystalline aluminum oxide, and are transparent and practically invisible.

They are one of the most innovative brackets in traditional cosmetic and orthodontic dentistry. Resistant to pigmentation, Advance Pure Ice-S brackets maintain their transparency until the end of the treatment with no deterioration or color change.



COLOR-CODED

- Color coding on the tie wings help ensure correct placement of the bracket during bonding and color can be easily removed

COMFORT AND AESTHETICS

- Made of 100% polycrystalline aluminum oxide
- Excellent aesthetic results
- Low-profile brackets with rounded edges and spherical ball hooks
- Resistant to wear, pigmentation and color changes

BONDING

- Anatomically-designed base helps create a perfect fit to the tooth surface
- Retaining edges on the base prevent leakage and facilitate the removal of excess adhesive
- Patented laser texturing of the bracket's base provides 30% greater retention

BASE WITH RETAINING EDGES

- Retention pins on the bracket help increase the bonding area

ACCURACY AND PERFORMANCE

- Anatomically-designed spherical hooks on the cuspids and premolars are not only strong, but also easier to use with elastics and accessories
- Carefully designed for more rotational control and less friction



ROTH 0.022"



PKG/5 **\$57.95**
10+ **\$53.95**

CASE



(605-2620) Roth (.022")
PKG/20 **\$189.95**
10+ **\$174.95**

Pure Ice-S

ROTH 0.022" | SLOT 0.022"x 0.028" (0.56 x 0.71mm)

	Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
	Angulation	0°	0°	11°	9°	5°	5°	9°	11°	0°	0°
	Hook	✓	✓	✓	-	-	-	-	✓	✓	✓
	Width (mm)	3.2	3.2	3.7	3.2	3.7	3.7	3.2	3.7	3.2	3.2
	MWD # ▶	605-2624	605-2624	605-2623	605-2622	605-2621	605-2625	605-2626	605-2627	605-2628	605-2628
	Advance # ▶	APISRUR4/SH	APISRUR4/SH	APISRUR3H	APISRUR2	APISRUR1	APISRUL1	APISRUL2	APISRUL3H	APISRUL4/SH	APISRUL4/SH
		5 15 U4SR	4 14 U4SR	3 13 U3R	2 12 U2R	1 11 U1R	1 21 U1L	2 22 U2L	3 23 U3L	4 24 U4SL	5 25 U4SL
Maxillary											
		LSR	L4R	L3R	L12RL	L12RL	L12RL	L12RL	L3L	L4L	LSL
		45 5	44 4	43 3	42 2	41 1	31 1	32 2	33 3	34 4	35 5
Mandibular											
		LSR	L4R	L3R	L12RL	L12RL	L12RL	L12RL	L3L	L4L	LSL
		45 5	44 4	43 3	42 2	41 1	31 1	32 2	33 3	34 4	35 5
	Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
	Angulation	0°	0°	5°	0°	0°	0°	0°	5°	0°	0°
	Hook	✓	✓	✓	-	-	-	-	✓	✓	✓
	Width (mm)	3.2	3.2	3.7	2.9	2.9	2.9	2.9	3.7	3.2	3.2
	MWD # ▶	605-2635	605-2634	605-2633	605-2629	605-2629	605-2629	605-2629	605-2630	605-2631	605-2632
	Advance # ▶	APISRLR5H	APISRLR4H	APISRLR3H	APISRL12UN	APISRL12UN	APISRL12UN	APISRL12UN	APISRL13H	APISRL14H	APISRL15H

MBT 0.022"



PKG/5 **\$57.95**
10+ **\$53.95**

CASE



(605-2640) MBT (.022")
PKG/20 **\$189.95**
10+ **\$174.95**

MBT 0.022" | SLOT 0.022"x 0.028" (0.56 x 0.71mm)

	Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
	Angulation	0°	0°	8°	8°	4°	4°	8°	8°	0°	0°
	Hook	✓	✓	✓	-	-	-	-	✓	✓	✓
	Width (mm)	3.2	3.2	3.7	3.2	3.7	3.7	3.2	3.7	3.2	3.2
	MWD # ▶	605-2644	605-2644	605-2643	605-2642	605-2641	605-2645	605-2646	605-2647	605-2648	605-2648
	Advance # ▶	APISMUR4/SH	APISMUR4/SH	APISMUR3H	APISMUR2	APISMUR1	APISMUL1	APISMUL2	APISMUL3H	APISMUL4/SH	APISMUL4/SH
		5 15 U4SR	4 14 U4SR	3 13 U3R	2 12 U2R	1 11 U1R	1 21 U1L	2 22 U2L	3 23 U3L	4 24 U4SL	5 25 U4SL
Maxillary											
		LSR	L4R	L3R	L12RL	L12RL	L12RL	L12RL	L3L	L4L	LSL
		45 5	44 4	43 3	42 2	41 1	31 1	32 2	33 3	34 4	35 5
Mandibular											
		LSR	L4R	L3R	L12RL	L12RL	L12RL	L12RL	L3L	L4L	LSL
		45 5	44 4	43 3	42 2	41 1	31 1	32 2	33 3	34 4	35 5
	Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
	Angulation	2°	2°	3°	0°	0°	0°	0°	3°	2°	2°
	Hook	✓	✓	✓	-	-	-	-	✓	✓	✓
	Width (mm)	3.2	3.2	3.7	2.9	2.9	2.9	2.9	3.7	3.2	3.2
	MWD # ▶	605-2655	605-2654	605-2653	605-2649	605-2649	605-2649	605-2649	605-2650	605-2651	605-2652
	Advance # ▶	APISMLR5H	APISMLR4H	APISMLR3H	APISML12UN	APISML12UN	APISML12UN	APISML12UN	APISML13H	APISML14H	APISML15H

SLB Precision Series

INTERACTIVE, DURABLE, AND COMFORTABLE

Advance SLB Precision Series Brackets have a nickel-titanium clip that is designed for opening and closing without the use of special instruments. The elastic nature and memory of the clip allows it to return to its original position without suffering permanent deformations. Together, the clip and wire exert light, continuous force, helping ensure the best control of tooth movements.

PRECISION MANUFACTURED FOR RESILIENCE

Manufactured using modern metal injection molding with high-strength stainless steel, this bracket is guaranteed to have the necessary resilience with an excellent finish.

ROUNDED, COMFORTABLE DESIGN

The edge-free design is rounded and polished to provide maximum comfort for the patient.

EASY-TO-HANDLE NITI CLIP


The clip is manufactured using a nickel-titanium alloy that prevents permanent deformations through its elastic memory to ensure greater durability and resistance. It can also be opened and closed using only an explorer or ligature applicator.

EXCELLENT SLIDING MECHANICS

Chamfers at the slot exits help to reduce friction and minimize the effects of binding and notching, which are harmful during orthodontic treatment.

MIDLINE MARKING

Permanent marking on the midline allows for easy bracket positioning during the bonding process.




BETTER FIT

An anatomically-designed base with mesiodistal and occlusal gingival curvatures makes for a great fit to the tooth surface.

BETTER RETENTION

The base is designed with pins, and undergoes a laser texturing process, which gives the brackets up to 30% greater retention, ensuring more efficient adhesion and safety to the orthodontist. The retention edges prevent the resin from leaking and facilitate the removal of excess adhesive, keeping it from obstructing the bracket clip.




Intraoral Archwire

Super-Elastic


NiTi

ADVANCE SUPER-ELASTIC NiTi

Preformed archwires designed to provide alignment, rotation and leveling with great safety and comfort for your patient. Advance Super-Elastic NiTi brings excellent shape memory and consistent forces to provide predictable treatment results. Due to their high flexibility coefficient, these NiTi archwires correct crowding without the need for bending. Advance Super-Elastic rectangular archwires are very efficient for aligning and leveling the brackets.



NATURAL



NATURAL 3

ADVANCE SUPER-ELASTIC NiTi ARCHWIRE – NATURAL

DIMENSIONS	UPPER	LOWER	PRICES
○ 0.012"	605-2488	605-2493	PKG/10 \$23.95 5+ \$21.95
○ 0.014"	605-2489	605-2494	
○ 0.016"	605-2490	605-2495	
○ 0.018"	605-2491	605-2496	
○ 0.020"	605-2492	605-2497	
□ 0.016"x0.016"	605-2498	605-2504	PKG/10 \$28.95 5+ \$26.95
□ 0.016"x0.022"	605-2499	605-2505	
□ 0.017"x0.025"	605-2500	605-2506	
□ 0.018"x0.025"	605-2501	605-2507	
□ 0.019"x0.025"	605-2502	605-2508	
□ 0.021"x0.025"	605-2503	605-2509	

○ ROUND □ SQUARE □ RECTANGULAR

ADVANCE SUPER-ELASTIC NiTi ARCHWIRE – NATURAL 3

DIMENSIONS	UPPER	LOWER	PRICES
○ 0.012"	605-2522	605-2527	PKG/10 \$23.95 5+ \$21.95
○ 0.014"	605-2523	605-2528	
○ 0.016"	605-2524	605-2529	
○ 0.018"	605-2525	605-2530	
○ 0.020"	605-2526	605-2531	
□ 0.016"x0.016"	605-2510	605-2516	PKG/10 \$28.95 5+ \$26.95
□ 0.016"x0.022"	605-2511	605-2517	
□ 0.017"x0.025"	605-2512	605-2518	
□ 0.018"x0.025"	605-2513	605-2519	
□ 0.019"x0.025"	605-2514	605-2520	
□ 0.021"x0.025"	605-2515	605-2521	

○ ROUND □ SQUARE □ RECTANGULAR

Intraoral Archwire

Thermo NiTi 35°

ADVANCE THERMO NiTi 35°

Advance Thermo NiTi 35° Archwires are thermally active at 35°C, which provides greater flexibility with lighter and more continuous force as compared to conventional NiTi alloys. Excellent sliding mechanics for predictable tooth movement.

This activation, which is above room temperature and close to body temperature, makes it easier for the archwire to adapt to the oral cavity.



ADVANCE THERMO NiTi 35° ARCHWIRE – NATURAL

DIMENSIONS	UPPER	LOWER	PRICES
○ 0.012"	605-2444	605-2449	PKG/10 \$28.95 5+ \$26.95
○ 0.014"	605-2445	605-2450	
○ 0.016"	605-2446	605-2451	
○ 0.018"	605-2447	605-2452	
○ 0.020"	605-2448	605-2453	
□ 0.016"x0.016"	605-2454	605-2460	PKG/10 \$34.95 5+ \$31.95
□ 0.016"x0.022"	605-2455	605-2461	
□ 0.014"x0.025"	605-2761	605-2760	
□ 0.017"x0.025"	605-2456	605-2462	
□ 0.018"x0.025"	605-2457	605-2463	
□ 0.019"x0.025"	605-2458	605-2464	
□ 0.021"x0.025"	605-2459	605-2465	

○ ROUND □ SQUARE □ RECTANGULAR

ADVANCE THERMO NiTi 35° ARCHWIRE – NATURAL 3

DIMENSIONS	UPPER	LOWER	PRICES
○ 0.012"	605-2422	605-2427	PKG/10 \$28.95 5+ \$26.95
○ 0.014"	605-2423	605-2428	
○ 0.016"	605-2424	605-2429	
○ 0.018"	605-2425	605-2430	
○ 0.020"	605-2426	605-2431	
□ 0.016"x0.016"	605-2432	605-2438	PKG/10 \$34.95 5+ \$31.95
□ 0.016"x0.022"	605-2433	605-2439	
□ 0.017"x0.025"	605-2434	605-2440	
□ 0.018"x0.025"	605-2435	605-2441	
□ 0.019"x0.025"	605-2436	605-2442	
□ 0.021"x0.025"	605-2437	605-2443	

○ ROUND □ SQUARE □ RECTANGULAR

Intraoral Archwire

Reverse

SE NiTi

ADVANCE REVERSE CURVE SE NiTi

Due to their differentiated shape, Advance Reverse Curve SE NiTi Intraoral Archwires are ideal for bite opening, vertical occlusion and torque control. They feature much greater resilience than stainless steel archwires and exert smooth, even forces for rapid tooth movement.



ADVANCE REVERSE CURVE SE NITI ARCHWIRE - NATURAL

DIMENSIONS	UPPER	LOWER	PRICES
○ 0.012"	605-2466	605-2471	PKG/10 \$42.95 5+ \$38.95
○ 0.014"	605-2467	605-2472	
○ 0.016"	605-2468	605-2473	
○ 0.018"	605-2469	605-2474	
○ 0.020"	605-2470	605-2475	
□ 0.016"x0.016"	605-2476	605-2482	PKG/10 \$49.95 5+ \$44.95
□ 0.016"x0.022"	605-2477	605-2483	
□ 0.017"x0.025"	605-2478	605-2484	
□ 0.018"x0.025"	605-2479	605-2485	
□ 0.019"x0.025"	605-2480	605-2486	
□ 0.021"x0.025"	605-2481	605-2487	

○ ROUND □ SQUARE □ RECTANGULAR

Intraoral Archwire

Copper Thermal NiTi

Copper 6%

ADVANCE COPPER THERMAL NiTi

Thermo-activated Advance Copper Thermal NiTi Archwires are 6% copper in composition and have an activation temperature of 35° C.

BENEFITS OF COPPER

- Promotes precise shape memory and greater resistance to permanent deformation compared to conventional alloys
- Enhances the thermal properties of the archwire and highlights its elastic capacity, providing greater performance and consistency in tooth movement



ADVANCE COPPER THERMAL NITI ARCHWIRE - NATURAL

DIMENSIONS	UPPER	LOWER	PRICES
○ 0.012"	605-2400	605-2405	PKG/10 \$54.95 5+ \$49.95
○ 0.014"	605-2401	605-2406	
○ 0.016"	605-2402	605-2407	
○ 0.018"	605-2403	605-2408	
□ 0.016"x0.022"	605-2411	605-2417	
□ 0.014"x0.025"	605-2410	605-2416	PKG/10 \$69.95 5+ \$62.95
□ 0.016"x0.025"	605-2412	605-2418	
□ 0.017"x0.025"	605-2413	605-2419	
□ 0.018"x0.025"	605-2414	605-2420	
□ 0.019"x0.025"	605-2415	605-2421	

○ ROUND □ RECTANGULAR

Intraoral Archwire

Stainless Steel

304 VAR

ADVANCE STAINLESS STEEL 304 VAR

Advance Stainless Steel 304 VAR archwires are manufactured with a superior Nickel-Chrome alloy (Stainless Steel) that provides excellent surface finish and polish, ensuring the high quality of the product. The archwires have the characteristics of surface gloss and good formability, weldability and strength.



ADVANCE STAINLESS STEEL ARCHWIRE - NATURAL

DIMENSIONS	UPPER	LOWER	PRICES
○ 0.012"	605-2544	605-2549	PKG/10 \$7.95 5+ \$6.95
○ 0.014"	605-2545	605-2550	
○ 0.016"	605-2546	605-2551	
○ 0.018"	605-2547	605-2552	
○ 0.020"	605-2548	605-2553	
□ 0.016"x0.016"	605-2532	605-2538	PKG/10 \$9.95 5+ \$8.95
□ 0.016"x0.022"	605-2533	605-2539	
□ 0.017"x0.025"	605-2534	605-2540	
□ 0.018"x0.025"	605-2535	605-2541	
□ 0.019"x0.025"	605-2536	605-2542	
□ 0.021"x0.025"	605-2537	605-2543	

○ ROUND □ SQUARE □ RECTANGULAR

Intraoral Archwire

Stainless Steel

304 VAR Braided

ADVANCE BRAIDED STAINLESS STEEL

Braided archwires are made of eight braided stainless-steel filaments, which provide greater flexibility and produce forces with a magnitude lower than those produced by solid steel archwires.



ADVANCE BRAIDED INTRAORAL ARCHWIRE - NATURAL

DIMENSIONS	UPPER	LOWER	PRICES
□ 0.016"x0.022"	605-2554	605-2559	PKG/10 \$34.95 5+ \$31.95
□ 0.017"x0.025"	605-2555	605-2560	
□ 0.018"x0.025"	605-2556	605-2561	
□ 0.019"x0.025"	605-2557	605-2562	
□ 0.021"x0.025"	605-2558	605-2563	

□ RECTANGULAR

ADVANCE PRECISION SERIES BONDING TUBES

Advance Precision Series Bonding Tubes are produced using the MIM (metal-injection molding) process. The tubes have an anatomically-designed, laser-etched base with mesiodistal and occlusogingival curvatures, buccal groove and smooth edges. Top notch for improved tube positioning, and a tapered tube entry for easier archwire insertion. Ball hook for patient comfort and adding elastics.



ROTH 0.022" SINGLE NON-CONVERTIBLE TUBE					BONDING	MBT 0.022" SINGLE NON-CONVERTIBLE TUBE					BONDING
USE	LEN.	TORQUE	DISTAL		0.022x0.028"	USE	LEN.	TORQUE	DISTAL		0.022x0.028"
Upper	7/6	4.0mm	-10°	8°	605-2800	Upper	7/6	4.0mm	-14°	10°	605-2804
	6/7	4.0mm	-10°	8°	605-2801		6/7	4.0mm	-14°	10°	605-2805
Lower	7/6	4.0mm	-25°	4°	605-2802	Lower	7/6	4.0mm	-20°	0°	605-2807
	6/7	4.0mm	-25°	4°	605-2803		6/7	4.0mm	-20°	0°	605-2806

PKG/10 **\$59.95**
5+ **\$53.95**

ADVANCE PREMIUM PRECISION SERIES BONDING TUBES

Advance Premium Precision Series Bonding Tubes are made of 17/4 PH surgical steel through an innovative laser-welding process. They have a low-profile design for greater comfort, and an anatomically-designed 80 gauge mesh base.



ROTH 0.022" SINGLE NON-CONVERTIBLE TUBE					BONDING	MBT 0.022" SINGLE NON-CONVERTIBLE TUBE					BONDING
USE	LEN.	TORQUE	DISTAL		0.022x0.028"	USE	LEN.	TORQUE	DISTAL		0.022x0.028"
Upper	7/6	3.2mm	-10°	10°	605-2820	Upper	7/6	3.2mm	-14°	10°	605-2824
	7/6	3.2mm	-10°	10°	605-2821		7/6	3.2mm	-14°	10°	605-2825
Lower	7/6	3.2mm	-25°	4°	605-2822	Lower	7/6	3.2mm	-25°	4°	605-2826
	7/6	3.2mm	-25°	4°	605-2823		7/6	3.2mm	-25°	4°	605-2827

ROTH 0.022" SINGLE CONVERTIBLE TUBE					BONDING	MBT 0.022" SINGLE CONVERTIBLE TUBE					BONDING
USE	LEN.	TORQUE	DISTAL		0.022x0.028"	USE	LEN.	TORQUE	DISTAL		0.022x0.028"
Upper	6	4 mm	-10°	10°	605-2850	Upper	6	4 mm	-14°	10°	605-2854
	6	4 mm	-10°	10°	605-2851		6	4 mm	-14°	10°	605-2855
Lower	6	4 mm	-25°	4°	605-2852	Lower	6	4 mm	-25°	4°	605-2856
	6	4 mm	-25°	4°	605-2853		6	4 mm	-25°	4°	605-2857

PKG/10 **\$64.95**
5+ **\$58.95**

ADVANCE PURE ICE BONDING TUBES

Designed for use with clear aligners only.

Advance Pure Ice Bonding Tubes have a large offset base to get close to the gingival surface. Due to their small size, the tubes allow bonding in tight areas such as partially-erupted second molars where only the mesiovestibular cusp is visible. They're made with the purest polycrystalline aluminum oxide, which ensures excellent aesthetic results. Advance Pure Ice Tubes are anatomically-designed and have comfortable, rounded edges. Designed for sectional treatment used in an aligner system.

0.022" ADVANCE PURE ICE TUBES					BONDING
PRESCRIPTION	LEN.	TORQUE	DISTAL		0.022x0.028"
Upper	UNIVERSAL	3.0 mm	-7°	0°	605-0051
	UNIVERSAL	3.0 mm	-7°	0°	605-0051
Lower	UNIVERSAL	3.0 mm	-7°	0°	605-0051
	UNIVERSAL	3.0 mm	-7°	0°	605-0051



PKG/2 **\$10.95**
5+ **\$9.95**

ADVANCE PRECISION SERIES BONDING BUTTONS

Stainless steel buttons used to retain elastics, ligatures and ligature wires. Laser-textured bases with retaining pins.

Pkg/10 \$26.95, 5+ \$24.95, 10+ \$22.95

(605-0065) Straight/Flat

(605-0066) Convex

(605-0064) Concave/Curved

STRAIGHT/FLAT



CONVEX



CONCAVE/CURVE



ADVANCE PURE ICE CERAMIC LINGUAL BUTTONS

Made with polycrystalline aluminum oxide for aesthetics, strength and durability. Universal design for both the upper and lower arches. Can also be used with conventional orthodontic mechanics or with clear aligners.



Pkg/2 \$34.95, 5+ \$32.95, 10+ \$30.95

(605-0050) Concave

ADVANCE PRECISION SERIES STRAIGHT BALL TIP HOOKS

For anchoring elastics, springs and metal ligature wires. Stainless steel. Rounded edges.



Pkg/10 \$29.95

(605-0068) Straight Ball Tip Hook

ADVANCE PRECISION SERIES CLOSED COIL SPRINGS

These nickel titanium springs are flexible, resilient, biocompatible, have great shape memory and longer activation intervals.



Closed Coil Springs

Pkg/4 \$19.95

(605-0055) 7mm

(605-0056) 9mm

(605-0057) 12mm

(605-0058) 15mm

Closed Coil Springs for Mini Screw

Pkg/4 \$21.95

(605-0060) 7mm

(605-0061) 9mm

(605-0062) 12mm

(605-0063) 15mm

ADVANCE PRECISION SERIES OPEN COIL SPRINGS

These nickel titanium springs are flexible, resilient, biocompatible and have great shape memory.



Open Coil Springs

Pkg/4 \$21.95

(605-0059) .75mm X 178mm

ADVANCE PRECISION SERIES ORTHO STOPS

Advance Precision Series stainless steel stops are crimpable and available in open (C-stop) and closed (tube) styles. Open stops are dead soft and easily fit up to a .018mm x .025mm archwire. Closed stops are 2mm long, fit up to a .21mm x .25mm archwire and slide into place using Weingart utility pliers for customized placement.



Pkg/10 \$29.95

(605-0069) Open Stop



Pkg/25 \$45.95

(605-0070) Closed Stop

ADVANCE CLEAR DISTALIZER CLASS II AESTHETIC CORRECTOR

Aesthetic, comfortable design that can be used for right and left dental arches, and with fixed braces or aligners. The canine base is made of polymer and the molar base is made of metal.

Pkg/1 pair \$189.95

23mm | (605-0052)

25mm | (605-0053)

27mm | (605-0054)

Contains: 1 pair of distalizers and 1 ruler.



Select the size you need using the included flexible millimeter ruler. Position the ruler between the central buccal groove of the first molar and the center of the buccal face of the canine.



23 25 27(mm)

Clear Distalizer
Rev. A 02.03.04.33 - 04/24





ADVANCE ELASTICS WILDLIFE OF YELLOWSTONE PARK

Introduce a little wild into your patient's next ortho visit with Advance Elastics featuring the Wildlife of Yellowstone Park! Memorable animal shapes printed on the resealable bags help staff quickly identify which elastic size and strength is needed. Advance Elastics have excellent memory, calibrated forces and a mild mint scent.

- Indicated for intercuspation controls, Classes II and III dental correction, and midline correction
- Calibrated forces and excellent shape memory
- Resealable bags have a printed odontogram as an usage guide for the patient and a free plastic inserter
- Elastics have a mild mint aroma

ELASTICS | WILDLIFE OF YELLOWSTONE PARK

DIMENSIONS	1/8" 3.2mm	3/16" 4.8mm	1/4" 6.3mm	5/16" 7.9mm	3/8" 9.5mm
FORCE	Red Fox	Sandhill Crane	Coyote	Osprey	Otter
LIGHT 2.5 OZ 71g					
ITEM NUMBER	605-3772	605-3775	605-3778	605-3781	605-3784
FORCE	Bobcat	Wolf	Wolverine	Bald Eagle	Mountain Goat
MEDIUM 4.5 OZ 128g					
ITEM NUMBER	605-3773	605-3776	605-3779	605-3782	605-3785
FORCE	Grizzly	Bison	Moose	Bighorn Sheep	Elk
HEAVY 6.5 OZ 184g					
ITEM NUMBER	605-3774	605-3777	605-3780	605-3783	605-3786

Contains latex.




EACH **\$43.95**
5+ **\$40.95**
10+ **\$38.95**

Each box contains 50, 100-count bags (For a total of 5,000 elastics!)



ADVANCE CHAINS

Advance Chains have excellent memory, color stability and are non-toxic. Chains are available in a variety of different colors, including aesthetically-pleasing crystal, pearl and pearl blue.



ROLL/5 FT **\$13.95**, 5+ **\$12.95**, 10+ **\$11.95**

ROLL/15 FT **\$23.95**, 5+ **\$21.95**, 10+ **\$20.95**


COLORS	CONT. (5ft)	SHORT (5ft)	CONT. (15ft)	SHORT (15ft)
BLUE	605-3790	605-3813	605-3700	605-3723
RED	605-3791	605-3814	605-3701	605-3724
PINK	605-3792	605-3815	605-3702	605-3725
BLACK	605-3793	605-3816	605-3703	605-3726
SEA GREEN	605-3794	605-3817	605-3704	605-3727
CRYSTAL GREEN	605-3795	605-3818	605-3705	605-3728
GREY	605-3796	605-3819	605-3706	605-3729
CRYSTAL	605-3797	605-3820	605-3707	605-3730
CRYSTAL BLUE	605-3798	605-3821	605-3708	605-3731
NAVY BLUE	605-3799	605-3822	605-3709	605-3732
PURPLE	605-3800	605-3823	605-3710	605-3733
WHITE	605-3801	605-3824	605-3711	605-3734

COLORS	CONT. (5ft)	SHORT (5ft)	CONT. (15ft)	SHORT (15ft)
LIGHT PURPLE	605-3802	605-3825	605-3712	605-3735
LEMON YELLOW	605-3803	605-3826	605-3713	605-3736
SILVER	605-3804	605-3827	605-3714	605-3737
BABY BLUE	605-3805	605-3828	605-3715	605-3738
BABY PINK	605-3806	605-3829	605-3716	605-3739
CRYSTAL RED	605-3807	605-3830	605-3717	605-3740
PEARL	605-3808	605-3831	605-3718	605-3741
PEARL BLUE	605-3809	605-3832	605-3719	605-3742
PEARL GREEN	605-3810	605-3833	605-3720	605-3743
CRYSTAL PURPLE	605-3811	605-3834	605-3721	605-3744
ORANGE	605-3812	605-3835	605-3722	605-3745

*Due to printing variations, colors may vary.

ADVANCE PATIENT TIES

Advance Patient Ties have excellent elasticity and memory, color stability and are non-toxic. Ties are available in a variety of different colors. Latex free.



COLORS	22 Ties (pkg/46)	11 Ties (pkg/92)	COLORS	22 Ties (pkg/46)	11 Ties (pkg/92)
BLUE	605-3746	605-3850	SILVER	605-3762	605-3864
RED	605-3747	605-3851	BABY BLUE	605-3764	605-3865
PINK	605-3748	605-3852	BABY PINK	605-3763	605-3866
BLACK	605-3749	605-3853	CRYSTAL RED	605-3765	605-3867
SEA GREEN	605-3750	605-3854	PEARL	605-3766	605-3868
CRYSTAL GREEN	605-3751	605-3855	PEARL BLUE	605-3767	605-3869
GREY	605-3752	605-3856	ASSORTED	605-3768	605-3870
CRYSTAL	605-3753	605-3857	PEARL GREEN	605-3769	605-3871
CRYSTAL BLUE	605-3754	605-3858	CRYSTAL PURPLE	605-3770	605-3872
NAVY BLUE	605-3755	605-3859	ORANGE	605-3771	605-3873
PURPLE	605-3756	605-3860	CORAL	605-3840	605-3874
WHITE	605-3757	605-3861	MAGENTA	605-3841	605-3875
LIGHT PURPLE	605-3760	605-3862	BEACHSIDE BLUE	605-3842	605-3876
LEMON YELLOW	605-3761	605-3863	NEON GREEN	605-3843	605-3877

*Due to printing variations, colors may vary.

PKG/46 STICKS OF 22 TIES
(For a total of 1,012 ties!)

EACH **\$26.95**
5+ **\$24.95**
10+ **\$22.95**

PKG/92 STICKS OF 11 TIES
(For a total of 1,012 ties!)

EACH **\$26.95**
5+ **\$24.95**
10+ **\$22.95**

ADVANCE SEPARATORS

Separating elastics used to expand interdental spaces. Radiopaque.

Pkg/1,000 **18.95**
5+ **17.95**
10+ **16.95**

(605-3836) Green (3/16")
(605-3837) Blue (5/32")



*Prices are subject to change without notice.



Introduce a little *wild* into your patient's
next ortho visit with **Advance Elastics**
featuring the Wildlife of Yellowstone Park!

Check out our
patient elastics,
chains and ties
on pages 18-19!



Advance Orthodontic Products from **MidwestDental®**
EQUIPMENT • SERVICES • SUPPLIES

(800) 766-2025 | MidwestDental.com

Not responsible for printing and typographical errors. Prices are subject to change without notice. Version 6_2024