

Dental Solutions from

**VOCO**  
THE DENTALISTS


**MidwestDental®**  
EQUIPMENT • SERVICES • SUPPLIES

**CORE** SEE PAGE 6

TOOTH-LIKE  
SURFACE  
HARDNESS!



**VARNISH** SEE BACK COVER



**BONDING** SEE PAGE 2



SEE PAGE 3 **CEMENTS**



**COMPOSITES** SEE PAGE 4



SEE PAGE 6 **CROWN & BRIDGE**

Order Today! (800) 766-2025 | [MidwestDental.com](http://MidwestDental.com)

**VOCO** Reference Catalog

Log in to  
see your  
pricing!

## Acrylics

### Hard Reline Material

#### Ufi Gel® hard C (VOCO)

Biocompatible direct hard denture relining material ideal for permanent chairside relining, total or partial and direct or indirect relining, and more! MMA-free. Bubble-free automix application. Esthetic and color stable.



#### (284-0004) Cartridge Intro Kit

Contains: 1 80g cartridge, 1 10ml adhesive, 1 Type 2 dispenser and 10 Type 8 mixing tips.

#### (284-0001) Cartridge Kit

Contains: 1 80g cartridge, 1 10ml bottle adhesive, and 10 Type 8 mixing tips.

#### (284-0002) Cartridge Refill

Contains: 1 80g cartridge and 10 Type 8 mixing tips.

#### (284-0003) Adhesive Refill

Contains: 1 10ml bottle adhesive.

#### Type 8 Mixing Tips (Pink)

(168-0049) Pkg/50

Non-stock

### Soft Reline Material

#### Ufi Gel® SC (VOCO)

Self-curing, silicone-based permanent soft reline material for total and partial dentures. Ideal for relieving pressure, damping the palatal vibrating line, cushioning sharp-edged alveolar processes, and supporting healing in implantology. Biocompatible, tasteless, and odorless. 1-minute processing of adhesive.



#### (284-0007) Intro Kit

Contains: 1 50ml cartridge, 1 10ml bottle adhesive, 1 10ml glazing base, 1 10ml glazing catalyst, 1 brush holder, 25 micro tims, 10 Type 8 mixing tips, 4 mixing pads, and 1 polisher disc.

#### (284-0006) Cartridge Refill

Contains: 1 50ml cartridge and 10 Type 8 mixing tips.

#### (284-0005) Adhesive Refill

Contains: 1 10ml bottle adhesive.

#### (284-0008) Base/Catalyst Refill

Contains: 1 10ml each of glazing base and catalyst.

#### (199-1895) Single Tim, pkg/100

#### Type 8 Mixing Tips (Pink)

(168-0049) Pkg/50

Non-stock

Non-stock

## Articulating

### Occlusal & Pressure Indicators

#### Fit Test® C&B (VOCO)

A fast diagnosis aid for all restorative applications. Ideal for the precise control and fit of marginal gaps, as well as adaptations of crowns, bridges, inlays/onlays, cast metal cores, metal denture bases and attachments. 45-second set time. Extremely tear resistant.



#### (284-0009) Syringe Kit

Contains: 1 5ml syringe and 10 Type 9 mixing tips.

#### Type 9 Mixing Tips, Short (Brown/White)

(160-0015) Pkg/50

## Bond and Etch

#### Futurabond® DC (VOCO)

Self-etch, dual-cured adhesive for all light-, dual-, and self-cured resin materials. Reinforced with nano particles for increased bonding value and marginal integrity. Moisture tolerant. Fast and easy one-coat application— apply, dry, cure in 35-seconds).



#### (080-0024) Liquid Set

Contains: 1 4ml Liquid 1, 1 4ml Liquid 2, 50 single tims, and 2 mixing pads.

#### (080-0025) .1ml SingleDose, pkg/50

#### (080-0026) .1ml SingleDose, pkg/200

#### (199-1895) Single Tims, pkg/100

#### Futurabond® M+ (VOCO)

A flexible and reliable solution that can be applied with or without phosphoric acid etchant that is suitable for direct and indirect restorations. Ideal for zirconium, aluminum oxide, and silicate ceramics without the need for additional primer. 35-second working time.



#### (080-0027) Bottle Kit

Contains: 1 5ml bottle, 50 single tims, and 2 mixing pads.

#### (080-0028) 5ml bottle, pkg/3

#### (080-0029) Dual Cure Activator, 2ml

#### (199-1895) Single Tims, pkg/100

#### Futurabond® NR (VOCO)

A nano-reinforced, single-step, self-etch adhesive. Indicated for direct restorations with all light-curing restorative materials for all classes of fillings. Features long-term marginal integrity, moisture tolerance, quick and easy application, and fluoride-release properties.




#### (080-0022) .1ml SingleDose, pkg/50

#### (080-0023) .1ml SingleDose, pkg/200

Non-stock

#### Futurabond® U (VOCO)

One universal adhesive for all your needs. Self-, selective-, or total-etch for all direct or indirect materials. Reinforced with nano particles for permanent bond strength. Bonds to all light-, dual-, and self-cure resin materials on metal, zirconia, aluminum oxide, or silicate ceramic without extra primer. Fast application with 35-second working time.



#### (080-0032) Trial Kit

Contains: 20 .1ml SingleDoses, 1 2ml syringe Vococid gel, and 10 Type 48 intra-oral tips.

#### (080-0030) .1ml SingleDose, pkg/50

#### (080-0031) .1ml SingleDose, pkg/200

#### Type 48 Intra-oral Tips

(160-0044) Pkg/50

Non-stock

#### Vococid Phosphoric Etchant Gel (VOCO)

35% orthophosphoric acid gel in a blue tint for visual control. Does not drip.



#### (080-0034) 2ml Syringe, pkg/5

Contains: 5 2ml syringes and 15 Type 48 intra-oral tips.

#### (080-0033) 5ml Syringe, each

Contains: 1 5ml syringe and 10 Type 42 intra-oral tips.

#### (199-1896) Type 42, pkg/30

#### (160-0044) Type 48, pkg/50

Non-stock

Non-stock

Non-stock



## Cements & Liners

### Cavity Liners & Bases

#### Calciur® (VOCO)

45% calcium hydroxide paste ideal for use under any lining or restorative material. Radiopaque with a high antimicrobial effect due to its high pH value of 12.5.

(144-0004) 2ml Syringe, each

Contains: 1 2ml syringe and 5 Type 47 intra-oral tips.

(160-0013) Type 47 Tips, pkg/30



#### Calcimol (VOCO)

Self-curing, radiopaque, 26% calcium hydroxide paste. Easy to apply, short setting time. Supports the formation of tertiary dentin.

(144-0005) Complete Kit

Contains: 1 13g base and 1 11g catalyst.



#### Calcimol LC (VOCO)

Light-curing, radiopaque, resin modified calcium ion-releasing base liner and pulp-capping material. High pH level of 9.5-11. Acid resistant.

(144-0006) Complete Kit

Contains: 2 2.5g syringes and 5 Type 41 intra-oral tips.

(520-2651) Type 41 Tips, pkg/100



#### Ionoseal® (VOCO)

Light-cured, glass ionomer composite liner with high fluoride release and biocompatibility. Radiopaque and acid resistant. Fast, 20-second light cure at 2mm.

(144-0007) Complete Kit

Contains: 3 2.5g syringes and Type 41 intra-oral tips.

(520-2651) Type 41 Tips, pkg/100



### Desensitizer

#### Admira® Protect (VOCO)

**CDT code is 9911 for desensitizing application!**

Light-cured, nano-filled desensitizer ideal for hypersensitive dentin and exposed cervical areas. High bond strength and wear resistance. No glutaraldehyde.

(144-0008) 4.5ml Bottle, each

Contains: 1 4.5ml bottle, 1 brush holder, and 50 micro tims.

(199-2029) Micro Tim, pkg/100



Non-stock

### Glass Ionomer Cements

#### Meron Glass Ionomer Luting Cement (VOCO)

Radiopaque with excellent wetting properties, low film thickness, high translucency, and ideal stability.

(160-0008) Standard Kit

Contains: 1 35g powder, 1 15ml liquid, and 1 green measuring spoon.

(160-0009) Mini Powder Kit

Contains: 1 15g powder, 1 7ml liquid, and 1 green measuring spoon.

(160-0010) Meron AC, pkg/50



Non-stock

#### Meron Plus QM Resin Modified Glass Ionomer Cement (VOCO)

Self-adhesive resin modified glass ionomer cement with twice the bond strength to dentin and twice the adhesive value to zirconium dioxide than other leading brands. Moisture and acid resistant with low film thickness for a precise marginal fit.

(160-0011) 8.5g Syringe, each

Contains: 1 8.5g syringe and 10 Type 10 mixing tips.

(160-0012) Pkg/3 8.5g syringes

Contains: 3 8.5g syringes and 30 Type 10 mixing tips.

Type 10 Mixing Tips, Long (Brown/White)

(160-0016) Pkg/50



### Permanent Cements

#### Bifix® QM (VOCO)

Universal, dual-cured, self-etch, radiopaque resin based adhesive cement with excellent adhesion to dentin and enamel. Ideal for inlays and onlays, crowns, veneers, pins and posts, and bridges.

(160-0020) Intro Kit

Contains: 1 10g QuickMix syringe in Universal, 1 4ml bottle each of FuturaBond DC liquid types 1 and 2, 5ml bottle Ceramic Bond, 10 Type 11 mixing tips, and 10 Type 4 intra-oral tips.

QuickMix Syringe Refills

Contains: 1 10g QuickMix syringe, 10 Type 11 mixing tips, and 10 Type 4 intra-oral tips.

(160-0021) Universal

(160-0022) Transparent

(160-0023) White Opaque

(160-0014) Type 4 Tips, pkg/50

(520-2655) Type 11 Tips, pkg/50



Non-stock

#### Bifix® SE (VOCO)

Universal, dual-cured, radiopaque self-adhesive resin cement compatible with common materials. Features high bond strength without an adhesive and high mechanical strength. Low solubility permits reliable cementation. Gel-phase cleanup, no refrigeration or hand mixing required.

(160-0024) Intro Kit

Contains: 3 5g QuickMix syringes in shades Universal, Transparent, and White Opaque, 15 Type 14 mixing tips, 15 Type 15 mixing tips and 15 Type 1 intra-oral tips.

QuickMix Syringe Refills

Contains: 1 10g QuickMix syringe, 10 Type 11 mixing tips, and 10 Type 4 intra-oral tips.

(160-0025) Universal

(160-0026) Transparent

(160-0027) White Opaque

Type 1 Intra-oral Tips

(160-0039) Pkg/50

Type 4 Intra-oral Tips

(160-0014) Pkg/50

Type 11 Mixing Tips

(520-2655) Pkg/50

Type 14 Mixing Tips (Brown/Orange)

(160-0018) Pkg/50

Type 15 Intra-oral Tips (Brown/Orange)

(160-0019) Pkg/50



Non-stock

Non-stock

Non-stock

## Temporary Cements

### Bifix® Temp (VOCO)

A dual-cured, transparent, composite-based temporary cement with strong adhesion properties. Easy to remove excess material. 5µm film thickness.

#### (168-0005) Universal Refill

Contains: 1 5ml QuickMix syringe and 10 Type 9 mixing tips.

#### Type 9 Mixing Tips Short (Brown/White)

(160-0015) Pkg/50



### Provicol® QM Aesthetic (VOCO)

Translucent zinc oxide, non-eugenol temporary cement with calcium hydroxide. Easy removal due to semi-elastic consistency. Low film thickness to enhance the fit and secure excellent retention. Ideal for the temporary cementation of temporary or permanent restorations.

#### (168-0010) QuickMix Syringe

Contains: 1 5ml QuickMix syringe and 10 Type 10 mixing tips.

#### Type 10 Mixing Tips, Long (Brown/White)

(160-0016) Pkg/50



### Provicol® QM Plus (VOCO)

Strong retention non-eugenol temporary cement with calcium hydroxide to reduce postoperative sensitivities. Increased strength for more retention and secure adhesion. Low film thickness allows easy placement for a perfect fit. Easy to clean up with peelable consistency. Non eugenol. Highly radiopaque.

#### (168-0009) QuickMix Syringe

Contains: 1 5ml QuickMix syringe and 10 Type 10 mixing tips.

#### Type 10 Mixing Tips Long (Brown/White)

(160-0016) Pkg/50



## Temporary Filling Materials

### Clip (VOCO)

Light-cured, provisional filling material that is non-sticky, flexible, and transparent. Easy to handle and remove without causing damage to preparation boundaries.

#### (168-0004) Caps Refill, pkg/25

Contains: 25 .25g caps.

#### (168-0002) Syringe Refill, pkg/2

Contains: 2 4g syringes.

#### (168-0003) Syringe Refill, pkg/3

Contains: 3 4g syringes.

### Clip Flow (VOCO)

Light-cured, flowable provisional filling material in Non-Dripping Technology (NDT) syringes for easy placement and minimal waste. Non-sticky, elastic, flowable, and durable. Easy to place and remove without causing damage to preparation boundaries.

#### (168-0001) Syringe Refill

Contains: 2 1.8g NDT syringes and 20 Type 46 intra-oral tips.

#### Type 46 Intra-oral Tips

(520-2653) Pkg/100



Non-stock

Non-stock



## Composites & Restoratives

### Universal Composites

### Admira® Fusion (VOCO)

Universal restorative material featuring Nano-ORMOCER® Technology (nano fillers, glass particles, and organically modified ceramics) to reduce shrinkage and increase biocompatibility. Non-sticky consistency for easy handling. High wear resistance.

#### Pkg/15 .2g Caps

(199-2000) A1	(199-2008) B2
(199-2001) A2	(199-2010) C2
(199-2002) A3	(199-2011) D3
(199-2004) A3.5	(199-2013) OA2
(199-2005) A4	(199-2016) Bleach
(199-2007) B1	



### Admira® Fusion x-tra Bulk Fill (VOCO)

A single shade, omni-chromatic restorative. Nano-ORMOCER® Technology (nano fillers, glass particles, and organically modified ceramics) reduces shrinkage and increases biocompatibility. BPA free.

#### (199-2020) 3g Syringe, each

(199-2021) Pkg/15 .2g Caps



### Grandio® Universal (VOCO)

A light-cured filling material featuring excellent radiopacity that cures under blue light. Contains inorganic fillers in a methacrylate matrix.

#### Pkg/20 .25g Caps

(199-1925) A1	(199-1938) B1
(199-1926) A2	(199-1933) C3
(199-1927) A3	(199-1935) D2
(199-1928) A3.5	(199-1937) Bleach
(199-1929) A4	



### Grandio®SO Universal (VOCO)

A universal, light-cured 89% filled nano-hybrid composite for anterior and posterior fillings of all classes. Outstanding wear resistance and enhanced color stability. Non-sticky. Excellent polishability and retention.

#### Pkg/16 .25g Caps

(199-1961) A1	(199-1968) B1
(199-1962) A2	(199-1969) B2
(199-1963) A3	(199-1971) C2
(199-1965) A3.5	(199-1975) OA2
(199-1966) A4	(199-1977) Bleach



### x-tra fil Bulk Fill (VOCO)

Quickly cured bulk-fill, packable posterior composite with low shrinkage due to multi-hybrid filler technology. Universal shade with high translucency and radiopacity. Ideal for fast core build-ups and Class I and II restorations.

#### (199-1903) Pkg/20 .25g Caps

### VisColor bulk (VOCO)

A highly-filled thermo-viscous bulk-fill composite. Initially flowable when warmed, then immediately sculptable as it cools. TVT polymer combined with super nano-hybrid fillers for fast viscosity change, flowability, low shrinkage, and excellent wear.

#### Caps Delivery

Contains: 16 .25g caps.

#### (199-2025) Universal

(199-2026) A1

(199-2027) A2

(160-0042) Capsule Warmer



Non-stock

## Flowable Composites

### Admira® Fusion Flow (VOCO)

A flowable, light-cured direct restorative material. Nano-ORMOCER® Technology (nano fillers, glass particles, and organically modified ceramics) reduces shrinkage up to 2.75% by volume. Excellent wettability.

#### Syringe Refill

Contains: 2 2g syringes and 20 Type 41 intra-oral tips.

(199-1866) A1 (199-1871) B1  
(199-1867) A2 (199-1872) B2  
(199-1868) A3 (199-1873) C2

(520-2651) Type 41 Tips, pkg/100



### Admira® Fusion x-base Bulk Fill (VOCO)

Flowable, universal-shade, bulk fill base restorative with Nano-ORMOCER® Technology (nano fillers, glass particles, and organically modified ceramics) to reduce shrinkage and inSmart self-leveling properties. Highly biocompatible and color stable.

#### (199-2022) Syringe Delivery

Contains: 2 2g syringes in universal shade and 20 Type 45 intra-oral tips.

(520-2652) Type 45 Tips, pkg/100



### Grandio® Flow (VOCO)

A light-curing, flowable restorative with innovative nano-technology. Radiopaque. Ideal for restorative therapy in anterior and posterior areas. Low viscosity for good wetting and adaptation to cavity walls.

#### Syringe Refill

Contains: 2 2g syringes and 20 Type 41 intra-oral tips.

(199-1801) A1 (199-1807) B1  
(199-1802) A2 (199-1808) OA2  
(199-1803) A3 (199-1809) BL

(520-2651) Type 41 Tips, pkg/100



### Grandio®SO Flow (VOCO)

Light-cured, flowable 81% filled nano-hybrid composite with outstanding strength, wear resistance, and color retention. Excellent wettability.

#### Syringe Refill

Contains: 2 2g syringes and 20 Type 41 intra-oral tips.

(199-1846) A1 (199-1849) A3.5  
(199-1847) A2 (199-1852) B1  
(199-1848) A3 (199-1858) White Opaque

#### (199-1860) Caps Delivery, A2

Contains: 16 .25g caps.

(520-2651) Type 41 Tips, pkg/100



### Grandio®SO Light Flow (VOCO)

Extremely flowable, yet highly stable nano-hybrid flowable composite with low viscosity, high performance and esthetics, and excellent handling.

#### (199-1818) Syringe Refill, A2

Contains: 2 2g syringes.

(520-2650) Type 44 Tips, pkg/100

(520-2652) Type 45 Tips, pkg/100



### Grandio®SO Heavy Flow (VOCO)

Highly viscous, light-cured, flowable 83% filled nano-hybrid composite. Features superior physical properties, strength, and wear resistance. Ideal for proximal and cervical areas. Non-Dripping Technology (NDT) syringe for easy placement and minimal waste.

#### Syringe Refill

Contains: 2 2g syringes and 20 Type 46 intra-oral tips.

(199-1826) A1 (199-1832) B1  
(199-1827) A2 (199-1836) Incisal  
(199-1828) A3 (199-1837) Bleach  
(199-1829) A3.5 (199-1838) White Opaque  
(199-1830) A4

#### Caps Delivery

Contains: 16 .25g caps.

(199-1839) A1 (199-1841) A3  
(199-1840) A2 (199-1842) A3.5

#### Type 46 Intra-oral Tips

(520-2653) Pkg/100



### x-tra base® Bulk Fill Flowable (VOCO)

Flowable bulk-fill posterior composite with low shrinkage, high strength and radiopacity, and smart self-leveling properties for excellent wall adaptation without slumping. Non-Dripping Technology (NDT) syringe for easy placement and minimal waste.

#### Syringe Refill

Contains: 2 2g syringes and 20 Type 46 intra-oral tips.

(199-1904) Universal  
(199-1905) A2

#### Caps Delivery

Contains: 16 .25g caps.

(199-1906) Universal  
(199-1907) A2

#### Type 46 Intra-oral Tips

(520-2653) Pkg/100



## Glass Ionomer Restoratives

### Ionolux® (VOCO)

Light-cured resin modified glass ionomer restorative material with composite-like esthetics, high fluoride release, and excellent translucency. Highly radiopaque. Ideal for Class III and V restorations, small Class I restorations, restorations of deciduous teeth, temporary restorations, core build-ups, and linings. No conditioner or adhesive necessary.

#### (199-1878) Intro Kit

Contains: 50 application capsules (10 shade A1, 30 shade A2, and 10 shade A3).

#### Capsule Refill, pkg/20

(199-1879) A1  
(199-1880) A2



### IonoStar® Plus (VOCO)

Fast-setting glass ionomer restorative material with high translucency, tooth-like fluorescence, and excellent physical properties for long-lasting restorations. Fluoride-releasing. Ideal for an abundance of restorations, as a base or liner, core build-up material, and extended fissure sealing. No conditioner required.

(199-1887) Capsule Refill A2, pkg/20



## Composites & Restoratives, continued

### Surface Sealant

#### Easy Glaze® (VOCO)

Nano-filled, light-cured liquid polish and surface sealant to finish temporary and permanent restorations in seconds, protecting against staining and discoloration. Features a natural gloss and tooth-like fluorescence. Can be used on all types of direct and indirect restorations. Works with LED, PAC, and halogen curing lights. Nano-filler for better wear resistance.

(080-0035) 5ml Bottle, each

Contains: 1 5ml bottle and 20 Easy Brushes.

(199-1893) Easy Brush, pkg/50



### Composite Dispenser

#### Capsule Dispenser (VOCO)

Low-force composite dispenser for restorative material in capsules.

(199-1900) Each



### Core Materials

#### Rebilda® DC (VOCO)

Dual-cure, core build-up composite in an automix delivery system. Results in tooth-like surface hardness. Features low heat development, even with large core volumes. High adhesion with system bond for fewer parapulpal pins. Highly radiopaque. Fluoride releasing.

##### Cartridge Delivery

Contains: 1 50g cartridge, 10 Type 12 mixing tips, and 10 Type 3 intra-oral tips.

(160-0029) Dentin

(160-0031) White

##### QuickMix Syringe Refills

Contains: 1 10g QuickMix syringe, 10 Type 11 mixing tips, and 10 Type 4 intra-oral tips.

(160-0033) Dentin

(160-0034) Blue

(160-0035) White

Type 3 Intra-oral Tips

(160-0040) Pkg/50

Type 4 Intra-oral Tips

(160-0014) Pkg/50

Type 12 Mixing Tips (Blue)

(160-0017) Pkg/50

Type 11 Mixing Tips

(520-2655) Pkg/50

Type 19, Blue, Bendable

(160-0041) Pkg/30



Non-stock

#### Grandio® Core Dual Cure (VOCO)

Flowable, dual-curing, nano-hybrid core build-up material with stackable consistency and dentin-like hardness. Highly radiopaque. Ideal for vital and non-vital teeth as well as the cementation of pins and posts. Automix syringe delivery.

##### QuickMix Syringe Refills

Contains: 1 10g QuickMix syringe, 10 Type 11 mixing tips, and 5 Type 4 intra-oral tips.

(160-0036) Dentin

(160-0037) Blue

(160-0038) White

Type 4 Intra-oral Tips

(160-0014) Pkg/50

Type 11 Mixing Tips

(520-2655) Pkg/50



## Temporary Crown & Bridge Materials

#### Structur® 2 SC (VOCO)

Fast setting temporary BIS-acrylic crown and bridge material that is strong and wear resistant. Sets fast in less than 1 minute. Provides long term color stability and composite-like polishability. Minimal shrinkage and low exothermic setting temperature. MMA-free.

##### Cartridge Delivery

Contains: 1 75g cartridge and 10 Type 6 mixing tips.

(168-0016) A1

(168-0022) A3.5

(168-0017) A2

(168-0019) B1

(168-0018) A3

(168-0020) BL

##### QuickMix Syringe Delivery

Contains: 1 8g syringe and 10 Type 10 mixing tips.

(168-0012) A2

(168-0013) A3

(160-0043) Type 6 Tips, pkg/50

(160-0016) Type 10 Tips, pkg/50



#### Structur® 3 (VOCO)

Nano-filled, quick setting, strong temporary crown and bridge material with "Wipe & Go" technology. Simply wipe with alcohol for a brilliant gloss—no polishing necessary. 45-second intra-oral setting time. Durable with high fracture resistance and compressive strength. Aesthetic with a natural gloss and tooth-like fluorescence. Color stable for long-term use.

##### Cartridge Delivery

Contains: 1 50ml cartridge and 10 Type 6 mixing tips.

(168-0034) A1

(168-0035) A2

(168-0036) A3

(160-0043) Type 6 Tips, pkg/50

#### Structur® Premium (VOCO)

Highly durable, fast setting temporary bis-acrylic crown and bridge material with added flexural strength. Ceramic-like esthetics with brilliant gloss and natural fluorescence. High-fracture resistance. Sets quickly (30-45 seconds).

##### Cartridge Delivery

Contains: 1 75g cartridge and 10 Type 6 mixing tips. Extruder gun sold separately.

(168-0026) A1

(168-0028) A3

(168-0027) A2

(168-0031) B1

(160-0043) Type 6 Tips, pkg/50





## Impression Materials

### Bite Registration Material

#### Registrado® X-tra (VOCO)

Extra-hard, fast, and flexible bite registration material with a non-dripping consistency. Self curing with a 40-second hard set and 90-second extreme final hardness. Thixotropic for easy application and no occlusion resistance.



##### (084-0004) Cartridge Delivery

Contains: 2 50ml cartridges, 10 Type 16 mixing tips, and 3 Type 2 intra-oral tips.

##### Type 2 Intra-oral Tips

(168-0047) Pkg/50

##### Type 16 Mixing Tips, Green

(168-0050) Pkg/50

Non-stock

### Vinyl Polysiloxane

#### V-Posil® VPS Impression System (VOCO)

VPS impression material with exceptional tear strength and elastic recovery. Up to 2 minutes of adjustable working time for control over the pace of the procedure. Abbreviated 2-minute thermo-activated set time improves patient comfort by saving time. Heavy Soft Fast material offers optimal viscosity for loading and remains soft after setting for easy removal.



##### 50ml Cartridge Refill

Contains: 2 50ml cartridges and 6 Type 16 mixing tips.

##### (528-0035) Heavy Soft Fast

(528-0040) Light Fast

(528-0042) X-Light Fast

##### 50ml 10-pack

Contains: 10 50ml cartridges.

##### (528-0036) Heavy Soft Fast

##### 380ml Dynamic Cartridge Refill

Contains: 1 380ml cartridge and 10 Type 31 mixing tips.

##### (528-0033) Heavy Soft Fast

##### 380ml Dynamic 2-pack

Contains: 2 380ml cartridges.

##### (528-0034) Heavy Soft Fast

(528-0032) Putty Fast

Non-stock

##### (528-0030) Putty Jar Fast

Contains: 1 450ml jar base, 1 450ml jar catalyst, 1 white spoon, and 1 blue spoon.

##### Type 3 Intra-oral Tips (for V-Posil Light Fast and X-Light Fast)

(160-0040) Pkg/50

##### Type 16 Mixing Tips, Green

(168-0050) Pkg/50

##### Type 31 Dynamic Mixing Tips with Bayonet Ring

(168-0052) Pkg/50

Non-stock

## Pins & Posts

### Posts

#### Rebilda® Post System (VOCO)

A complete post and core system with fiber reinforced resin posts. Highly translucent and radiopaque with dentin-like elasticity for gentle root treatment. High retention rate for strong bond and excellent light conduct and diffusion for use with dual-cured materials.



##### Post Refills, pkg/5

(648-2905) #12, 1.2mm, Green

(648-2906) #15, 1.5mm, Black

### Overdenture Products

#### Quick Up® (VOCO)

Tissue-colored implant pick-up system with new Quick Up method for faster procedure completion. Work without vent holes and excess material removal, under fill with 2.5 minute fast-setting, self-cured material, and finish with a light-cured component. Direct pick-up results in higher accuracy compared to indirect pick-up. Low-heat, self-curing material avoids tissue irritations. Tasteless and odorless. High bond strength to acrylic denture base.



##### (737-2800) System Kit

Contains: 1 7.5g QuickMix syringe of Quick Up, 1 2g syringe of Quick Up LC, 1 4ml adhesive, 1 5ml QuickMix syringe of Fit Test C&B, 15 Type 11 mixing tips, 10 Type 10 mixing tips, 15 Type 44 intra-oral tips, 15 Type 1 intra-oral tips, 2 mixing pads, and 20 single tips.

##### (737-2801) 7.5g QuickMix Syringe

Contains: 1 7.5g QuickMix syringe of Quick Up, and 10 Type 10 mixing tips.

##### (737-2802) 2g Syringe

Contains: 1 2g syringe of Quick Up LC, and 15 Type 44 intra-oral tips.

##### (737-2803) Adhesive, 4ml, each

##### Type 1 Intra-oral Mixing Tips

(160-0039) Pkg/50

##### Type 10 Mixing Tips, Long (Brown/White)

(160-0016) Pkg/50

##### Type 11 Mixing Tips

(520-2655) Pkg/50

##### Type 44 Intra-oral Mixing Tips


(520-2650) Pkg/100

(199-1895) Single Tim, pkg/100

## Retraction Systems

#### Retraction Paste (VOCO)

Astringent paste for the effective temporary retraction and moisture management of the gingival sulcus. Aluminum chloride-based. Initially flowable, then stable for optimal retraction. Turquoise blue color for high contrast. Easy and precise application.



(708-0001) Pkg/25 .3g Caps

(708-0002) Pkg/100 .3g Caps

## Preventives

### Fluoride Paste (Take-Home)

#### Remin Pro (VOCO)

**CDT code is 9910 for desensitizing medicament!**

Triple protective paste composed of fluoride, nano-hydroxyapatite, and Xylitol. No risk of milk protein allergies. Fast acting and great tasting.

Pkg/12 40g tubes

(452-0012) Mint

(452-0014) Assorted\*

\*Assorted contains: 4 40g tubes each of melon, mint, and strawberry.



### Fluoride Varnish

#### Profluorid® Varnish (VOCO)

**CDT code is 1206 for topical application of varnish, and CPT code is 99188 for medical claims!**

5% sodium fluoride varnish with rapid release in the first 15 minutes followed by continuous release for up to 4 hours. Thin, clear application and superior adherence to teeth.

(452-0017) .25ml SingleDose, Melon

Contains: 50 .25ml SingleDoses and 50 Easy Brushes.

Pkg/50 .4ml SingleDose Delivery

Contains: 50 .4ml SingleDoses and 50 Easy Brushes.

(452-0018) Melon (452-0036) Piña Colada

(452-0021) Caramel (452-0032) Bubblegum

(452-0022) Cherry (452-0034) Cola Lime

(452-0023) Mint

Pkg/200 .4ml SingleDose Delivery

Contains: 200 .4ml SingleDoses and 200 Easy Brushes.

(452-0019) Melon (452-0037) Piña Colada

(452-0024) Caramel (452-0033) Bubblegum

(452-0025) Cherry (452-0035) Cola Lime

(452-0026) Mint

(452-0020) Assorted SingleDose

Contains: 48 .4ml SingleDoses (12 each of cherry, caramel, melon and mint) and 48 Easy Brushes.

Pkg/10ml Tube, each

(452-0016) Melon (452-0029) Mint

(452-0027) Caramel (452-0031) Piña Colada

(452-0028) Cherry



## Pit & Fissure Sealant

#### Grandio® Seal (VOCO)

A 70% filled nano-hybrid pit and fissure sealant with low shrinkage, superb flow and wetting properties, and drip-free, bubble-free placement. Radiopaque. Non-Dripping Technology (NDT) syringe for easy placement and minimal waste.

(716-0001) Intro Kit

Contains: 5 2g NDT syringes of Grandio Seal, 1 5ml syringe of Vococid, 45 Type 45 intra-oral tips, and 10 Type 42 intra-oral tips.

(716-0002) NDT Syringe Refill

Contains: 2 2g NDT syringes and 15 Type 45 intra-oral tips.

Type 42 Mixing Tips

(199-1896) Pkg/30

Type 45 Intra-oral Tips

(520-2652) Pkg/100



## Desensitizer

#### Profluorid® L (VOCO)

**CDT code is 9910 for desensitizing medicament!**

Clear, 5% sodium fluoride liquid varnish composed of a synthetic resin matrix to ensure effective desensitizing during the bleaching process without interfering with results. Easy to apply. Free of glutaraldehyde, gluten, and methacrylate.

SingleDose Delivery

(080-0036) Pkg/50 with 50 single tims

SingleDose Delivery

(080-0037) Pkg/200 with 200 single tims

(080-0038) Bottle Delivery

Contains: 1 4.5ml bottle, 2 mixing palettes, and 20 Easy Brushes.

(199-1893) Easy Brush, pkg/50



# MidwestDental®

EQUIPMENT • SERVICES • SUPPLIES

(800) 766-2025 • [MidwestDental.com](http://MidwestDental.com)

Innovation  
Precision  
Excellence  
MADE IN GERMANY

LOG IN TO SEE YOUR PRICING

